

Ligation Protocol for Cloning with Instant Sticky-end Ligase Master Mix (M0370)

Protocol

1. Transfer master mix to ice prior to reaction set up. Mix tube by finger flicking before use.
2. Combine 20–100 ng of vector* with a 3-fold molar excess of insert and adjust volume to 5 µl with dH₂O.
3. Add 5 µl of Instant Sticky-end Ligase Master Mix, mix thoroughly by pipetting up and down 7-10 times, and place on ice. The sample is now ready to be used for transformation.
4. Use for transformation or store at -20°C.
5. **Do not heat inactivate.**

Heat inactivation dramatically reduces transformation efficiency.

* In-house testing has demonstrated that maximal transformation efficiency is achieved using between 20–100 ng of vector (sticky) and a corresponding 3-fold molar excess of the insert to be ligated into the vector.