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**A protocol of molecular detection of phytoplasmas and *Xylella* spp. in
post-entry quarantine for plants.**

supported by a Grant-in-Aid for “Regulatory research projects for food safety,
animal health and plant protection, no. 2508.” from the Ministry of Agriculture,
Forestry, and Fisheries of Japan.

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National Agriculture and Food Research Organization (NARO), Japan

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A protocol of molecular detection of phytoplasmas and *Xylella* spp. in post-entry quarantine for plants.

Version 1.2. (November 26th, 2018)

Takao Ito & Ryoji Nakaune
NARO, Japan

Introduction

This is a protocol for universal molecular detection of phytoplasmas (PP) and *Xylella* spp. (XL). The protocol shows TaqMan real-time and conventional PCR assays to detect various PP and XL strains from various plant samples. TaqMan real-time PCR can detect PP and/or XL with plant internal control (IC) simultaneously. We describe here the methods of Ito & Suzuki (2017) with some modifications in Japanese and English. This protocol is developed in the study supported by a Grant-in-Aid for “Regulatory research projects for food safety, animal health and plant protection, no. 2508.” from the Ministry of Agriculture, Forestry, and Fisheries of Japan.

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1. Preparation

Check the Appendix I for instruments, the Appendix II for reagents, and the Appendix III for the other information.

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2. Samples

Collect symptomatic leaves when pathogens are at high concentration. We describe here sampling from grapevines.

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2.1. Sampling

Collect some symptomatic leaves (Figs. 1 & 2) in late summer to early autumn. When no symptoms are observed, randomly collect some leaves attached to canes.

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Fig. 1. Examples of symptomatic leaves of grapevines infected with phytoplasmas

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Fig. 2. Examples of symptomatic leaves of grapevines infected with *Xylella fastidiosa*.

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2.2. Extraction

DNA can be extracted using commercial kits such as Nucleon Phytopure (GE Healthcare, Buckinghamshire, UK) and DNeasy Plant Mini Kit (Qiagen, Hilden, Germany), or using cetyltrimethylammonium bromide (CTAB)-based buffer (EPPO, 2016; IPPC, 2016). For simple diagnosis, we describe here a crude extraction modified from the method of Nakaune and Nakano (2006).

2.2.1. Crude extraction using the Multi-beads shocker (Yasui kikai, Osaka, Japan)

2.2.1.1. Put a metal bead and 1mL extraction buffer (Appendix II) into a 2mL tube, and cool them until use.

2.2.1.2. Collect leaf petioles (<50mg) of a sample and put them into the tube.

2.2.1.3. Crush samples in the Multi-beads shocker (2,500 rpm for 60 sec).

2.2.1.4. Centrifuge the tubes ($9,000 \times g$ for 10 min at 4°C), and transfer the supernatants to new tubes.

2.2.1.5. Next steps, or keep them in a freezer until use.

* The crude extracts can be used for detection of viruses and viroids (Nakaune and Nakano, 2006).

2.2.2. Isopropanol precipitation

2.2.2.1. Add an equal volume of cold isopropanol to the crude extract and mix (Some extracts can be mixed for a bulk assay).

2.2.2.2. Centrifuge the tubes ($20,000 \times g$ for 5 min at 4°C), and discard supernatants.

2.2.2.3. Dry pellets

2.2.2.4. Suspend the pellets in one-fifth volume of TE buffer.

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3. Molecular detection methods

3.1. TaqMan real-time PCR

The real-time PCR can detect targets without laborious post-amplification processing. The assay shows threshold cycle (Ct) values that are PCR cycle numbers at which a detectable amount of amplicon product has been generated. **Triplex real-time PCR assay to simultaneously detect PP, XL, and IC is described with red characters in this protocol.** **Duplex assay to simultaneously detect PP, and IC is described with blue characters.** **Duplex assay to simultaneously detect XL, and IC is described with green characters.** Singleplex assays to detect PP or XL only can be modified from the duplex assays. You can use the assays properly for different samples.

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Targets	Probes	Sequence (5' – 3') and modification
PP	UPH-Pb	TGACGGGACTCCGCACAAGCG (5' FAM-3' BHQ1)
IC	C18S-Pt	ACACACCGCCGTCGCTCC (5' TAMRA-3' BHQ2a)

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3.1.3. Primer stocks

You should utilize oligo-synthesizing companies supplying inosine (I)-containing primers at low prices. We use primers purified with reverse-phase cartridges. Primer stocks can be kept in a freezer.

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3.1.4. Probe stocks

Applied Biosystems (Foster City, CA, USA) and some others can supply TaqMan MGB probes (5' FAM-3' MGB/5' VIC-3' MGB). We have not tested 5' HEX-3' MGB primers which can be used instead of 5' VIC-3' MGB probes. We use probes purified with HPLC. Probe stocks should be kept in a freezer.

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3.1.6. Probe mixtures

We describe here probe mixtures from probe stocks of X, Y, and $100\mu\text{M}$ (pmol/μL). Keep the mixtures in a freezer until use.

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3.1.7. Setting of the StepOnePlus (Applied Biosystems)

3.1.7.1. Start the StepOnePlus software → [Advanced Setup] → [Quantitation–Standard Curve]/[TaqMan reagents]/[Fast] (default parameters) → [Run Method] → [Reaction Volume Per Well] should be (10μL), and thermal cycling parameters should be set as follows.

	Steps	Temp (°C)	Time	Cycles
1	Uracil–N–glycosylase incubation	50	2 min	1
2	Initial denaturation	95	20 sec	1
3	Denaturation	95	1 sec	50
	Primer annealing/extension (Data collection)	60	20 sec	

3.1.7.2. [Plate Setup] → [Define Targets and Samples] → [Define Targets]

3.1.7.2.1. Triplex real-time PCR

Target Name	Reporter	Quencher
PP	FAM	NFQ–MGB
XL	VIC	NFQ–MGB
IC	TAMRA	None

3.1.7.2.2. Duplex real-time PCR to detect PP/IC

Target Name	Reporter	Quencher
PP	FAM	None
IC	TAMRA	None

3.1.7.2.3. Duplex real-time PCR to detect XL/IC

Target Name	Reporter	Quencher
XL	FAM	NFQ–MGB
IC	TAMRA	None

3.1.7.3. [Define Samples] → Enter sample names.

3.1.7.4. [Assign Targets and Samples] → Assign sample names. Assign targets of PP/XL/IC for the triplex assay, targets of PP/IC for the duplex assays to detect PP/IC, and targets of XL/IC for the duplex assays to detect XL/IC.

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3.1.8. Reagent mixture

TaqMan FAST Advanced Master Mix (Applied Biosystems)

	Reagents	1 reaction	10 reactions
1	Sterile water	2.5	25
2	TaqMan FAST Advanced Master Mix	5.0	50
3	Primer mixture	1.0	10
4	Probe mixture	1.0	10
	Total	9.5 μ L	95 μ L

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3.1.9. TaqMan real-time PCR

3.1.9.1. Dispense 9.5 μ L of the reagent mixture to PCR tubes (wells of a plate).

3.1.9.2. Add 0.5 μ L of samples to the tubes (wells of a plate). Two reactions should be performed for each sample.

3.1.9.3. Set the tubes (the plate) in the StepOnePlus and [START RUN].

* Positive controls, healthy plant samples, and buffer samples should be tested for comparison in the same run.

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3.1.10. Data analysis(Figs. 3 & 4)

3.1.10.1. 「File」→「Export」→「Start Export」

3.1.10.2. In the exported file, Ct<45 can be considered positive for PP/XL and Ct<40 can be considered positive for IC.

3.1.10.3. PP (XL) should be positive in positive controls. PP (XL) and IC should be negative and positive in healthy plant samples, respectively. All the targets should be negative for buffer

samples.

3.1.10.4. No positives from plant samples show PCR errors.

3.1.10.5. PP (XL) positives in the two reactions can show PP (XL) in the samples. If the outcome is critical (e.g. post-entry quarantine samples, new host records, and new distributions) or (sub)species must be determined, direct sequencing should be done (3.1.10.7).

3.1.10.6. Only one PP (XL) positive in the two reactions may show irregular non-specific amplification. Additional two reactions should be performed. To confirm PP positives, a loop-mediated isothermal amplification (LAMP) assay with a Phytoplasma universal detection kit (Nippon Gene, Tokyo, Japan) can be also utilized.

3.1.10.7. Direct sequencing of the PCR products (approximately 170bp)

3.1.10.7.1. Purify PCR products using reagents such as QIAquick PCR purification kit (Qiagen) or ExoSAP-IT (Thermo Fisher Scientific, Waltham, MA, USA)

3.1.10.7.2. Directly sequence the products with each of the primers UPH-F and D-UPHr2 for PP, and each of the primers XrDf1 and D-XrDr2 for XL.

3.1.10.7.3. Search nucleotide sequence identities using BLASTN (<http://ddbj.nig.ac.jp/blast/blastn?lang=ja>, <http://blast.ncbi.nlm.nih.gov/Blast.cgi>) or [S1 and S2 Datasets] (Appendix III) of the products (~120bp) after trimming away primer regions. Products amplified from PP and XL can show the highest nucleotide sequence identities with PP and XL, respectively (Fig. 5).

* When you cannot obtain nucleotide sequence data from heterogeneous or low-amount products, use templates diluted, concentrated, or purely extracted, reduce or increase PCR cycles, perform singleplex PCR, *etc.*

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Fig. 3. Amplification plots of TaqMan triplex real-time PCR to detect PP/XL/IC

Well	Sample Name	Target Name	Reporter	Quencher	Cr	Threshold	Automatic Ct Threshold
G4	X-PD	IC	TAMRA	None	Undetermined	0.011243913	TRUE
G4	X-PD	PP	FAM	NFQ-MGB	Undetermined	0.024366653	TRUE
G4	X-PD	XL	VIC	NFQ-MGB	27.08343506	0.042214422	TRUE
H4	X-PD	IC	TAMRA	None	Undetermined	0.011243913	TRUE
H4	X-PD	PP	FAM	NFQ-MGB	Undetermined	0.024366653	TRUE
H4	X-PD	XL	VIC	NFQ-MGB	27.09623337	0.042214422	TRUE
G11	Tp2	IC	TAMRA	None	19.06248665	0.011243913	TRUE
G11	Tp2	PP	FAM	NFQ-MGB	23.83968163	0.024366653	TRUE
G11	Tp2	XL	VIC	NFQ-MGB	Undetermined	0.042214422	TRUE
H11	Tp2	IC	TAMRA	None	19.20702171	0.011243913	TRUE
H11	Tp2	PP	FAM	NFQ-MGB	23.78330994	0.024366653	TRUE
H11	Tp2	XL	VIC	NFQ-MGB	Undetermined	0.042214422	TRUE
H12	PX	IC	TAMRA	None	23.6335907	0.011243913	TRUE
H12	PX	PP	FAM	NFQ-MGB	20.76914597	0.024366653	TRUE
H12	PX	XL	VIC	NFQ-MGB	27.72828865	0.042214422	TRUE
G7	P2	IC	TAMRA	None	21.13648224	0.011243913	TRUE
G7	P2	PP	FAM	NFQ-MGB	Undetermined	0.024366653	TRUE
G7	P2	XL	VIC	NFQ-MGB	Undetermined	0.042214422	TRUE
H7	P2	IC	TAMRA	None	21.51067734	0.011243913	TRUE
H7	P2	PP	FAM	NFQ-MGB	Undetermined	0.024366653	TRUE
H7	P2	XL	VIC	NFQ-MGB	Undetermined	0.042214422	TRUE
A12	B	IC	TAMRA	None	Undetermined	0.011243913	TRUE
A12	B	PP	FAM	NFQ-MGB	Undetermined	0.024366653	TRUE
A12	B	XL	VIC	NFQ-MGB	Undetermined	0.042214422	TRUE

Fig. 4. An exported result.

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3.2. Conventional PCR

Conventional PCR can be used without real-time PCR instruments, although it necessitates gel electrophoresis of the amplified products. Each of PCR assays targeting PP and XL should be performed in a singleplex assay. Although they can be performed in a duplex assay, you cannot judge which of the targets is amplified.

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3.2.2. Primer stocks

You should utilize oligo-synthesizing companies supplying inosine (I)-containing primers at low prices. We use primers purified with reverse-phase cartridges. Keep primer stocks in a freezer until use.

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3.2.3. Primer mixtures

Mix primer stocks of PP or XL at 25 μ M (pmol/ μ L) each. Keep them in a freezer.

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3.2.5. PCR

3.2.5.1. Dispense 9.5μL of the reagent mixture to PCR tubes.

3.2.5.2. Add 0.5μL of samples to the tubes. If using a thermal cycler without a heated lid, overlay the reaction mix with mineral oil to prevent evaporation during thermal cycling. Cap the tubes and centrifuge briefly

3.2.5.3. Set the tubes on the PCR instrument.

3.2.5.4. Start thermal cycling shown below.

	Steps	Temp. (°C)	Time	Cycles
1	Initial denaturation	95	2 min	1
2	Denaturation	95	20 sec	43
	Primer annealing	55	20 sec	
	Primer extension	72	20 sec	
3	Final extension	72	5 min	1
4	Soak	4	∞	

* You can use any thermal cycling parameters which detect positive controls.

3.2.5.5. Electrophoreses the amplified fragments (with dye) in a 1.5–2% agarose gel and visualize them with ethidium bromide staining (Fig. 6).

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Fig. 6. 1.5% agarose gel electrophoresis of amplified fragments of XL.

M: DNA size marker, 1–7: other bacteria samples, 8–15: XL samples, 16: TE buffer

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4. Quick start guide

Triplex real-time PCR to detect phytoplasmas (PP), *Xylella* spp. (XL), and plant internal controls (IC) from crude extracts

1	Extraction	1.1. Crude extraction 1.1.1. Put leaf petioles (<50mg), a metal beads, and 1mL extraction buffer into a tube. 1.1.2. 2,500 rpm 60 sec. (the Multi-beads shocker) 1.1.3. 9,000 x g 10min 4°C 1.1.4. Transfer the supernatant to a new tube. Next steps, or keep it in a freezer. 1.2. Isopropanol precipitation 1.2.1. Add an equal volume of cold isopropanol to the crude extract and mix. 1.2.2. 20,000 x g 5 min 4°C 1.2.3. Discard supernatants and dry pellets. 1.2.4. Suspend the pellet in one-fifth volume of TE. Next steps, or keep it in a freezer.																																																																
2	Real-time PCR	2.1. Reagent mixture <table><tr><th>Reagents</th><th>1 reaction</th><th>10 reactions</th></tr><tr><td>Sterile water</td><td>2.5</td><td>25</td></tr><tr><td>TaqMan FAST Advanced Master Mix</td><td>5</td><td>50</td></tr><tr><td>Primer mixture</td><td>1</td><td>10</td></tr><tr><td>Probe mixture</td><td>1</td><td>10</td></tr><tr><td>Total (µL)</td><td>9.5</td><td>95</td></tr></table> 2.2. Dispense 9.5 µL of the reagent mixture to PCR tubes 2.3. Add 0.5 µL of the extract (1.2.4) to the tube. 2.4. Set the tubes and run the StepOnePlus with the following parameters: 50°C 2 min. → 95°C 20 sec. → 95°C 1 sec. → 60°C 20 sec. (50 cycles) * Targets (Reporter/Quencher) : PP (FAM/NFQ-MGB) 、XL (VIC/NFQ-MGB) 、IC (TAMRA/None)	Reagents	1 reaction	10 reactions	Sterile water	2.5	25	TaqMan FAST Advanced Master Mix	5	50	Primer mixture	1	10	Probe mixture	1	10	Total (µL)	9.5	95																																														
Reagents	1 reaction	10 reactions																																																																
Sterile water	2.5	25																																																																
TaqMan FAST Advanced Master Mix	5	50																																																																
Primer mixture	1	10																																																																
Probe mixture	1	10																																																																
Total (µL)	9.5	95																																																																
3	Data analysis	3.1. Export data 3.2. Consider positives of PP/XL at Ct<45 and IC at Ct<40. <table><tr><th>Well</th><th>Sample Name</th><th>Target Name</th><th>Ct</th></tr><tr><td>G4</td><td>X-PD</td><td>IC</td><td>Undetermined</td></tr><tr><td>G4</td><td>X-PD</td><td>PP</td><td>Undetermined</td></tr><tr><td>G4</td><td>X-PD</td><td>XL</td><td>27.08343506</td></tr><tr><td>H11</td><td>Tp2</td><td>IC</td><td>19.20702171</td></tr><tr><td>H11</td><td>Tp2</td><td>PP</td><td>23.78330994</td></tr><tr><td>H11</td><td>Tp2</td><td>XL</td><td>Undetermined</td></tr><tr><td>H12</td><td>PX</td><td>IC</td><td>23.6335907</td></tr><tr><td>H12</td><td>PX</td><td>PP</td><td>20.76914597</td></tr><tr><td>H12</td><td>PX</td><td>XL</td><td>27.72828865</td></tr><tr><td>G7</td><td>P2</td><td>IC</td><td>21.13648224</td></tr><tr><td>G7</td><td>P2</td><td>PP</td><td>Undetermined</td></tr><tr><td>G7</td><td>P2</td><td>XL</td><td>Undetermined</td></tr><tr><td>A12</td><td>B</td><td>IC</td><td>Undetermined</td></tr><tr><td>A12</td><td>B</td><td>PP</td><td>Undetermined</td></tr><tr><td>A12</td><td>B</td><td>XL</td><td>Undetermined</td></tr></table>	Well	Sample Name	Target Name	Ct	G4	X-PD	IC	Undetermined	G4	X-PD	PP	Undetermined	G4	X-PD	XL	27.08343506	H11	Tp2	IC	19.20702171	H11	Tp2	PP	23.78330994	H11	Tp2	XL	Undetermined	H12	PX	IC	23.6335907	H12	PX	PP	20.76914597	H12	PX	XL	27.72828865	G7	P2	IC	21.13648224	G7	P2	PP	Undetermined	G7	P2	XL	Undetermined	A12	B	IC	Undetermined	A12	B	PP	Undetermined	A12	B	XL	Undetermined
Well	Sample Name	Target Name	Ct																																																															
G4	X-PD	IC	Undetermined																																																															
G4	X-PD	PP	Undetermined																																																															
G4	X-PD	XL	27.08343506																																																															
H11	Tp2	IC	19.20702171																																																															
H11	Tp2	PP	23.78330994																																																															
H11	Tp2	XL	Undetermined																																																															
H12	PX	IC	23.6335907																																																															
H12	PX	PP	20.76914597																																																															
H12	PX	XL	27.72828865																																																															
G7	P2	IC	21.13648224																																																															
G7	P2	PP	Undetermined																																																															
G7	P2	XL	Undetermined																																																															
A12	B	IC	Undetermined																																																															
A12	B	PP	Undetermined																																																															
A12	B	XL	Undetermined																																																															

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6.1.1. Cell disruptors

We use the Multi-beads shocker (Yasui kikai, Osaka, Japan). You can use any other tools which can crush plant tissues (Nakaune & Nakano, 2006).

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6.1.2. Centrifuges

We use refrigerated centrifuges with maximum 20,000 × g (15,000 rpm).

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6.1.3. Real-time PCR instruments

We developed this protocol to use the StepOnePlus system (Applied Biosystems). You can use any other real-time PCR instruments, for which the protocol may have to be modified.

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6.1.4. Thermal cyclers (for conventional PCR)

We use GeneAmp 9700 (Applied Biosystems). You can use any other thermal cyclers.

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6.1.5. Electrophoresis apparatuses (for conventional PCR)

You can run PCR products in agarose gels in any electrophoresis apparatuses

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6.1.6. Detectors (for conventional PCR)

UV-transilluminator and photo equipment, *etc.*

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6.1.7. Tips

We use filter tips to avoid cross-contamination.

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6.1.8. Tubes

We use 2mL tubes (ST-0250F) / metal beads (MC-0218R(S)) (Yasui kikai) for the Multi-beads shocker, 1.5mL tubes, optical PCR tubes for real-time PCR, and PCR tubes for conventional

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6.1.9. Support bases for optical PCR tubes for real-time PCR.

We use support bases for optical PCR tubes (plates) to keep the tubes clean.

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6.2.2. Reagents for real-time PCR

- ・TaqMan FAST Advanced Master Mix (Applied Biosystems)
 - * We developed this protocol to use this reagent.
 - * Inosine-containing primers should not be used generally with proofreading PCR enzymes, as described below. Anyway, you can use other reagents, for which the protocol may have to be modified.

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6.2.3. Reagents for conventional PCR

•GoTaq Hot Start DNA polymerase (Promega)

* Inosine-containing primers should not be used generally with proofreading PCR enzymes having 3' → 5' exonuclease activity. However, not only the non-proofreading GoTaq polymerase but also a proofreading TaKaRa *Ex Taq* (Takara, Shiga, Japan) can work well in the assay. Preliminary tests should be done for any other polymerases.

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6.2.4. Primers and probes

[3. Molecular detection methods] section describes them in detail. Ito & Suzuki (2017) shows sequence alignments of partial 16S rDNA genes of Phytoplasmas and *Xylella* spp. with the primers and probes in S1 and S2 Figs. If new species or subspecies are found, nucleotide sequences of their 16S rDNA should be compared with those of the primers and probes in advance. Preliminary tests to detect target sequences of the new (sub)species, which can be synthesized using gBlocks of Integrated DNA Technologies (Coralville, IA, USA), are also recommended.

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6.2.5. Extraction buffer

You can use this extraction buffer modified from that of Ito & Suzuki (2017).

Reagents	Final conc.	Amount per 1L
Trisodium citrate dihydrate	50 mM	14.7 g
PVP K25 or K30	2 % (w/v)	20 g

- Dissolve them in 1L distilled water.
- Sterile the buffer with an autoclave, and keep it in a refrigerator for storage.
- Just before use, add DTT at final concentration of 10 mM.

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6.3.1. Datasets of nucleotide sequences of PP/XL.

You can use nucleotide sequence files (*.seq) of PP/XL in S1 and S2 Datasets (Ito & Suzuki, 2017) for nucleotide sequence identity searches of the PCR products after trimming away primer regions, instead of the BLASTN search. The Datasets must be updated by yourself.

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6.3.2. Positive control plants

The positive control should be approximately one-tenth of the amount of leaf tissue used per plant for the DNA extraction (IPPC, 2016).

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6.3.3. Positive control DNA of PP

You can use synthesized DNA as positive controls. For example, the nucleotide sequence (M30790–partial16SrDNA, 480bp) can be synthesized using gBlocks of Integrated DNA Technologies (Coralville, IA, USA). Synthesized DNA should be diluted at 20 copies (10 cells: one PP cell has two 16S rDNA copies)/0.5μL with healthy plant extracts.

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6.3.4. Positive control DNA of XL

You can use synthesized DNA as positive controls. For example, the nucleotide sequence (AF192343–partial16SrDNA, 480bp) can be synthesized using gBlocks of Integrated DNA Technologies. Synthesized DNA should be diluted at 20 copies (10 cells: one XL cell has two 16S rDNA copies)/0.5μL with healthy plant extracts.

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6.3.5. Plants

Collect symptomatic leaves when pathogens are at high concentration in sampling from any plants. Symptoms and concentration of the bacteria in a plant can depend upon environmental factors, strains, and the host plant species (EPPO, 2106; IPPC, 2016).

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