

Biohazard Garbage Protocol



Autoclave Biohazard Waste

1. Locate Biohazard garbage
2. If the waste container is $> \frac{3}{4}$ full, then it needs to be processed for autoclaving.
 - a. Autoclave the processed waste on gravity cycle.
 - b. Instructions on cycle times based on the amount of trash you are autoclaving are on the autoclave
3. Instructions on how to process **non-sharps** waste for autoclaving:
 - a. Transfer bag into a wide autoclavable tray.
 - i. Tray will have PP on bottom. Don't use HDPE! It will MELT!
 - b. Pour 100ml di water into bag containing trash.
 - c. Pour 100ml di water into tray (Yep! Between the container and the bag!)
 - d. Close bag loosely put a strip of autoclave tape on bag.
 - e. Write "Wilhelm" on the autoclave tape.
4. Instructions on how to process **sharps** waste for autoclaving:
 - a. Transfer sharps container into a wide autoclavable tray.
 - i. Tray will have PP on bottom. Don't use HDPE! It will MELT!
 - b. Pour 100 ml water into interior of the sharps container.
 - c. Close lid until it snaps shut, and put autoclave tape on lid.
 - d. Write "Wilhelm" on the autoclave tape.
5. After autoclaved waste has cooled, put the autoclaved waste into designated trash bin next to autoclave.