

## Cell Hashing protocol

For experiments involving cell hashing, we recommend using the [cost per cell calculator](#) from the Satija lab to plan experiments, determine number of hashes, number of cells to load, expected doublet rates (detected and undetected) and cost considerations:

### Cell staining for Drop-seq or 10x Genomics

- Obtain all single cell suspensions from different samples/conditions that will be multiplexed in the run. Keep samples in separate tubes until after cell hashing and shortly before loading cells into the single cell RNA-seq instrument. When aiming to super-load the same sample into one run, divide the sample up into equal proportions before staining with distinct cell hashing antibodies. Keep cell suspensions on ice (unless otherwise stated) at all times.
- Carefully count all cells to ensure accurate quantitation.
  - Make note of cell viability (>95%) and also include dead cells in the total cell count!
  - If you observe many dead cells, live cell enrichment (e.g. by FACS) is recommended!
- Resuspend ~1-2 million cells in 100 µl Staining buffer (2%BSA/0.02%Tween, PBS).
- Add 10 µl Fc Blocking reagent (FcX, BioLegend).
- Incubate for 10 minutes at 4°C.
- Add 1 µg of single cell hashing antibody to cells.
- Incubate for 30 minutes at 4°C.
- Wash cells 3 times with 1 mL Staining buffer (2%BSA/0.02%Tween, PBS), spin 5 minutes 400g at 4°C.
- Resuspend cells in PBS at appropriate concentration for downstream application.  
(e.g. for **10x** ~500 cells/µl; for **Drop-seq** [~200 cells/µl]; for **super-loading** ~1,500 cells/µl or higher).
- Filter cells through 40 µm strainers (e.g. Flowmi cell strainer).
- **Verify cell concentration by counting on hemocytometer after filtration.**
- **Pool all different samples/conditions at desired proportions and immediately proceed to next step.**

Run **Drop-seq** ([Macosko et al., 2015](#)) or **10x Genomics single cell 3' v2** assay as described until before cDNA amplification.

At cDNA amplification step:

**Add “additive” primer to cDNA PCR to increase yield of HTO products:**

**HTO PCR additive primer (2 µM): 1 µl (for 10x Genomics) or 0.4 µl (for Drop-seq)**

Subtract the total volume of additive primer from the water added to the PCR reaction.

**After cDNA amplification: Separate HTO-derived cDNAs (<180bp) and mRNA-derived cDNAs (>300bp).**

- Perform SPRI selection to separate mRNA-derived and antibody-oligo-derived cDNAs.
- **DO NOT DISCARD SUPERNATANT FROM 0.6X SPRI. THIS CONTAINS THE HASHTAGS.**
  - Add 0.6X SPRI to cDNA reaction as described in 10x Genomics or Drop-seq protocol.
  - Incubate 5 minutes and place on magnet.
  - **Supernatant contains hashtags.**
  - Beads contain full length mRNA-derived cDNAs.
- **mRNA-derived cDNA >300bp (beads fraction).**
  - Proceed with standard 10x or Drop-seq protocol for cDNA sequencing library preparation.
- **Hashtags <180bp (supernatant fraction).**
  - Purify Hashtags using two 2X SPRI purifications per manufacturer protocol:
    - Add 1.4X SPRI to supernatant to obtain a final SPRI volume of 2X SPRI.
    - Transfer entire volume into a low-bind 1.5mL tube.
    - Incubate 10 minutes at room temperature.
    - Place tube on magnet and wait ~2 minutes until solution is clear.
    - Carefully remove and discard the supernatant.


- Add 400  $\mu$ l 80% Ethanol to the tube without disturbing the pellet and stand for 30 seconds (only one Ethanol wash).
- Carefully remove and discard the ethanol wash.
- Centrifuge tube briefly and return it to magnet.
- Remove and discard any remaining ethanol.
- Resuspend in beads in 50  $\mu$ l water.
- Perform another round of 2X SPRI purification by adding 100  $\mu$ l SPRI reagent directly onto resuspended beads.
- Mix by pipetting, and incubate 10 minutes at room temperature.
- Place tube on magnet and wait ~2 minutes until solution is clear.
- Carefully remove and discard the supernatant.
- Add 200  $\mu$ l 80% Ethanol to the tube without disturbing the pellet and stand for 30 seconds (1<sup>st</sup> Ethanol wash).
- Carefully remove and discard the ethanol wash.
- Add 200  $\mu$ l 80% Ethanol to the tube without disturbing the pellet and stand for 30 seconds (2<sup>nd</sup> Ethanol wash).
- Carefully remove and discard the ethanol wash.
- Centrifuge tube briefly and return it to magnet.
- Remove and discard any remaining ethanol and allow the beads to air dry for 2 minutes (do not over dry beads).
- Resuspend beads in 45  $\mu$ l water.
- Pipette mix vigorously and incubate at room temperature for 5 minutes.
- Place tube on magnet and transfer clear supernatant into two PCR tubes.
- Amplify HTO sequencing library:
  - Prepare 100uL PCR reaction with purified small fraction:
    - 45  $\mu$ l purified Hashtag fraction
    - 50  $\mu$ l 2x KAPA Hifi PCR Master Mix.
    - 2.5  $\mu$ l TruSeq DNA D7xx\_s primer (containing i7 index) 10  $\mu$ M.
    - 2.5  $\mu$ l P5 oligo at 10  $\mu$ M depending on application:
      - For Drop-seq use P5-SMART-PCR hybrid oligo.
      - For 10x use SI PCR oligo.
    - Cycling conditions:
 

|      |        |  |               |
|------|--------|--|---------------|
| 95°C | 3 min  |  |               |
| 95°C | 20 sec |  |               |
| 64°C | 30 sec |  | ~ 8-12 cycles |
| 72°C | 20 sec |  |               |
| 72°C | 5 min  |  |               |
- Purify PCR product using 1.6X SPRI purification by adding 160  $\mu$ l SPRI reagent.
  - Incubate 5 minutes at room temperature.
  - Place tube on magnet and wait 1 minute until solution is clear.
  - Carefully remove and discard the supernatant.
  - Add 200  $\mu$ l 80% Ethanol to the tube without disturbing the pellet and stand for 30 seconds (first Ethanol wash).
  - Carefully remove and discard the ethanol wash.
  - Add 200  $\mu$ l 80% Ethanol to the tube without disturbing the pellet and stand for 30 seconds (second Ethanol wash).
  - Carefully remove and discard the ethanol wash.
  - Centrifuge tube briefly and return it to magnet.
  - Remove and discard any remaining ethanol and allow the beads to air dry for 2 minutes.
  - Resuspend beads in 20  $\mu$ l water.
  - Pipette mix vigorously and incubate at room temperature for 5 minutes.
  - Place tube on magnet and transfer clear supernatant to PCR tube.
- Hashtag libraries are now ready to be sequenced.
  - Quantify library by standard methods (QuBit, BioAnalyzer, qPCR).
  - Hashtag library will be around 180 bp (Figure 1).

## Sequencing Cell Hashing libraries:

- We estimate that an average of 100 molecules of HTO per cell is sufficient to achieve useful information. The number of reads required to obtain 100 molecules depends on the complexity of the sequencing library (e.g. duplication rate). HTO and cDNA sequencing libraries can be pooled at desired proportions. To obtain sufficient read coverage for both libraries we typically sequence HTO libraries in 5-10% of a lane and cDNA library fraction at 90% of a lane (HiSeq2500 Rapid Run Mode Flowcell).

### HTO library structure:



### Read 1:


Cell Barcode UMI TTTTTTTTTT ...

### Read 2:

Antibody Barcode B AAAAAAAAAAAAAAAAAA ...

## Figures

Figure 1.



**Figure 1. Hashtag library verification.** (a) A TSO-RT-oligo product (~140 bp) can be amplified during the HTO PCR by carryover primers from cDNA amplification. The product will not cluster but will interfere with quantification. This example figure shows ADT libraries. Sequential 2X SPRI purification of the HTO fraction after cDNA amplification reduces carryover of primers from cDNA amplification, and minimizes the amplification of this product during HTO-library amplification. To further enrich for HTO specific product the purified HTO library can be reamplified for ~3 additional cycles with HTO specific primer sets or P5/P7 generic primers. (b) A clean HTO library will contain a predominant single peak at around 180 bp.

## Oligonucleotide sequences:

### Hashtag oligos (HTOs):

These contain standard TruSeq DNA read 2 sequences and can be amplified using truncated versions of Illumina's TruSeq DNA primer sets (see example D701\_s below). See example below with a 12nt barcode:

5' GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTAGGACCATCCAABAAAAAAAAAAAAAAAAAAAAAAAAAAAAA\*AA

### Oligos required for HTO library amplification:

- Drop-seq P5-SMART-PCR hybrid primer (for Drop-seq only)  
5' AATGATACGGCGACCACCGAGATCTACACGCTGTCCGCGGAAGCAGTGGTATCAACGCAGAGT\*A\*C
- 10x Genomics SI-PCR primer (for 10x Single Cell Version 2 only)  
5' AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGC\*T\*C
- HTO cDNA PCR additive primer  
5' GTGACTGGAGTTCAGACGTGTGC\*T\*C
- Illumina TruSeq D701\_s primer (for HTO amplification; i7 index 1, shorter than the original Illumina sequence)  
5' CAAGCAGAAGACGGCATACGAGATCGAGTAATGTGACTGGAGTTCAGACGTGT\*G\*C

\* Phosphorothioate bond  
B C or G or T; not A nucleotide

**Materials and Kits needed:**

- FC blocking reagent (e.g. BioLegend FcX)
- Desalting columns (e.g. Bio-Rad, Micro Bio-Spin 6 Columns. #732-6221)
- 8-strip PCR tubes, **emulsion safe (!)** (e.g. TempAssure PCR 8-strips, USA Scientific)
- Bioanalyzer chips and reagents (DNA High Sensitivity and small RNA kit, Agilent)
- SPRIselect reagent (GE Healthcare, B23317)
- E-gel 4% (Invitrogen, USA)
- Low-bind 1.5 mL tubes (Common lab suppliers)
- PCR Thermocycler (e.g. Bio-Rad, T100)
- Magnetic tube rack (e.g. Invitrogen, USA)
- Qubit (Invitrogen, USA)
- Hemocytometer (e.g. Fuchs Rosenthal)
- DMSO (Common lab suppliers).
- PBS (Common lab suppliers)
- Tween20 (Common lab suppliers)
- Biotin (Common lab suppliers)
- TE pH 8.0 (Common lab suppliers)
- BSA (Common lab suppliers)