

1. Objective

The isolation of lymphocytes and pan-mononuclear cells from human whole blood for scRNAseq analysis.

2. Materials

50mL Centrifuge Tube (Fisher Scientific, Cat. No.: 12-565-271)
5mL Falcon™ Round-Bottom Polypropylene Tubes (Fisher Scientific, Cat. No.: 14-959-11A)
Biotin anti-human CD235ab Antibody (100µg) (Biolegend, Cat. No.: 306618)
Biotin anti-human CD66b Antibody (100µg) (Biolegend, Cat. No.: 305120)
BioMag®Plus Streptavidin (Bangs Laboratories, Cat. No.: BP628)
Cryogenic Vials (Fisher Scientific, Cat. No.: 09-761-71)
Cryostor CS10 (Fisher Scientific, Cat. No.: NC9930384)
Dulbecco's Phosphate Buffered Saline (DPBS) (Fisher Scientific, Cat. No.: 14-190-144)
EDTA 0.5M pH 8.0 (Fisher Scientific, Cat. No.: 15-575-020)
Fetal Bovine Serum (FBS) (Fisher Scientific, Cat. No.: 10-099-14)
Ficoll-Paque™ PLUS Media (Fisher Scientific, Cat. No.: 45-001-749)
Human TruStain FcX (Biolegend, Cat. No.: 422302)
NC-Slide A8 (Chemometec, Cat. No.: 942-0003)
Penicillin-Streptomycin-Glutamine (100X) (Fisher Scientific, Cat. No.: 10-378-016)
Solution 13 AO/DAPI (Chemometec, Cat. No.: 910-3013)
Dead Cell Removal Kit (Miltenyi, Cat. No.: 130-090-101)
 Dead Cell Removal Microbeads
 20x Binding Buffer Stock Solution
MS Columns (Miltenyi, Cat. No.: 130-042-101)
RosetteSep Granulocyte Depletion Cocktail (Stemcell Technologies, Cat. No.: 15624)

3. Equipment

Centrifuge
Cell Counter - NC-3000
EasyEights™ EasySep™ Magnet (Stemcell Technologies, Cat. No.: 18103)
MACS Multistand (Miltenyi, Cat. No.: 130-108-934)

4. Protocol

4.1.1. Preparing Buffer

- 4.1.2. Create the following **DPBS Solution-EDTA** in a bottle of DPBS by using the table below:

<i>Component</i>	<i>Volume (mL)</i>	<i>Starting Conc.</i>	<i>Final Conc.</i>
DPBS	474	-	-
FBS	25	100%	5%
EDTA	1	0.5M	1mM

Table 1.

4.2. Biosafety Notes

- All materials required for sample processing are to be prepared and in the biosafety cabinet before handling of blood samples.
- All blood sample manipulation takes place in a biosafety cabinet unless specifically stated.
- All centrifugation steps must take place in capped containers. Upon completion of a centrifugation step, return the capped container to the biosafety cabinet remove all tubes, and spray and wipe them with >70% ethanol before continuing to the next step.

4.3. Preparation of Blood

- 4.3.1. Place sample box into biosafety cabinet.
- 4.3.2. Remove samples from box and from the containment bags, discard bags, spray sample containers with >70% ethanol and wipe down.
- 4.3.3. Record the total volume of whole blood to be processed below:
_____ mL
- 4.3.4. Spin the whole blood 400 x g for 10 minutes, 20°C in the anti-coagulant tubes and remove the plasma layer to cryovials. Be sure to not remove any red blood cells with the plasma. **Add the same volume to all cryovials.**
Record number of vials: _____ and the volume per vial _____ mL.
NOTE: Ensure that cryovials are decontaminated prior to removal from the biosafety cabinet.
- 4.3.5. Add the same volume of **DPBS Solution-EDTA** as removed plasma back to the blood tube. Pipette to mix.

4.4. RosetteSep and Ficoll-Paque

- 4.4.1. Add 50µL of RosetteSep Cocktail/1mL of blood directly to the blood tube. Pipette to mix
- 4.4.2. Incubate samples at room temperature for 20 minutes.
- 4.4.3. Transfer sample to a 50mL tube and dilute up to 25ml with **DPBS Solution-EDTA solution.**
- 4.4.4. Aliquot 15mL of Ficoll-Paque Media PLUS to a 50mL tube.
- 4.4.5. Using the **slow** setting on the pipette gun, **gently** layer the blood/DPBS Solution-EDTA mixture on top of the 15mL of Ficoll-Paque Media PLUS. **Take extra care not to disturb the blood-ficoll**

interface while layering. Disturbing the interface excessively prevents the mononuclear cells from becoming a clean layer.

- 4.4.6. Spin for 20 minutes, 1200 x g at 20°C with **no brake**, 4 acceleration.

Note: Centrifuge should be pre-warmed to 20°C.

- 4.4.7. Remove the mononuclear cell layer (PBMC) from each tube and transfer to a new 50mL tube.

Take extra to care to avoid pulling cells from the ficoll layer (underneath the mononuclear cell layer) as this typically contains a lot of granulocytes. Pulling from the plasma layer is not an issue.

- 4.4.8. Top PBMC with cold DPBS Solution-EDTA to 40mL (ensure at least 2-3 volumes are added) and centrifuge the cell suspension(s) for 10 minutes at 400 x g, 4°C.

- 4.4.9. (Platelet Spin) Discard the supernatant, top tube to 40mL with cold DPBS Solution-EDTA, and centrifuge the cell suspension for 10 minutes at 120 x g, 4°C.

- 4.4.10. Remove the supernatant (**caution: pellet may be loose**), and re-suspend the cell pellet in 4mL **Dulbecco's Phosphate Buffered Saline (DPBS).**

4.5. Cell Counting of COVID Samples

- 4.5.1. Add 0.050mL of sample, 0.050mL of DPBS, and 0.005mL of Solution 13 to a 1.5mL centrifuge tube, incubate for 2 minutes at room temperature.

- 4.5.2. Add 0.100mL of BD Cytofix Fixation Buffer to the samples and incubate 30 minutes, 4°C, and protect from light.

- 4.5.3. Aliquot 0.010mL of sample to the well of a NC-Slide A8 and count on the NC-3000.

Record number and viability below, calculate total cells:

cell number: _____ cells/mL, _____ % viable

4.6. Division of Sample for Analysis and Freeze-down

- 4.6.1. Aliquot up to 1×10^7 cells to a 5mL Falcon Round-Bottom tube and place on ice for subsequent sample clean-up (Step 4.9).

- 4.6.2. Freeze down up to 1×10^8 cells in approximately 1×10^7 aliquots (1mL each) using Cryostor CS10 Medium, a Mr. Frosty, and a -80°C freezer. Record the number of vials frozen: _____ and the cells per cryovial frozen: _____.

4.7. Sample Clean Up for scRNAseq – CD66b and CD235ab removal

- 4.7.1. Centrifuge the single cell suspension for 5 minutes at $400 \times g$, 4°C .
- 4.7.2. Discard the supernatant and resuspend the cell pellet in 50L of DPBS Solution-EDTA; add 10 μL of Human TruStain FcX to the single cell suspension and incubate for 10 minutes, 4°C .
- 4.7.3. Add 10 μL of biotinylated anti-CD66b and biotinylated anti-CD235ab to the sample and incubate for 30 minutes, 4°C .
- 4.7.4. While the single cell suspension is incubating add 0.200mL of BioMag Plus Streptavidin Beads to a 5mL Falcon Round-Bottom tube.
- 4.7.5. Add 2mL of DPBS Solution-EDTA to the BioMag Plus Streptavidin Beads and place on a magnet for 5 minutes.
- 4.7.6. Remove all the supernatant from the BioMag Plus Streptavidin Beads, remove from the magnet and resuspend the beads in 0.100mL of DPBS Solution-EDTA.
- 4.7.7. Once step 4.7.3 is complete, add 3mL of DPBS Solution-EDTA to the single cell suspension and centrifuge for 5 minutes at $400 \times g$, 4°C .
- 4.7.8. Resuspend the single cell suspension in the BioMag Plus Streptavidin Beads from step 4.7.6, and incubate at room temperature for 5 minutes.
- 4.7.9. Add 3mL of **DPBS** to the tube and place on a magnet for 5 minutes.
- 4.7.10. Remove supernatant from tube and transfer to a separate 5mL Falcon Round Bottom tube.

4.8. Sample Clean Up for scRNAseq – Dead Cell Removal

- 4.8.1. Centrifuge the single cell suspension for 5 minutes at $400 \times g$, 4°C , and discard supernatant.
- 4.8.2. Resuspend cell pellet in 0.100mL of Dead Cell Removal Microbeads, mix well, incubate, room temperature, 15 minutes.
- 4.8.3. While the cell suspension is incubating, place an MS Column onto the MACS Multistand and rinse with 0.500mL 1x Binding Buffer Solution.
- 4.8.4. Post incubation, apply cell suspension to the MS Column and capture the flow through in a 5mL Falcon Round Bottom tube. Rinse with 1.5mL of 1x Binding Buffer and capture in the same tube.
- 4.8.5. Centrifuge the single cell suspension for 5 minutes at $400 \times g$, 4°C , and discard supernatant.
- 4.8.6. Resuspend cell pellet in 0.500mL DPBS, and count cells.

4.9. Cell Counting of COVID Samples (10x)

- 4.9.1. Add 0.050mL of sample, 0.050mL of DPBS, and 0.005mL of Solution 13 to a 1.5mL centrifuge tube, incubate for 2 minutes at room temperature.
- 4.9.2. Add 0.100mL of BD Cytofix Fixation Buffer to the samples and incubate 30 minutes, room temperature, and protect from light.

4.9.3. Aliquot 0.010mL of sample to the well of a NC-Slide A8 and count on the NC-3000.

4.9.4. Record number and viability below, calculate total cells:

cell number: _____ cells/mL, _____ % viable

4.10. **10X Encapsulation**

4.10.1. Follow the appropriate 10X protocol (Chromium Next GEM Single Cell 3' Reagent Kits v3.1 User Guide – Rev D) for encapsulation of cells from the airway sample.