

Genomic DNA from environmental biofilm

Protocol at a glance

From Genomic DNA from soil - NucleoSpin® Soil

		2mL of biofilm 18,000g, 30m, +4°C Remove supernatant containing ethanol Suspend the pellet with 700µL of SL1 Transfert the mixed biofilm pellet / SL1 to a NucleoSpin® Bead Tube Type A												
1	Prepare sample													
2	Adjust lysis conditions	150 µL Enhancer SX												
3	Sample lysis	Horizontally vortex 5 min at RT or use other homogenizers according to manufacturers protocol												
4	Precipitate contaminants	  11,000 x g, 2 min 150 µL SL3 Vortex 5 s 0–4 °C, 5 min 11,000 x g, 1 min												
5	Filter lysate	  Load supernatant on NucleoSpin® Inhibitor Removal Column (red ring) 11,000 x g, 1 min												
6	Adjust binding conditions	 250 µL SB Vortex 5 s												
7	Bind DNA	  Load 550 µL sample on NucleoSpin® Soil Column (green ring) 11,000 x g, 1 min Load remaining sample 11,000 x g, 1 min												
8	Wash silica membrane	  <table border="1"> <tbody> <tr> <td>1st</td> <td>500 µL SB</td> <td>11,000 x g, 30 s</td> </tr> <tr> <td>2nd</td> <td>550 µL SW1</td> <td>11,000 x g, 30 s</td> </tr> <tr> <td>3rd</td> <td>650 µL SW2 </td> <td>11,000 x g, 30 s</td> </tr> <tr> <td>4th</td> <td>650 µL SW2 </td> <td>11,000 x g, 30 s</td> </tr> </tbody> </table>	1 st	500 µL SB	11,000 x g, 30 s	2 nd	550 µL SW1	11,000 x g, 30 s	3 rd	650 µL SW2 	11,000 x g, 30 s	4 th	650 µL SW2 	11,000 x g, 30 s
1 st	500 µL SB	11,000 x g, 30 s												
2 nd	550 µL SW1	11,000 x g, 30 s												
3 rd	650 µL SW2 	11,000 x g, 30 s												
4 th	650 µL SW2 	11,000 x g, 30 s												
9	Dry silica membrane	  11,000 x g, 2 min												
10	Elute DNA	  Final elution in 30µL of SE heated at 50°C RT, 1 min 11,000 x g, 30 s												