

HIFI-Barcode SOP –

Filling reference gaps via assembling DNA barcodes using high-throughput sequencing

Version 1 - 201710

Shanlin Liu^{1, 2, 4, #} Chentao Yang^{2, #}, Chengran Zhou^{2, 5 #}, Xin Zhou^{1, 3, *}

1. Beijing Advanced Innovation Center for Food Nutrition and Human Health,
College of Plant Protection, China Agricultural University, Beijing 100193,
People's Republic of China

2. BGI-Shenzhen, Shenzhen, 518083, China

3. National Engineering Research Center for Fruit and Vegetable Processing,
China Agricultural University, Beijing 100193, People's Republic of China

4. Centre for GeoGenetics, Natural History Museum of Denmark, University of
Copenhagen, Copenhagen, Denmark

5. Key Laboratory of Bio-Resources and Eco-Environment, Ministry of
Education, College of Life Sciences, Sichuan University, Chengdu 610000,
China.

[#]Equal contribution

^{*}Corresponding to: Xin Zhou

Table of Contents

1 DESCRIPTION	1
2 DNA PREPARATION.....	3
3 DNA AMPLIFICATION	4
3.1 PRIMER SETS	4
3.2 PCR.....	7
3.2.1 <i>PCR reaction mixture</i>	7
3.2.2 <i>PCR program</i>	7
3.3 PCR GELS.....	8
4 SEQUENCING	9
4.1 HISEQ OR BGISEQ500	9
4.2 PACBIO (OPTIONAL)	9
5 BIOINFOMATICS DATA ANALYSIS.....	10
5.1 HIFI-BARCODE-HISEQ ASSEMBLY	10
5.2 HIFI-BARCODE-PACBIO ASSEMBLY	12
6 REFERENCES	12
7 CITATION	22

1 Description

Over the past decade, biodiversity scientists have dedicated tremendous efforts in constructing DNA reference barcodes for rapid species registration and identification. Although analytical cost for standard DNA barcoding has been significantly reduced since early 2,000, further dramatic reduction on barcoding costs is unlikely because the Sanger sequencing is approaching its limits in throughput and chemistry cost. Constraints in barcoding cost not only led to unbalanced barcoding efforts around the globe, but also refrained High-Throughput-Sequencing (HTS) based taxonomic identification from applying binomial species names, which provide crucial linkages to biological knowledge. We developed an Illumina-based pipeline, HIFI-Barcode, to produce full-length COI barcodes from pooled PCR amplicons generated by individual specimens. The new pipeline generated accurate barcode sequences that were comparable to Sanger standards, even for different haplotypes of the same species that were only a few nucleotides different from each other. Additionally, the new pipeline was much more sensitive in recovering amplicons at low quantity. The HIFI-Barcode pipeline successfully recovered barcodes from over 78% of the PCR reactions that didn't show clear bands on the electrophoresis gel. Moreover, sequencing results based on the single molecular sequencing platform, Pacbio, confirmed the accuracy of the HIFI-Barcode results. Altogether, the new pipeline can provide an improved solution to produce full-length reference barcodes at about 1/10 of the current cost, enabling construction of comprehensive barcode libraries for local fauna, leading to a feasible direction for DNA barcoding global biomes.




Figure 1. Schematic illustration of HIFI-Barcode pipeline.

2 DNA preparation

Individual Genomic DNA could be extracted using the Glass Fiber Plate method following manufacturer's protocol or other existing method.

3 DNA amplification

3.1 Primer sets

Ninety-six pairs of different tags were added to both ends of a common COI barcode primer set (LCO1490 and HCO2198) with each tag containing 5 bps allowing for ≥ 2 bp differences from each other.

Table 1. Indexed Primer sequences

Primer	5'to3'	Primer	5'to3'
Rev001	AAAGCTAACTTCAGGGTGACCAAAAAATCA	For001	AAAGCGGTCAACAAATCATAAAGATATTGG
Rev002	AACAGTAACTTCAGGGTGACCAAAAAATCA	For002	AACAGGGTCAACAAATCATAAAGATATTGG
Rev003	AACCTTAACTTCAGGGTGACCAAAAAATCA	For003	AACCTGGTCAACAAATCATAAAGATATTGG
Rev004	AACTCTAACTTCAGGGTGACCAAAAAATCA	For004	AACTCGGTCAACAAATCATAAAGATATTGG
Rev005	AAGCATAAACTTCAGGGTGACCAAAAAATCA	For005	AAGCAGGTCAACAAATCATAAAGATATTGG
Rev006	AAGGTTAACTTCAGGGTGACCAAAAAATCA	For006	AAGGTGGTCAACAAATCATAAAGATATTGG
Rev007	AAGTGTAACCTTCAGGGTGACCAAAAAATCA	For007	AAGTGGGTCAACAAATCATAAAGATATTGG
Rev008	AATGGTAACTTCAGGGTGACCAAAAAATCA	For008	AATGGGGTCAACAAATCATAAAGATATTGG
Rev009	AACTTAACTTCAGGGTGACCAAAAAATCA	For009	AACTGGTCAACAAATCATAAAGATATTGG
Rev010	ACAGATAAACTTCAGGGTGACCAAAAAATCA	For010	ACAGAGGTCAACAAATCATAAAGATATTGG
Rev011	ACCATTAACTTCAGGGTGACCAAAAAATCA	For011	ACCATGGTCAACAAATCATAAAGATATTGG
Rev012	ACCTATAAACTTCAGGGTGACCAAAAAATCA	For012	ACCTAGGTCAACAAATCATAAAGATATTGG
Rev013	ACGAATAAACTTCAGGGTGACCAAAAAATCA	For013	ACGAAGGTCAACAAATCATAAAGATATTGG
Rev014	ACGTTTAACTTCAGGGTGACCAAAAAATCA	For014	ACGTTGGTCAACAAATCATAAAGATATTGG
Rev015	ACTGTAACTTCAGGGTGACCAAAAAATCA	For015	ACTGTGGTCAACAAATCATAAAGATATTGG
Rev016	AGAACTAACTTCAGGGTGACCAAAAAATCA	For016	AGAACGGTCAACAAATCATAAAGATATTGG
Rev017	AGACATAAACTTCAGGGTGACCAAAAAATCA	For017	AGACAGGTCAACAAATCATAAAGATATTGG
Rev018	AGAGTTAACTTCAGGGTGACCAAAAAATCA	For018	AGAGTGGTCAACAAATCATAAAGATATTGG
Rev019	AGATGTAACTTCAGGGTGACCAAAAAATCA	For019	AGATGGGTCAACAAATCATAAAGATATTGG

Rev020	AGCTTTAACTTCAGGGTGACCAAAAAATCA	For020	AGCTTGGTCAACAAATCATAAAGATATTGG
Rev021	AGGATTAACTTCAGGGTGACCAAAAAATCA	For021	AGGATGGTCAACAAATCATAAAGATATTGG
Rev022	AGTAGTAACTTCAGGGTGACCAAAAAATCA	For022	AGTAGGGTCAACAAATCATAAAGATATTGG
Rev023	AGTCTTAACTTCAGGGTGACCAAAAAATCA	For023	AGTCTGGTCAACAAATCATAAAGATATTGG
Rev024	AGTGATAAACTTCAGGGTGACCAAAAAATCA	For024	AGTGAGGTCAACAAATCATAAAGATATTGG
Rev025	AGTTCTAACTTCAGGGTGACCAAAAAATCA	For025	AGTTCGGTCAACAAATCATAAAGATATTGG
Rev026	ATACCTAACTTCAGGGTGACCAAAAAATCA	For026	ATACCGGTCAACAAATCATAAAGATATTGG
Rev027	ATCACTAACTTCAGGGTGACCAAAAAATCA	For027	ATCACGGTCAACAAATCATAAAGATATTGG
Rev028	CAAAGTAACTTCAGGGTGACCAAAAAATCA	For028	CAAAGGGTCAACAAATCATAAAGATATTGG
Rev029	CAACTTAACTTCAGGGTGACCAAAAAATCA	For029	CAACTGGTCAACAAATCATAAAGATATTGG
Rev030	CAATCTAACTTCAGGGTGACCAAAAAATCA	For030	CAATCGGTCAACAAATCATAAAGATATTGG
Rev031	CAGAATAAACTTCAGGGTGACCAAAAAATCA	For031	CAGAAGGTCAACAAATCATAAAGATATTGG
Rev032	CATACTAACTTCAGGGTGACCAAAAAATCA	For032	CATACGGTCAACAAATCATAAAGATATTGG
Rev033	CATCATAAACTTCAGGGTGACCAAAAAATCA	For033	CATCAGGTCAACAAATCATAAAGATATTGG
Rev034	CCAATTAACTTCAGGGTGACCAAAAAATCA	For034	CCAATGGTCAACAAATCATAAAGATATTGG
Rev035	CGATTTAACTTCAGGGTGACCAAAAAATCA	For035	CGATTGGTCAACAAATCATAAAGATATTGG
Rev036	CGTATTAACCTTCAGGGTGACCAAAAAATCA	For036	CGTATGGTCAACAAATCATAAAGATATTGG
Rev037	CGTTATAAACTTCAGGGTGACCAAAAAATCA	For037	CGTTAGGTCAACAAATCATAAAGATATTGG
Rev038	CTAACTAACTTCAGGGTGACCAAAAAATCA	For038	CTAACGGTCAACAAATCATAAAGATATTGG
Rev039	CTACATAAACTTCAGGGTGACCAAAAAATCA	For039	CTACAGGTCAACAAATCATAAAGATATTGG
Rev040	CTATGTAACTTCAGGGTGACCAAAAAATCA	For040	CTATGGGTCAACAAATCATAAAGATATTGG
Rev041	CTCAATAAACTTCAGGGTGACCAAAAAATCA	For041	CTCAAGGTCAACAAATCATAAAGATATTGG
Rev042	CTGATTAACTTCAGGGTGACCAAAAAATCA	For042	CTGATGGTCAACAAATCATAAAGATATTGG
Rev043	CTGTATAAACTTCAGGGTGACCAAAAAATCA	For043	CTGTAGGTCAACAAATCATAAAGATATTGG
Rev044	CTTAGTAACTTCAGGGTGACCAAAAAATCA	For044	CTTAGGGTCAACAAATCATAAAGATATTGG
Rev045	CTTCTTAACTTCAGGGTGACCAAAAAATCA	For045	CTTCTGGTCAACAAATCATAAAGATATTGG
Rev046	GAAACTAACTTCAGGGTGACCAAAAAATCA	For046	GAAACGGTCAACAAATCATAAAGATATTGG
Rev047	GAACATAAACTTCAGGGTGACCAAAAAATCA	For047	GAACAGGTCAACAAATCATAAAGATATTGG
Rev048	GAATGTAACTTCAGGGTGACCAAAAAATCA	For048	GAATGGGTCAACAAATCATAAAGATATTGG
Rev049	GACTTTAACTTCAGGGTGACCAAAAAATCA	For049	GACTTGGTCAACAAATCATAAAGATATTGG

Rev050	GAGATTAACTTCAGGGTGACCAAAAAATCA	For050	GAGATGGTCAACAAATCATAAAGATATTGG
Rev051	GAGTATAAACTTCAGGGTGACCAAAAAATCA	For051	GAGTAGGTCAACAAATCATAAAGATATTGG
Rev052	GATAGTAACTTCAGGGTGACCAAAAAATCA	For052	GATAGGGTCAACAAATCATAAAGATATTGG
Rev053	GATCTTAACTTCAGGGTGACCAAAAAATCA	For053	GATCTGGTCAACAAATCATAAAGATATTGG
Rev054	GATGATAAACTTCAGGGTGACCAAAAAATCA	For054	GATGAGGTCAACAAATCATAAAGATATTGG
Rev055	GATTCTAACTTCAGGGTGACCAAAAAATCA	For055	GATTGGTCAACAAATCATAAAGATATTGG
Rev056	GCAAATAAACTTCAGGGTGACCAAAAAATCA	For056	GCAAAGGTCAACAAATCATAAAGATATTGG
Rev057	GCTATTAACTTCAGGGTGACCAAAAAATCA	For057	GCTATGGTCAACAAATCATAAAGATATTGG
Rev058	GCTTATAAACTTCAGGGTGACCAAAAAATCA	For058	GCTTAGGTCAACAAATCATAAAGATATTGG
Rev059	GGAATTAACCTTCAGGGTGACCAAAAAATCA	For059	GGAATGGTCAACAAATCATAAAGATATTGG
Rev060	GGATATAAACTTCAGGGTGACCAAAAAATCA	For060	GGATAGGTCAACAAATCATAAAGATATTGG
Rev061	GGTTTTAACTTCAGGGTGACCAAAAAATCA	For061	GGTTTGGTCAACAAATCATAAAGATATTGG
Rev062	GTAGATAAACTTCAGGGTGACCAAAAAATCA	For062	GTAGAGGTCAACAAATCATAAAGATATTGG
Rev063	GTCATTAACTTCAGGGTGACCAAAAAATCA	For063	GTCATGGTCAACAAATCATAAAGATATTGG
Rev064	GTGAATAAACTTCAGGGTGACCAAAAAATCA	For064	GTGAAGGTCAACAAATCATAAAGATATTGG
Rev065	GTGTTTAACTTCAGGGTGACCAAAAAATCA	For065	GTGTTGGTCAACAAATCATAAAGATATTGG
Rev066	GTTACTAACTTCAGGGTGACCAAAAAATCA	For066	GTTACGGTCAACAAATCATAAAGATATTGG
Rev067	G TTCATAAACTTCAGGGTGACCAAAAAATCA	For067	G TTCAGGTCAACAAATCATAAAGATATTGG
Rev068	TAAGGTAACTTCAGGGTGACCAAAAAATCA	For068	TAAGGGGTCAACAAATCATAAAGATATTGG
Rev069	TACTGTAACTTCAGGGTGACCAAAAAATCA	For069	TACTGGGTCAACAAATCATAAAGATATTGG
Rev070	TAGGATAAACTTCAGGGTGACCAAAAAATCA	For070	TAGGAGGTCAACAAATCATAAAGATATTGG
Rev071	TAGTCTAACTTCAGGGTGACCAAAAAATCA	For071	TAGTCGGTCAACAAATCATAAAGATATTGG
Rev072	TATCGTAACTTCAGGGTGACCAAAAAATCA	For072	TATCGGGTCAACAAATCATAAAGATATTGG
Rev073	TATGCTAACTTCAGGGTGACCAAAAAATCA	For073	TATGCGGTCAACAAATCATAAAGATATTGG
Rev074	TCACATAAACTTCAGGGTGACCAAAAAATCA	For074	TCACAGGTCAACAAATCATAAAGATATTGG
Rev075	TCAGTTAACTTCAGGGTGACCAAAAAATCA	For075	TCAGTGGTCAACAAATCATAAAGATATTGG
Rev076	TCATGTAACTTCAGGGTGACCAAAAAATCA	For076	TCATGGGTCAACAAATCATAAAGATATTGG
Rev077	TCCAATAAACTTCAGGGTGACCAAAAAATCA	For077	TCCAAGGTCAACAAATCATAAAGATATTGG
Rev078	TCCTTTAACTTCAGGGTGACCAAAAAATCA	For078	TCCTTGGTCAACAAATCATAAAGATATTGG
Rev079	TCGATTAACTTCAGGGTGACCAAAAAATCA	For079	TCGATGGTCAACAAATCATAAAGATATTGG

Rev080	TCGTATAAACTTCAGGGTGACCAAAAAATCA	For080	TCGTAGGTCAACAAATCATAAAGATATTGG
Rev081	TCTCTTAACTTCAGGGTGACCAAAAAATCA	For081	TCTCTGGTCAACAAATCATAAAGATATTGG
Rev082	TGAAGTAACTTCAGGGTGACCAAAAAATCA	For082	TGAAGGGTCAACAAATCATAAAGATATTGG
Rev083	TGACTTAACTTCAGGGTGACCAAAAAATCA	For083	TGACTGGTCAACAAATCATAAAGATATTGG
Rev084	TGAGATAAACTTCAGGGTGACCAAAAAATCA	For084	TGAGAGGTCAACAAATCATAAAGATATTGG
Rev085	TGCTATAAACTTCAGGGTGACCAAAAAATCA	For085	TGCTAGGTCAACAAATCATAAAGATATTGG
Rev086	TGGAATAAACTTCAGGGTGACCAAAAAATCA	For086	TGGAAGGTCAACAAATCATAAAGATATTGG
Rev087	TGTACTAACTTCAGGGTGACCAAAAAATCA	For087	TGTACGGTCAACAAATCATAAAGATATTGG
Rev088	TGTCATAAACTTCAGGGTGACCAAAAAATCA	For088	TGTCAGGTCAACAAATCATAAAGATATTGG
Rev089	TGTGTAACTTCAGGGTGACCAAAAAATCA	For089	TGTGTGGTCAACAAATCATAAAGATATTGG
Rev090	TTACGTAACTTCAGGGTGACCAAAAAATCA	For090	TTACGGGTCAACAAATCATAAAGATATTGG
Rev091	TTAGCTAACTTCAGGGTGACCAAAAAATCA	For091	TTAGCGGTCAACAAATCATAAAGATATTGG
Rev092	TTCTCTAACTTCAGGGTGACCAAAAAATCA	For092	TTCTCGGTCAACAAATCATAAAGATATTGG
Rev093	TTGACTAACTTCAGGGTGACCAAAAAATCA	For093	TTGACGGTCAACAAATCATAAAGATATTGG
Rev094	TTGCATAAACTTCAGGGTGACCAAAAAATCA	For094	TTGCAGGTCAACAAATCATAAAGATATTGG
Rev095	TTGGTTAACTTCAGGGTGACCAAAAAATCA	For095	TTGGTGGTCAACAAATCATAAAGATATTGG
Rev096	TTTCCTAACTTCAGGGTGACCAAAAAATCA	For096	TTTCCGGTCAACAAATCATAAAGATATTGG

3.2 PCR

3.2.1 PCR reaction mixture

Each PCR reaction should contain 1 µl of DNA template, 16.2 µl of molecular biology grade water, 3 µl of 10X reaction buffer (Mg²⁺ plus), 2.5 µl of dNTPs mix (10 mM), 1 µl of forward and reverse primers (10 mM), and 0.3 µl of TaKaRa Ex Taq polymerase (5 U/µl).

3.2.2 PCR program

The amplification program includes a thermocycling profile of 94°C for 1 min, 5 cycles of 94°C for 30 sec, 45°C for 40 sec, and an extension at 72°C for 1 min,

followed by 35 cycles of 94°C for 30 sec, 51°C for 40 sec, and 72°C for 1 min, with a final extension at 72°C for 10 min, and finally holding at 12°C.

3.3 PCR gels

All amplicons could be visualized on a 1.2% 96 Agarose E-gel (Biowest Agarose).

4 SEQUENCING

4.1 Hiseq or BGISEQ500

- All PCR products from each plate should be pooled using 1 µl per sample resulting in two 96 µl mixtures.
- PCR amplicons should be fragmented to construct library of an insert-size of 250 bp and sequenced with a strategy of 150 PE.

4.2 Pacbio (optional)

- All PCR products from each plate should be pooled using 1 µl per sample resulting in two 96 µl mixtures.
- Sequencing using PacBio RS II.

5 BIOINFOMATICS DATA ANALYSIS

5.1 HIFI-Barcode-Hiseq assembly

HIFI-Barcode-Hiseq could be downloaded from <https://github.com/comery/HIFI-barcode-hiseq>.

DESCRIPTION

HIFIBarcode is used to produce full-length COI barcodes from pooled PCR amplicons generated by individual specimens.

INSTALLATION

- Clone from github

```
$ git clone https://github.com/comery/HIFI-barcode-hiseq.git
```

```
$ tar -zxvf HIFIBarcode.v1.0.tar.gz
```

Source code

(1) Wrapper script: HIFIBarcode.v1.0.pl: the wrapper script to run other Perl scripts to do the work you choose.

(2) Main scripts: Seven Perl scripts in folder HIFIBarcode.V1.0/bin/, including:

- 1_split_extract.pl
- 2_uniqu_sort_cluster.Pro.pl
- 3_sep_extract_overlap.pl
- 4_cluster_fromend.pl
- 5_forgap_filling.pl
- 6_rename_kmer.pl
- 7_final.pl

(3) Published softwares:

VSEARCH v2.4.4. Rognes, Torbjørn, et al. "VSEARCH: a versatile open source tool for metagenomics." *PeerJ* 4 (2016): e258. COPE CMR v1.0.3. Liu, Binghang, et al. "COPE: an accurate k-mer-based pair-end reads connection tool to facilitate genome assembly." *Bioinformatics* 28.22 (2012): 2870-2874. SOAPBarcode. Liu, Shanlin, et al. "SOAPBarcode: revealing arthropod biodiversity through assembly of Illumina shotgun sequences of PCR amplicons." *Methods in Ecology and Evolution* 4.12 (2013): 1142-1150.

Pre-requisites

PERL v5

EXAMPLES

- 1: run wrapper script to get shell text and then sh runHIFIBarcode.sh to run HIFIBarcode

```
perl HIFIBarcode.v1.0.pl --fq1 test_1.fq --fq2 test_2.fq --index  
index_primer.txt --length 5 --cpunum 10 --outdir test --outpre testout
```

- 2: run shell

```
sh test/runHIFIBarcode.sh
```

NOTE

The proposed method assemble short-read Illumina sequences based on k-mer sequence matches and such misassembly was not observed in their real data, but it's still possible. The pipeline also provides an additional note file with a suffix of "note.txt" with notes alerting users about the possibilities which have similar or same scores comparing to their best alternative.

LATEST RELEASE

Version 1.1 201709

5.2 HIFI-Barcode-Pacbio assembly

HIFI-Barcode-Pacbio could be downloaded from <https://github.com/comery/HIFI-barcode-pacbio>.

DESCRIPTION

HIFIBarcode is used to produce full-length COI barcodes from pooled PCR amplicons generated by individual specimens.

INSTALLATION

- Clone from github

```
$ git clone https://github.com/comery/HIFI-barcode-pacbio.git
```

- Or go to website <https://github.com/comery/HIFI-barcode-pacbio> and click 'Download ZIP'

Requirements

(1) software

- PicBio smrtanalysis download from <http://www.pacb.com/products-and-services/analytical-software/smrt-analysis/>

(2) programming language

- standard perl
- standard python(python2 is ok)

(3) perl module

- Bio::Perl(exactly Bio::Seq)

(4) main perl and python scripts in bin/

- 1.primer_like_extract.pl
- 2.cluster_count_passes_length.pl

- `ccs_passes.py` | see more information to go to [https://github.com/PacificBiosciences/Bioinformatics-Training/wiki/Extracting-Reads-of-Insert-\(CCS\)-number-of-passes](https://github.com/PacificBiosciences/Bioinformatics-Training/wiki/Extracting-Reads-of-Insert-(CCS)-number-of-passes) or directly download from https://raw.githubusercontent.com/PacificBiosciences/Bioinformatics-Training/master/scripts/ccs_passes.py

- `fish_ccs.pl`

DATA requirements:

(1) pacbio original H5 file input

- `01.data/*.h5`

(2) primers list

- `experiment_data/primer.lst`

for GGTCAACAAATCATAAAGATATTGG

rev TAAACTTCAGGGTGACCAAAAAATCA

(3) index(barcodes for identifying samples) list

- `experiment_data/index.xls`

001 AAAGC

002 AACAG

003 AACCT

004 AACTC

005 AAGCA

... ..

(4) samples_location.tab

- samples name and corresponding location in 96-cell plate

1 A01

2 B01

3 C01

4 D01

5 E01

. ...

Overview of steps

If you installed PacBio smrtanalysis, I suppose you get the file path of setup.sh, more about Pacbio Data : <http://www.pacb.com/wp-content/uploads/SMRT-Link-User-Guide-v4.0.0.pdf>

e.g: setup_path='/path/PicBio/smrtanalysis/current/etc/setup.sh'

step 1 extract CCS from h5 files

Input:

my_inputs.fofn (*my_inputs.fofn contains files list of Pacbio H5 file in 01.data/, like this-> ./01.data/m170506_092957_42199_c101149142550000001823255607191735_s1_p0.1.bax.h5*)

Output:

log

data *.ccs.fasta *.ccs.fastq *.ccs.h5 reads_of_insert.fasta
reads_of_insert.fastq slots.pickle

workflow

results

Run:

```
$ source /path/PicBio/smrtanalysis/current/etc/setup.sh
```

```
$ fofnToSmrtpipeInput.py my_inputs.fofn > my_inputs.xml
```

```
$ smrtpipe.py --params=settings.xml xml:input.xml
```

step 2 extract passes number from CCS h5 files

Input:

/data/*.ccs.h5

Output:

ccs_passes.lst

Run:

*Note: Before run the scripts, please in sure that you have run source
/path/PicBio/smrtanalysis/current/etc/setup.sh*

```
$ python bin/ccs_passes.py data/*.ccs.h5 >ccs_passes.lst
```

step 3 filtering CCS by passes number (>15)**Input:**

ccs_passes.lst

data/reads_of_insert.fasta

Output:

ccs_passes_15.fa

Run:

```
$ awk '$2>=15{print $1}' ccs_passes.lst >ccs_passes_15.lst
```

```
$ perl ./bin/fish_ccs.pl ccs_passes_15.lst  
data/reads_of_insert.fasta >ccs_passes_15.fa
```

step 4 assigning CCS to samples by index**Input:**

experiment_data/primer.fa

experiment_data/index.xls

ccs_passes_15.fa

Output:

Note: "outdir" name is up to you, here default value is "02.assignment"

02.assignment/

assign.log.txt

ccs.successfully_assigned.fa

check.ccs_passes_15.fa.log

Run:

```
$ perl ./bin/1.primer_like_extract.pl -p experiment_data/primer.fa -index  
experiment_data/index.xls -fa ccs_passes_15.fa -cm 2 -cg 1
```

step 5 clustering CCS of each sample to find best one

Input:

ccs.successfully_assigned.fa

check.ccs_passes_15.fa.log

ccs_passes.lst

Output:

cluster.top1.fas

cluster.id.txt

cluster.all.fa

Run:

```
$ cd 02.assignment/  
  
$ perl ../bin/2.cluster_count_passes_length.pl -ccs  
ccs.successfully_assigned.fa -pattern check.ccs_passes_15.fa.log -  
passes ../ccs_passes.lst
```

```
$ perl ../bin/change_name-location.pl cluster.top1.fas >hifi-barcode-  
pacbio.cluster.top1.fa
```

ALL DONE!

So, "hifi-barcode-pacbio.cluster.top1.fa" is final result!

LATEST RELEASE

Version 1.0 201707

6 REFERENCES

1. Hebert PD, Hollingsworth PM, Hajibabaei M: **From writing to reading the encyclopedia of life.** *Phil Trans R Soc B* 2016, **371**(1702):20150321.
2. Yu DW, Ji YQ, Emerson BC, Wang XY, Ye CX, Yang CY, Ding ZL: **Biodiversity Soup: Metabarcoding of arthropods for rapid biodiversity assessment and biomonitoring.** *Methods in Ecology and Evolution* 2012, **3**(4):613-623.
3. Zhou X, Li Y, Liu S, Yang Q, Su X, Zhou L, Tang M, Fu R, Li J, Huang Q: **Ultra-deep sequencing enables high-fidelity recovery of biodiversity for bulk arthropod samples without PCR amplification.** *GigaScience* 2013, **2**(1):4.
4. Taberlet P, Coissac E, Hajibabaei M, Rieseberg LH: **Environmental DNA.** *Molecular Ecology* 2012, **21**:1789–1793.
5. Calvignac-Spencer S, Merkel K, Kutzner N, Kühl H, Boesch C, Kappeler PM, Metzger S, Schubert G, Leendertz FH: **Carion fly-derived DNA as a tool for comprehensive and cost-effective assessment of mammalian biodiversity.** *Molecular Ecology* 2013, **22**(4):915-924.
6. Schnell IB, Thomsen PF, Wilkinson N, Rasmussen M, Jensen LR, Willerslev E, Bertelsen MF, Gilbert MTP: **Screening mammal biodiversity using DNA from leeches.** *Current biology* 2012, **22**(8):R262-R263.
7. Ratnasingham S, Hebert PDN: **BOLD: The Barcode of Life Data System (<http://www.barcodinglife.org>).** *Molecular ecology notes* 2007, **7**(3):355-364.
8. Zhou X, Frandsen PB, Holzenthal RW, Beet CR, Bennett KR, Blahnik RJ, Bonada N, Cartwright D, Chuluunbat S, Cocks GV: **The Trichoptera barcode initiative: a strategy for generating a species-level Tree of Life.** *Phil Trans R Soc B* 2016, **371**(1702):20160025.
9. Miller SE, Hausmann A, Hallwachs W, Janzen DH: **Advancing taxonomy and bioinventories with DNA barcodes.** *Phil Trans R Soc B* 2016, **371**(1702):20150339.
10. Zhou X, Robinson JL, Geraci CJ, Parker CR, Flint Jr OS, Etnier DA, Ruiter D, DeWalt RE, Jacobus LM, Hebert PD: **Accelerated construction of a regional DNA-barcode reference library: caddisflies (Trichoptera) in the Great Smoky Mountains National Park.** *Journal of the North American Benthological Society* 2011, **30**(1):131-162.
11. Zhou X, Adamowicz SJ, Jacobus LM, DeWalt RE, Hebert PD: **Towards a**

comprehensive barcode library for arctic life-Ephemeroptera, Plecoptera, and Trichoptera of Churchill, Manitoba, Canada. *Frontiers in zoology* 2009, **6**(1):1.

12. Telfer AC, Young MR, Quinn J, Perez K, Sobel CN, Sones JE, Levesque-Beaudin V, Derbyshire R, Fernandez-Triana J, Rougerie R: **Biodiversity inventories in high gear: DNA barcoding facilitates a rapid biotic survey of a temperate nature reserve.** *Biodiversity data journal* 2015(3).

13. Valentini A, Taberlet P, Miaud C, Civade R, Herder J, Thomsen PF, Bellemain E, Besnard A, Coissac E, Boyer F: **Next-generation monitoring of aquatic biodiversity using environmental DNA metabarcoding.** *Molecular Ecology* 2016.

14. Hajibabaei M, Shokralla S, Zhou X, Singer GAC, Baird DJ: **Environmental barcoding: a next-generation sequencing approach for biomonitoring applications using river benthos.** *PloS one* 2011, **6**(4):e17497.

15. Coissac E, Riaz T, Puillandre N: **Bioinformatic challenges for DNA metabarcoding of plants and animals.** *Molecular Ecology* 2012, **21**:1834–1847.

16. Bienert F, Dedanieli S, Miquel C, Coissac E, Poillot C, Brun JJ, Taberlet P: **Tracking earthworm communities from soil DNA.** *Molecular Ecology* 2012, **21**(8):2017-2030.

17. Baird DJ, Hajibabaei M: **Biomonitoring 2.0: a new paradigm in ecosystem assessment made possible by next-generation DNA sequencing..** *Molecular Ecology* 2012, **21**(8):2039-2044.

18. Gillett CP, Crampton-Platt A, Timmermans MJ, Jordal B, Emerson BC, Vogler AP: **Bulk de novo mitogenome assembly from pooled total DNA elucidates the phylogeny of weevils (Coleoptera: Curculionoidea).** *Molecular biology and evolution* 2014:msu154.

19. Krulwich R, LANDER E: **Cracking the Code of Life.** In.: *Public Broadcasting Service*; 2001.

20. Shokralla S, Gibson JF, Nikbakht H, Janzen DH, Hallwachs W, Hajibabaei M: Next-generation DNA barcoding: **using next-generation sequencing to enhance and accelerate DNA barcode capture from single specimens.** *Molecular Ecology Resources* 2014, **14**(5):892-901.

21. Shokralla S, Porter TM, Gibson JF, Dobosz R, Janzen DH, Hallwachs W, Golding GB, Hajibabaei M: **Massively parallel multiplex DNA sequencing for specimen identification using an Illumina MiSeq platform.** *Scientific reports* 2015, **5**.

22. Meier R, Wong W, Srivathsan A, Foo M: **\$1 DNA barcodes for reconstructing**

complex phenomes and finding rare species in specimen-rich samples. *Cladistics* 2016, 32(1):100-110.

23. Liu S, Li Y, Lu J, Su X, Tang M, Zhang R, Zhou L, Zhou C, Yang Q, Ji Y: **SOAPBarcode: revealing arthropod biodiversity through assembly of Illumina shotgun sequences of PCR amplicons.** *Methods in Ecology and Evolution* 2013, 4(12):1142-1150.

24. Ivanova NV, Dewaard JR, Hebert PD: **An inexpensive, automation-friendly protocol for recovering high-quality DNA.** *Molecular Ecology Resources* 2006, 6(4):998-1002.

25. Vrijenhoek R: **DNA primers for amplification of mitochondrial cytochrome c oxidase subunit I from diverse metazoan invertebrates.** *Molecular marine biology and biotechnology* 1994, 3(5):294-299.

26. Rognes T, Flouri T, Nichols B, Quince C, Mahé F: **VSEARCH: a versatile open source tool for metagenomics.** *PeerJ* 2016, 4:e2584.

27. Liu B, Yuan J, Yiu SM, Li Z, Xie Y, Chen Y, Shi Y, Zhang H, Li Y, Lam TW: **COPE: an accurate k-mer-based pair-end reads connection tool to facilitate genome assembly.** *Bioinformatics* 2012, 28(22):2870-2874.

28. Letunic I, Bork P: **Interactive tree of life (iTOL) v3: an online tool for the display and annotation of phylogenetic and other trees.** *Nucleic Acids Research* 2016:gkw290.

29. Li H, Durbin R: **Fast and accurate short read alignment with Burrows–Wheeler transform.** *Bioinformatics* 2009, 25(14):1754-1760.

30. Eid J, Fehr A, Gray J, Luong K, Lyle J, Otto G, Peluso P, Rank D, Baybayan P, Bettman B: **Real-time DNA sequencing from single polymerase molecules.** *Science* 2009, 323(5910):133-138.

31. Pimm SL, Jenkins CN, Abell R, Brooks TM, Gittleman JL, Joppa LN, Raven PH, Roberts CM, Sexton JO: **The biodiversity of species and their rates of extinction, distribution, and protection.** *Science* 2014, 344(6187):1246752.

32. Díaz S, Demissew S, Carabias J, Joly C, Lonsdale M, Ash N, Larigauderie A, Adhikari JR, Arico S, Báldi A: **The IPBES Conceptual Framework—connecting nature and people.** *Current Opinion in Environmental Sustainability* 2015, 14:1-16.

33. Luo A, Lan H, Ling C, Zhang A, Shi L, Ho SY, Zhu C: **A simulation study of sample size for DNA barcoding.** *Ecology and Evolution* 2015, 5(24):5869-5879.

34. Tang M, Tan M, Meng G, Yang S, Su X, Liu S, Song W, Li Y, Wu Q, Zhang A: **Multiplex sequencing of pooled mitochondrial genomes—a crucial step toward biodiversity**

analysis using mito-metagenomics. *Nucleic Acids Research* 2014, **42**(22):e166-e166.

35. Tan M, Zhang R, Xiang C, Zhou X: **The complete mitochondrial genome of the pumpkin fruit fly, *Bactrocera tau* (Diptera: Tephritidae).** *Mitochondrial DNA Part A* 2016, **27**(4):2502-2503.

36. Timmermans MJ, Dodsworth S, Culverwell C, Bocak L, Ahrens D, Littlewood DT, Pons J, Vogler AP: **Why barcode? High-throughput multiplex sequencing of mitochondrial genomes for molecular systematics.** *Nucleic acids research* 2010, **38**(21):e197-e197.

7 CITATION

Liu, Yang, Zhou, Zhou. (2017). Filling reference gaps via assembling DNA barcodes using high-throughput sequencing - moving toward to barcoding the world.