

TABLE I. Dissociation Medium

Stock solution	Component	Stock solution concentration (g/100 ml)	Final concentration in medium (mM)
1	$\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$	1.47	5
	$\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$	0.49	1
	KCL	4.18	2.8
	TES ^a	4.95	11
2	Na_2HPO_4	0.19	0.67
	KH_2PO_4	0.12	0.44
3	Na Pyruvate	1.10	5
	Na_3 Citrate $\cdot 2\text{H}_2\text{O}$	2.94	5

^aTES = N-tris (hydroxymethyl) methyl-2-aminomethanesulfonic acid (Sigma). HEPES can be substituted for TES. Make up stock solutions (which are 20 $\times$ concentrated in distilled water), filter sterilize (0.22 μ or 0.45 μ Millipore; sterilization by autoclaving causes precipitation), and store frozen. In making up stock solution 1 it is important that the ingredients be added and dissolved one at a time to prevent precipitation. To make the dissociation medium, add 5 ml of each stock solution to 85 ml distilled water, and then adjust the pH to 6.9–7.0 with 2 N NaOH (usually 0.1 ml). Filter sterilized, the medium can be stored at 4°C.