

iPSC editing with TALENs
Protocol by Ruilin Tian, Jason Hong and Martin Kampmann
9/2019

Reference:

Tian et al (2019). CRISPR Interference-Based Platform for Multimodal Genetic Screens in Human iPSC-Derived Neurons. Neuron pii: S0896-6273(19)30640-3.
[Epub ahead of print] PubMed PMID: 31422865.

Reagents

Component	Cat. No	Manufacturer
Lipofectamine Stem	STEM00008	Invitrogen (Life Tech/Fisher)
E8 Medium	A1517001	Gibco (Life Tech/Fisher)
Stemflex Medium	A3349401	Gibco (Life Tech/Fisher)
Matrigel	356231	Corning (Discovery Labware)
Knockout DMEM	10829-018	Gibco (Life Tech/Fisher)
DPBS	14190144	D8537-24X500mL (Sigma)
Accutase	A11105-01	Gibco (Life Tech/Fisher)
Opti-Mem Reduced Serum	31985070	Gibco (Life Tech/Fisher)
Rock inhibitor (Y-27632)	1254/10	R&D Systems (Tocris/Biotechne)
TALEN-L (e.g. for CLYBL: pZT-C13-L1)	#62196	Addgene
TALEN-R (e.g. for CLYBL: pZT-C13-R1)	#62197	Addgene
BCL-XL plasmid	pEF1-BCL-XL-wpre-polyA P1102	Gift from Xiaobing Zhang, described in PMID: 30239926

1. Pre-coat 6 well plate
 - a. Pre-coat 6 well plate with Matrigel (diluted 1:50 with Knockout DMEM) by adding 1mL to each well and let it sit in incubator for at least 20 minutes.
2. Passaging and Plating iPSC Cells
 - a. Grow iPSCs until they are ~85% confluent in one well of a 6 well plate or larger format.

- b. Remove old medium and wash them with DPBS 1x and add Accutase for 3-5 min to lift the cells.
- c. Singularize the cells by gently pipetting them up and down several times.
- d. In a 15mL conical tube, add DPBS and then add the lifted cells in Accutase.
- e. Spin down at 200g for 5 min.
- f. Aspirate supernatant and resuspend cells in Stemflex.
- g. Perform cell count.
- h. Remove the Matrigel in the pre-coated plate and add appropriate amount of Stemflex medium with Rock inhibitor (1000x)
- i. Re-seed 0.5M cells into a 6 well plate so that it can reach ~60% confluency the next day.
- j. Put plate in incubator and culture overnight.

3. Transfection

- a. The following day, change the media with new E8 media.
- b. Prepare the following mixes:
 - i. Tube 1
 1. Opti-Mem - 100uL
 2. Lipofectamine Stem Reagent - 10uL
 - ii. Tube 2
 1. Opti-Mem - 100uL
 2. DNA mix (0.5-5.0ug total)
 - a. Optimized TALENS ratio
 - i. Your Plasmid of choice - 1.5ug
 - ii. TALENS L - 0.75ug
 - iii. TALENS R - 0.75ug
 - iv. BCL-XL - 0.3ug
- c. Add tube 2 to tube 1 and mix well.
- d. Incubate the mixture for 10 minutes.
- e. Add the whole mixture to your cells in the 6 well plate.
- f. Incubate overnight and change the media to Stemflex the next day.
- g. If confluent, passage and expand them into 10cm plate.
- h. When confluent in 10cm plate, continue to selection protocol.