

Protocol for nuclei isolation from mouse lung for single nucleus RNASeq

Version 1.0 (3/13/20)

Jeffrey R Koenitzer and Benjamin D Humphreys

This protocol is adapted directly from Joshi et al. <https://www.protocols.io/view/single-nucleus-isolation-from-frozen-human-lung-ti-zu8f6zw>, with adjustments to RNase inhibitor concentrations and removal of FACS sorting steps.

Materials	Company	Cat. No.	Storage
Reagent for buffer and stock			
1x DPBS	Gibco	14190144	4 °C
Nuclei EZ Lysis Buffer (from kit)	Sigma	N 3408 (from NUC-101)	4 °C
cComplete ULTRA Tablets, Mini, EDTA-free, EASYpack	Roche	05 892 791 001	4 °C
RNasin Plus Ribonuclease Inhibitors	Promega	N2615	-20°C
SUPERaseIN RNase Inhibitor	Thermo Fisher Scientific	AM2696	-20°C
RNase free H ₂ O	Thermo Scientific	AM9938	RT
RNaseZap	Ambion	AM9780	RT
Bovine serum albumin, 10% solution, nuclease free	Millipore Sigma	1187358001	-20°C
Tools			
GentleMacs C tube	Miltenyi Biotech	130-093-237	
GentleMacs Dissociator	Miltenyi Biotech	130-093-235	
pluriStrainer 40 µm	pluriSelect	43-50040	
pluriStrainer 5 µm	pluriSelect	43-50005	
Fuchs-Rosenthal disposable hemocytometer	INCYTO	DHC-F015	
Falcon 15 mL Polystyrene Conical Tube	Fisher Scientific	352095	
RNase-free 50 ml Conical Tubes	Ambion	AM12502	
TPP 60mm Tissue Culture Dishes	MIDSCI	TP93060	

Buffers and Solutions

Working Solutions

cOmplete stock (10x)

- 1 mL Nuclei EZ Prep lysis buffer
- 1 tablet of cOmplete ULTRA tablets

Lysis Buffer - 2 ml per $< 6 \text{ mm}^3$ tissue

- 200 μl cOmplete stock
- 1.775 ml of Nuclei EZ Prep lysis buffer
- 12.5 μl of RNasin Plus
- 12.5 μl of SUPERaseIN

Wash Buffer - 4 ml per $< 6 \text{ mm}^3$ tissue

- 3.575 ml of dPBS
- 400 μl of 10% BSA
- 12.5 μl of RNasin Plus
- 12.5 μl of SUPERaseIN


Resuspension Buffer – 1 ml

- 0.987 ml of dPBS
- 10 μl of 10% BSA
- 6.25 μl of RNasin Plus
- 6.25 μl of SUPERaseIN

Procedure

1. Precool all instruments (including centrifuges), buffers, and tubes.
Note: All steps are performed on ice or in cold room (4°C) to minimize RNA degradation.
2. Remove mouse lung sample from -80°C freezer and trim a ~6 mm³ piece. Thaw on small plastic weighing boat until able to insert 28G needle to tissue, then inject 1 ml of Lysis Buffer to 'inflate' the tissue. Add the remaining Lysis Buffer and chop to the smallest pieces possible with scissors (60 s).
3. Transfer the minced tissue and buffer to a GentleMacs C tube. Close, invert, and transfer directly to MACS Tissue dissociator.
4. Run m_lung_01 program and m_lung_02 program in sequence, stopping the latter after 20 seconds. Place tube on ice.
5. Reduce foam by centrifugation (1 min at 750 x g). Use a wide-bore pipette tip to pass lysate through 40 µm strainer to 50 mL conical tube. Wash strainer with 4 mL Wash Buffer.
6. Pass suspension through a 5 µm strainer into 50 ml conical tube.
7. Centrifuge at 600 x g for 5 min at 4°C.
8. Resuspend in 500 µl Resuspension Buffer. Count nuclei by hemocytometer and dilute to desired concentration (e.g. 10,000 per µl). *Note: it is easier to resuspend in lower volumes and dilute than to concentrate nuclei via further centrifugation.*
9. Proceed to 10x Chromium.

Expected Results



40x

Nuclear suspension with scant debris



Foam after GentleMacs (L) and after centrifugation (R)