

QUALITATIVE AND QUANTITATIVE COMPOSITION

One aerosol aluminium spray can contains micronized zeolite 97.5%, ionic silver 2.5%, n-butane (propellant)

DEVICE DESCRIPTION

ARGOGEN is a Medical Device consisting in an absorbent mineral spray powder containing 2.5% ionic silver. The mineral powder is white and composed essentially of a main ingredient micronized zeolite, a silicate mineral with a porous structure and high absorption capacity. Ionic silver, a medicinal substance, is chemically-incorporated into the natural porous structure of zeolite and have an antimicrobial function. The device, when sprayed on wounds and skin lesions, absorbs wound exudates, creating a dry and clean environment that facilitates the natural healing process. Ionic silver, at the same time, provides an effective barrier to microbial penetration sources preventing further contamination and keeps the area of application protected from microbial proliferation.

PACKAGING

The device is packaged in aerosol aluminium spray can. N-butane is used as propellant.

The final content can be 10 g powder/125 mL or 4 g powder/50 mL.

INTENDED USE

ARGOGEN is intended to treat skin lesions and wounds of different origin, when an absorption of exudates is needed, such as burn wounds, pressure ulcers, abrasions. Product remains on the treated area, does not rub off, and is easily removed with the use of saline spray.

INDICATIONS

This Medical Device is indicated for adult patients in case of:

- pressure ulcers (bedsores);
- first and second degree burn wounds;
- skin lesions, such as wounds, abrasions, skin cuts, minor burns.

CONTRAINDICATIONS

This Medical Device is contraindicated for:

- People with known allergy to silver;
- Purposes other than its intended use.

It should not be used:

- In patients with severe heart, lung and other organ disease;
- In patients with not-treated diabetes mellitus or diabetic foot syndrome;
- In patients with peripheral artery disease;
- Thrombocytopenic patients;
- In case of concomitant infection of another site or osteomyelitis or erysipelas or phlegmon;
- In case of complicated deep tissue infection requiring systemic antimicrobial therapy;
- In case of active viral hepatitis (A/B/C) or active HIV infection or active syphilis;
- In presence of neoplastic growth in the ulcer;
- In case of history or clinical signs of impairment of the cochlea or vestibular system;
- In patients with neuromuscular diseases (ie Myasthenia gravis, Parkinson's disease);
- In case of aminoglycoside treatment or other antibiotic treatments for chronic venous ulcers administered less than four weeks before.

In case of Paediatric population or Pregnant or breastfeeding women, please consult your Doctor before use or use the product under Medical Supervision.

INSTRUCTION FOR USE

Wounds should be cleaned before the application with sterile isotonic solution or sterile water. Strongly shake aerosol can before use and remove the cap. Be sure to direct the dispenser erogator in the direction of the lesion to be treated, keeping the can in a vertical position at a distance of about 10 cm from the surface to be treated. Spray in short intervals, lightly pressing the dispenser button and pointing the spray on the area to be treated. Do not disperse in the air, do not inhale, and avoid any possible flammable risks during the use. In case of accidental ingestion or inhalation of the dust, consult your doctor. The application must continue without interruption once or twice a day until complete healing, and, when appropriate, ensure the product on the wound using a sterile gauze or a sterile breathable bandage. At each renewal of the application, thoroughly cleanse the lesion with with sterile isotonic solution or sterile water. The powder contained in will tend to stick to the wound (due to the presence of exudate) or to the macerated skin, while on the surrounding healthy and dry skin it will tend to remove spontaneously. It is recommended not to use the product for more than 4 weeks. Consult your doctor in case of incomplete healing or if unwanted symptoms arise. The powder dispensed exposed to light may tend to become dark or brown, due the oxidation of the silver ions present in the formulation. However this phenomenon does not compromise product effectiveness.

SPECIAL WARNINGS AND PRECAUTIONS FOR USE

EXTERNAL USE ONLY. Extremely flammable aerosol. Pressurized container: may burst if heated. Keep away from heat, sparks, open flames, hot surfaces. No smoking. Do not pierce or burn, even after use. Protect from sunlight. Do no expose to temperatures exceeding 50°C/122°F. Do

not spray on an open flame or other ignition source. Keep out of reach of children.

INTERACTION

Typically applied medicinal products such as proteolytic enzymes medicinal products at the wound site must be removed and the wound cleansed prior application as per instruction for use previously specified. Remaining antibacterial medicinal products reduce the activity of the product by decreasing its efficacy.

UNDESIRABLE COLLATERAL EFFECTS

None known and not available in the published literature. The product is biocompatible. There are no reported side effects. In case of overdose, the absorption of massive doses of silver can cause argyria.

SERIOUS INCIDENT REPORTING:

The user and/or patient are recommended to report any serious incident that has occurred in relation to the device to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

SHELF-LIFE

36 months referred to the product in intact packaging correctly stored.

TRANSPORTATION AND STORAGE

Extremely flammable aerosol. Keep away from heat, sparks, open flames, hot surfaces. Protect from sunlight. Do no expose to temperatures exceeding 50°C/122°F.

The expiration date refers to the product in its intact packaging and correctly stored. Do not use the Medical Device after the expiration date.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP): the document is available on the European portal Eudamed at the following URL: <https://ec.europa.eu/tools/eudamed> or alternatively in SIGNORINI MEDICALE S.r.l. website (<https://signorinimedicae.it>)

MANUFACTURED by
SIGNORINI MEDICALE S.r.l.
via dell'artigianato, 14
Albaredo d'Adige (VR)



DISTRIBUTED by

Instructions for Use Edition: 02/2022

SYMBOLS USED ON THE PRODUCT LABELS					
	Medical Device according to MDR 2017/745		Manufacturer		Batch/Lot Number
	Expiry Date		Keep away from light and heat sources		Keep Dry
	Consult instructions for use		Conformity to Directive 75/324/CEE (aerosol dispenser)		Nominal content (volume)
	Nominal total capacity of the container (according to Directive 2007/45/EC)		Flammable (according to CLP – EU Regulation n° 1272/2008)		