

Prepare buffers (Volumes above are indicated per sample)

10x stocks of Wash buffer base (no spermidine/Tween) and Labeling buffer base (no spermidine/SAM/Hia5)

2 mL Wash buffer -> take **488 uL** and add Digitonin + 50XProtInh to make **Lysis buffer**

50 uL **labeling buffer** / 250.000 cells (typically ~ 1E6 cells -> **200 uL**)

Make 25 mM **spermidine stocks fresh monthly** and store at -20deg.

Add 2uL spermidine (liquid, ~6.38M) to 350 uL 0.1M HCl and 160 uL H₂O – **check pH with strips!**

Protocol 6mA labeling and HMW DNA extraction

- 1) Place the Dounce homogenizer and pestles on ice, chill the centrifuge to 4 °C and preheat the ThermoMixer to 37°C..
- 2) **Carefully** transfer 3-4 (2mm diameter) tissue punch biopsies (~25mg) to the chilled Dounce homogenizer. Keep the Dounce homogenizer on ice during the entire disruption process
- 3) Add 500 uL of the Fiber-Seq Lysis buffer and let the tissue thaw for 1 minute.
- 4) Gently homogenize the tissue 10X with pestle A and 10X with B.
 - Push the tissue with the pestle firmly into the bottom of the Dounce chamber with each stroke (Down + Up = 1X).
 - Keep the tissue between tip of pestle and the bottom of the Dounce chamber for thorough homogenization.
 - Homogenate may become foamy, but this is not a cause for concern. In the next step, transfer any foam that forms.
- 5) Incubate on ice for 5min before adding 1000uL of Wash buffer.
- 6) Transfer the lysate to a 2mL Protein LoBind microcentrifuge tube.
- 7) Rinse the pestles and homogenizer with the remaining 500 uL Wash buffer and add to the sample.
- 8) Pellet homogenate by centrifuging at 700 x g and 4 °C for 5 min. Discard supernatant.
- 9) Resuspend the pellet in 200uL of Hia5 labeling buffer – use a 1mL or wide bore tip.
- 10) Incubate on a ThermoMixer at 37deg and 900 rpm for 30 min.
- 11a) **Continue with Circulomics CBB Tissue protocol from step 8 onwards**
- 11b) **Continue with NEB Monarch HMW**

Protocol Nanopore sequencing (LSK-110, PromethION)

According to ONT protocol *Genomic DNA by Ligation (SQK-LSK110)* with the following modifications:

- 1) Start with 3-4ug of HMW DNA in 150 uL and shear 25x with a 26G needle or in Megaruptor to 35kb.
- 2) Adjust volumes end-prep and FFPE repair step accordingly (*i.e.* vol x 3), omit control strand (CS)
- 3) Extend end-prep and FFPE repair steps from 5 to 30min (*i.e.* 30min at 20deg and 30min at 65deg)
- 4) Extend ligation step to 1h at RT
- 5) Elute the AMPure cleanups for 10 minutes and 20 minutes after the end-prep and ligation steps.
- 6) This should yield a ~3x library, aim to load near the high end of the 5-50fmol range, typically ~8uL

Note 1: library prep yield is typically 30-50%

Note 2: the amount loaded can be reduced during subsequent flushes to balance seq yield with # flushes

Fiber-seq brain tissue (BA24/MB) ASAP

10X Wash/Lysis base (10X-WLB) (200uL / sample)

	Stock	Final for 10X	V (uL)
Tris-HCl pH 7.4	1000 mM	100 mM	100
NaCl	5000 mM	100 mM	20
Water			880
<i>TotVol</i>			1000

10X Hia5 labeling base (10X-H5B) (20uL / sample)

	Stock	Final for 10X	V (uL)
Tris-HCl pH 8	1000 mM	150 mM	150
NaCl	5000 mM	150 mM	30
KCl	1000 mM	600 mM	600
EDTA pH 8	500 mM	10 mM	20
EGTA pH 8	250 mM	5 mM	20
Water			180
<i>TotVol</i>			1000

1X Wash buffer

	Stock	Final	V (uL)	X5
10X WLB	10 X	1 X	230	1150
BSA	10 %	0.1 %	23	115
Spermidine pH 7.4	25 mM	0.5 mM	46	230
Tween-20	10 %	0.1 %	23	115
Water			1978	9890
<i>TotVol</i>			2300	11500

1X Lysis buffer

	Stock	Final	V (uL)	X5
Wash buffer			488	2440
Digitonin	5 %	0.02 %	2	10
ProteaseInh	50 X	1 X	10	50
<i>TotVol</i>			500	2500

1X Hia5 labeling buffer

	Stock	Final	V (uL)	X5
10X-H5B	10 X	1 X	20	100
BSA	10	0.1	2	10
Spermidine pH 7.4	25	0.5	4	20
SAM	32	0.8	5	25
Hia5 enzyme	250	5	4	20
Water			165	825
			200	1000