

WP3 CAPACITY BUILDING WORKSHOP

University of Nairobi
25th-26th October 2022

Detection of mycotoxigenic species in plant samples with the portable qPCR Q3 device

Marco Cereda

STMicroelectronics, Italy



life.augmented

STMicroelectronics (ST): company presentation

October 2022

We are creators and makers of technology



One of the world's largest semiconductor companies



48,000 employees of which
8,400 in R&D



Over **80** sales & marketing
offices serving over **200,000**
customers across the globe



14 manufacturing sites



Signatory of the United Nations Global Compact (UNGC)
Member of the Responsible Business Alliance (RBA)

Our strategy stems from key long-term enablers

Smart Mobility



Helping car manufacturers make driving safer, greener, and more connected for everyone

Power & Energy



Enabling industries to increase energy efficiency everywhere and the use of renewable energy

Internet of Things & Connectivity



Supporting the proliferation of smart, connected IoT devices with products, solutions & ecosystems

We offer quality, flexibility, and supply security



Sweden
Norrköping

France
Crolles
Rousset
Tours
Rennes

Italy
Agrate
Catania
Marcianise



Morocco
Bouskoura

Malta
Kirkop

- Front-End (Wafer fabrication)
- Back-End (Assembly & Test)

China
Shenzhen

Philippines
Calamba



Malaysia
Muar

Singapore

We are drivers of your innovation

Advanced R&D centers around the world for close collaboration with operations, customers, and partners



~**8,400** people working in R&D and product design

13.5% of revenues invested in R&D in 2021


187 active R&D partnerships

~**18,500** patents
~**550** new filings in 2021

Open innovation with startups in **15** Proof of Concept centers
50 startups engaged in our programs

- T3.2: “Developing and adopting simple, rapid and inexpensive new ICT-based tools for detection and quantification of plant pests and pathogens”, led by ST
- ST has developed a **molecular analysis** workflow to analyze crops of interest, exploiting its “**Q3**” portable device

- Target: **DNA/RNA**
- Plant pests and pathogens can be detected by molecular analyses
- Most used technique: real-time Polymerase Chain Reaction (**qPCR**)
- **Pros: sensitivity; specificity; accuracy**
- **Cons: quite complex; expensive and bulky instrument → tackled in T3.2**
- Workflow:



Q3 (BIO2BIT®) platform: overview

- **Compact, Point-of-Care qPCR platform by ST:**
 1. Single-use Lab-on-a-Chip cartridge
 2. Portable instrument (14 x 7 x 8.5 cm; 300 g)
 3. Software
- Advantages vs standard qPCR equipment: **portable, less expensive, easier to use, faster**



- Background: general-purpose qPCR platform
- Objectives in EWA-BELT:
 - A. Task 3.2A: development and optimization** of one target application: molecular detection of pathogens in plants of interest for West/East Africa
 - Easy DNA extraction procedure
 - qPCR on Q3 (targets, reagents, protocol)
 - EWA-BELT-dedicated Q3 software, to easily run the qPCR test, whose output files can be uploaded to PLANTHEAD portal
 - B. Task 3.2B: training to end-users and field validation in Africa**

- Development and optimization of different protocols for DNA extraction from **maize flour** (samples provided by UNISS)
- Study of DNA sequences for qPCR detection of **fungal species** of interest, that **can produce mycotoxins**:
 - *Aspergillus flavus*, *A. ochraceus*, *A. parasiticus*, *Fusarium culmorum*, *F. graminearum*, *F. proliferatum*, *F. verticillioides*, *Penicillium verrucosum*
 - **But some strains** exist that are **genetically unable to produce mycotoxins** (e.g. *A. flavus* strains in Aflasafe®)...
- Setup and optimization of qPCR reaction mixtures
 - Triplex (3 targets): fungus 1, fungus 2, maize DNA sequence as control
- Porting of qPCR onto Q3



- Optimization of protocols for DNA extraction from **flour of different crops** (maize, peanuts, wheat, sorghum, coffee samples by UNISS)
- Study of **fungal** DNA sequences of **mycotoxin biosynthesis genes**:
 - **Ochratoxin A** (OTA)
 - **Aflatoxins** (Afla)
 - **Fumonisin** (Fum)
 - **Trichothecenes** (Tri)
 - **Nivalenol** trichothecene (Niv)
 - 3-acetyldeoxynivalenol trichothecene (**3ADON**)
 - 15-acetyldeoxynivalenol trichothecene (**15ADON**)
 - **Zearalenone** (Zea)
- Setup, optimization, and porting to Q3 of qPCR reaction mixtures:
 - Triplex: mycotoxin 1, mycotoxin 2, pan-plant DNA sequence as control







If present in sample, there is a contamination by some **fungus** (which exact species is not important) that for sure is **genetically able to produce** that **toxin(s)**!

Porting to Q3 on phase 2 (1/3)


- **qPCR** reaction mixtures ported on **Q3** (time-to-result = **45 min**)
- Development and testing of **EWA-BELT-dedicated Q3 software**:
 - **Very easy-to-use**: step-by-step writings to guide the user; all parameters pre-defined
 - Embedded algorithm for raw data interpretation, showing **results easily understandable** by non-expert users
 - **Output files uploadable to PLANTHEAD portal**



Porting to Q3 on phase 2 (2/3)



Porting to Q3 on phase 2 (3/3)




DNA extraction demo

- Allows to **purify DNA** from all the other components of a sample
- **Needed before qPCR**, because DNA must be released from cells to be accessible to qPCR reagents; also, qPCR enzyme is inhibited by many natural substances present in samples
- **Plants are a tough matrix to purify DNA from**: solid & complex matrix, cells have a resistant cell wall, many qPCR inhibitors (polyphenols, humic acids, polysaccharides)
- Different methods: purely chemical (laborious, toxic reagents), spin column (needs centrifuge), **magnetic beads** (typically needs robot)

Chosen extraction method: magnetic beads

- Aim: simplest DNA protocol, with very low equipment needed
- Material needed for developed protocol: only two pipettes and one small magnetic rack
- Core technology: paramagnetic, positively-charged microbeads with pores → high surface where DNA can bind
 - Diameter: 0.5 to 8 μm



*Scanning electron micrograph
of magnetic beads*

DNA extraction: steps

1. **Lysis** → to break cell walls and cell/nucleus membranes, releasing DNA (and all other cellular components) in solution
2. **Binding** (to beads) → appropriate salt/pH conditions make DNA bind to beads
3. **Washes** → to remove all other components, leaving DNA bound to beads
4. **Elution** → to detach DNA from beads, releasing it into solution



Developed DNA extraction protocol

0. Sample aliquoting and dissolution:

- a. Weigh **200 mg of sample** (flour) into a tube.
- b. Add **300 µL water**.

1. Lysis:

- a. Add **250 µL CTAB**. Incubation 5'.
- b. Add **250 µL TLA**. Incubation 5'.
- c. Transfer 400 µL supernatant (lysed sample) to a new tube for further processing. Discard the original tube.

2. DNA binding to beads:

- a. Add **30 µL NucleoMag C-Beads** + **400 µL MC2** (Binding Buffer). Magnet 30". Discard liquid. Remove tube from magnet.

3. Washes:

- a. Add **400 µL MC3** (Wash Buffer 1). Magnet 30". Discard liquid. Remove tube from magnet.
- b. Add **400 µL MC4** (Wash Buffer 2). Magnet 30". Discard liquid. Remove tube from magnet.
- c. Add **400 µL Ethanol 80%**. Magnet 30". Discard liquid. Keep tube close to magnet.
- d. Gently add 400 µL MC5 (Wash Buffer 3). Do not pipette up/down. Incubation 1' (not more). Discard liquid. Remove tube from magnet.

4. Elution:

- a. Add **100 µL MC6** (Elution Buffer). Incubation 5'. Magnet 1'. Transfer liquid, containing purified DNA, to a new tube.
- b. Store purified DNA at +4°C for up to 24 hours, or at -20°C for longer storage.

qPCR demo

Q3 cartridge layout

**OTA-
Afla**


**Fum-
Tri**

**3ADON-
15ADON**

**Niv-
Zea**

- Polymerase chain reaction (PCR) is a lab technique for rapidly producing (amplifying) millions-billions of copies of a specific DNA segment, which can then be detected/studied in greater detail
- Short synthetic DNA fragments called primers are designed to select a **target DNA segment** to be **amplified**
- With multiple rounds of DNA synthesis, an enzyme (DNA polymerase) exponentially amplifies that segment
 - $C_n = C_0 * 2^n$, where n is the number of rounds
- Every round consists in **thermal cycle** between $\approx 95^\circ\text{C}$ (denaturation) and $\approx 60^\circ\text{C}$ (primers annealing to target + extension of primers (amplification) by enzyme)
- RNA can also be targeted, adding a second enzyme



PCR: mechanism



- In traditional PCR, reaction is analyzed at the end point, when a plateau has typically been reached → No information on initial quantity of target DNA/RNA sequence in the sample, only qualitative answer (yes/no)
- Real-time PCR, or quantitative PCR (**qPCR**) combines PCR with a way to measure amplification while it is occurring
- Measurement during exponential phase provides **quantitative information**. Typically, optical (fluorescence) detection
- Technological requirements for a qPCR platform:
 1. Precise and accurate **thermal control** (alternate heating/cooling)
 2. **Fluorescence measurement**, typically at different wavelengths

Quantification in qPCR

- A fluorescence threshold is defined in the exponential phase; the moment when fluorescence overcomes this threshold is called **C_q** (or C_t) value for a certain sample → main output of qPCR
- C_q is inversely proportional to the Log of target DNA/RNA initial quantity



Q3 (BIO2BIT®) platform: overview

- **Compact, Point-of-Care qPCR platform by ST:**
 1. Single-use Lab-on-a-Chip cartridge
 2. Portable instrument (14 x 7 x 8.5 cm; 300 g)
 3. Software
- Advantages vs standard qPCR equipment: **portable, less expensive, easier to use, faster**





Q3 platform elements: cartridge

- Six reaction chambers, 5 μ L each
- Key element is a **silicon chip**, manufactured by ST Italy, integrating thermo-resistive heater and temperature sensor for precise and accurate **thermal control**
- Silicon chip assembled with plastics, wax and RFID tag to obtain cartridge
- Special **wax** avoids evaporation of reagents during thermal cycling, and seals reagents at the end
- **RFID tag**, writable/readable by Q3 instrument and all NFC-enabled devices, allows robust traceability and anti-fraud protection



Cartridge silicon chip: metal side

- Chips are produced in parallel onto an 8" silicon wafer



Distributed metal resistor
acting as **heater**


Concentrated metal
resistor acting as
temperature sensor

Contacts to Q3 instrument


- Silicon is a good thermal conductor

Cartridge silicon chip: reagents side


- Hydrophilic regions: silicon oxide
- Hydrophobic regions: polymer patterned through photolithography




Well cross sections



Solid wax + liquid qPCR mix before analysis start



Thermal cycling
 $\geq 60^{\circ}\text{C}$




Melted wax + liquid qPCR mix after analysis start


Chip temperature distribution

Simulation of temperature distribution @ 85°C

Characterization with thermo-sensitive liquid crystal @ 60°C



$$\Delta T = 0.1^{\circ}\text{C}$$



$$\Delta T = 0.13^{\circ}\text{C}$$

Q3 platform elements: instrument



1. Drives the cartridge allowing the thermal cycling required by qPCR
 2. Acquires **fluorescence signals** coming from cartridge during reaction, on 3 different optical bands
 3. Allows communication with host through USB or Bluetooth
- Only 14 x 7 x 8.5 cm and weighs barely 300 g
 - Internal memory for backup storage of up to 350 tests
 - RFID tag reading/writing and cryptographic scheme for data protection and chip cloning prevention
 - Motorised lid and touch button with multi-color status LED






Instrument: optics

Compact fluorescence microscope,
with no moving parts:

- 4 LEDs illuminating the cartridge
- 4 ST CMOS image sensors, closely spaced, looking at the same area
- Each sensor and LED is behind a different optical filter



Images obtained
from each optical
channel:



Q3 platform elements: software


1. Background: **general-purpose software**

- Full protocol flexibility
- General data analysis (Cq calculation)





2. Developed for **EWA-BELT**: “closed”, dedicated **software**

- Pre-defined protocol
- Writings guide the user step-by-step
- Dedicated data analysis (positive/negative for each mycotoxin)




Results of phase 1 (1/2)

- **Simple DNA extraction protocol** developed and tested on maize: only two pipettes and one magnet needed; **≈ 40 min**
- **Two qPCR reaction mixtures optimized**, targeting *F. culmorum*, *F. graminearum*, *A. ochraceus*, *P. verrucosum*
- Results on gold standard qPCR equipment:



Results of phase 1 (2/2)

- qPCR reaction mixtures ported on **Q3** (time-to-result = **45 min**)




- Fungi clearly detected
- Good consistency of replicates

Preliminary results of phase 2 (1/3)


- Simple DNA extraction protocols developed and tested on maize and peanut: only two pipettes and one magnet needed; **≈ 40 min**
- **Four qPCR mixes developed**, targeting genes for 8 mycotoxins
- Results on gold standard qPCR equipment:

A. flavus-contaminated Kenyan peanut



OTA-Afla mix

F. proliferatum-contaminated maize




Fum-Tri mix


Preliminary results of phase 2 (2/3)

ADON mix

F. culmorum-contaminated maize




F. graminearum-contaminated maize



F. graminearum-contaminated maize


Niv-Zea mix



Preliminary results of phase 2 (3/3)

- **qPCR** reaction mixtures ported on **Q3** (time-to-result = **45 min**)
- Tests run with **EWA-BELT-dedicated Q3 software**

F. proliferatum-contaminated maize



- Optimization of DNA extraction protocols for different plants (wheat, sorghum, coffee)
- Further optimization of qPCR detection of some mycotoxins
- Fine-tuning of algorithm embedded in EWA-BELT software for raw data interpretation

Possible applications of Q3

Infectious Diseases

COVID-19 ⁽¹⁾,
Sepsis ⁽¹⁾,
Malaria ⁽²⁾,
Ebola ⁽³⁾



Cardiology

Antiplatelet
drug
metabolism ^(4,5)



Pediatrics

Meningitis,
Pertussis ⁽⁶⁾



Gynecology

HPV,
Thrombophilia



Oncology

Pharmaco-
genetics



Agro-food

Parasites in
mollusks ⁽⁷⁾,
meat
identification,
EWA-BELT



(1) <https://alifax.com/assets/Uploads/STM-BRO001-MOLECULAR-MOUSE-03-06-17062022-EN.pdf>

(2) <http://dx.doi.org/10.1016/j.jmoldx.2019.04.008>

(3) <http://dx.doi.org/10.1016/j.jviromet.2017.11.009>

(4) <http://dx.doi.org/10.1016/j.jacc.2018.02.029>

(5) <http://dx.doi.org/10.1016/j.cca.2015.10.003>

(6) <http://dx.doi.org/10.1109/SBMicro.2016.7731331>

(7) <https://doi.org/10.1007/s13197-019-03972-7>

Thank you!
