

DNA isolation from blood using DNeasy Blood & Tissue Kit (QIAGEN 69504)

Blood 100 μ L in 1.5 mL tube	↓
↓	6,000 x g, 1', RT
+ 10 μ L of Lysozyme (200 mg/mL in DW)	↓
↓	column placed in a new 2-ml tube
+ 10 μ L of Triton X-100 (10% (v/v) in DW)	↓
↓	+ 500 μ L of Buffer AW1
37°C, 30 min	↓
↓	6,000 x g, 1', RT
+ 20 μ L Proteinase K (supplied)	↓
↓	column placed in a new 2-mL tube
+ 80 μ L of PBS	↓
↓	+ 500 μ L of Buffer AW2
+ 4 μ L of RNase A (100 mg/mL in DW)	↓
↓	20,000 x g, 3', RT
RT, 2 min	↓
↓	column placed in a new 1.5-mL tube
+ 200 μ L of Buffer AL & Vortex mixing	↓
↓	+ 200 μ L of Buffer AE
56°C, 15 min	↓
↓	RT, 1'
cool down to RT	↓
↓	6,000 x g, 1'
+ 200 μ L of EtOH & Vortex mixing	↓
↓	eluent
apply onto column placed in a new 2-mL tube	