

Catalase activity assessment

Material:

1. Sample
2. Catalase buffer (potassium phosphate buffer + Tween + EDTA)
3. H_2O_2 10mM

Procedure:

1. The H_2O_2 decomposition rate is measured by spectrophotometry set at 240 nm and 25 °C, recording within 15-s intervals over 150 s.
2. Blank the instrument versus a quartz cuvette containing 600 μL of potassium phosphate buffer (50 mM; pH 7.0) with Tween 0.002 % and EDTA 1 %.
3. Pipette 20 μL of H_2O_2 10 mM into a quartz cuvette. Run only one test at a time. Read until 45s.
4. Add the supernatant (3 μL) to H_2O_2 , reading until 150s.
5. The results were normalized by protein content and expressed as a percentage of the placebo group.