

ELISA IL-1 β (DuoSet – R&D Systems, Minneapolis, MN, EUA)

1. Add 100 μ L of capture antibody (Work concentration= 4 μ g/mL) in 96-well plate.

To reach the work concentration, dilute 50 μ L from the capture antibody aliquot at the concentration of 800 μ g/mL (with 40 μ g) in 10mL of Carbonate-Bicarbonate Buffer to obtain a concentration of 4 μ g/mL.

2. Incubate the plate overnight in a wet chamber at 4°C
3. Wash the plate 3x with PBS-T (250 μ L/well)
4. Block with 200 μ L of PBS-T 5% skimmed milk powder
5. Incubate 3h at room temperature
6. Wash the plate 3x with PBS-T
7. To prepare the curve, dilute 25 μ L from the recombinant aliquot (2500pg/mL) in 225 μ L of Reagent Diluent (RD) 1% BSA-PBS to obtain the concentration of 250pg/mL

First point of the curve: 250pg/mL.

Add 100 μ L of IL-1 Rec on the 1st and 2nd well.

Add 100 μ L of RD 1% BSA-PBS on 2nd, 3rd, 4th, 5th, 6th, 7th and 8th well.

Serial dilution from the 2nd well to the last well.

8. Add 100 μ L of samples
9. Incubate for 3h at room temperature
10. Wash the plate 3x with PBS-T
11. Add 100 μ L of detection antibody (Work concentration =200ng/mL)

To reach the work concentration, dilute 50 μ L from the detection antibody aliquot at the concentration of 40 μ g/mL (with 2 μ g) in 10mL of RD to obtain a concentration of 200ng/mL

12. Incubate for 2h at room temperature
13. Wash the plate 3x with PBS-T
14. Add 100 μ L of Streptavidin 1:200
15. Incubate in dark chamber for 30 minutes at room temperature
16. Wash the plate 3x with PBS-T
17. Add 100 μ L of substrate

1 tablet of OPD (SigmaFast OPD, Sigma, Life Science, EUA) + 1 buffer + 20 μ L of distilled water → to do 2 plates

18. Incubate in dark chamber for 20 minutes at room temperature
19. Block with 50 μ L of H₂SO₄ 2N
20. Epoch reading at 492 nm (Test-dependent wavelength)(Epoch-BioTek, Winooski, VT, EUA)

ELISA IL-6 (DuoSet – R&D Systems (Minneapolis, MN, EUA))

1. Add 100µL of capture antibody (Work concentration= 2µg/mL) in 96-well plate.

To reach the work concentration, dilute 50µL from the capture antibody aliquot at the concentration of 800µg/mL (with 40µg) in 10mL of Carbonate-Bicarbonate Buffer to obtain a concentration of 2µg/mL.

2. Incubate the plate overnight in a wet chamber at 4°C
3. Wash the plate 3x with PBS-T (250µL/well)
4. Block with 200µL of PBS-T 5% skimmed milk powder
5. Incubate 3h at room temperature
6. Wash the plate 3x with PBS-T
7. To prepare the curve, dilute 25µL from the recombinant aliquot (6000pg/mL) in 225µL of Reagent Diluent (RD) 1% BSA-PBS to obtain the concentration of 600pg/mL

First point of the curve: 600pg/mL.

Add 100µL of IL-1 Rec on the 1st and 2nd well.

Add 100µL of RD 1% BSA-PBS on 2nd, 3rd, 4th, 5th, 6th, 7th and 8th well.

Serial dilution from the 2nd well to the last well.

8. Add 100µL of samples
9. Incubate for 3h at room temperature
10. Wash the plate 3x with PBS-T
11. Add 100µL of detection antibody (Work concentration =50ng/mL)

To reach the work concentration, dilute 50µL from the detection antibody aliquot at the concentration of 10µg/mL (with 0,5µg) in 10mL of RD to obtain a concentration of 50ng/mL

12. Incubate for 2h at room temperature
13. Wash the plate 3x with PBS-T
14. Add 100µL of Streptavidin 1:200
15. Incubate in dark chamber for 30 minutes at room temperature
16. Wash the plate 3x with PBS-T
17. Add 100µL of substrate

1 tablet of OPD (SigmaFast OPD, Sigma, Life Science, EUA) + 1 buffer + 20µL of distilled water → to do 2 plates

18. Incubate in dark chamber for 20 minutes at room temperature
19. Block with 50µL of H₂SO₄ 2N
20. Epoch reading at 492 nm (Test-dependent wavelength)(Epoch-BioTek, Winooski, VT, EUA)

ELISA IL-8 (DuoSet – R&D Systems, (Minneapolis, MN, EUA))

1. Add 100µL of capture antibody (Work concentration= 4µg/mL) in 96-well plate.

To reach the work concentration, dilute 50µL from the capture antibody aliquot at the concentration of 800µg/mL (with 40µg) in 10mL of Carbonate-Bicarbonate Buffer to obtain a concentration of 4µg/mL.

2. Incubate the plate overnight in a wet chamber at 4°C
3. Wash the plate 3x with PBS-T (250µL/well)
4. Block with 200µL of PBS-T 5% skimmed milk powder
5. Incubate 3h at room temperature
6. Wash the plate 3x with PBS-T
7. To prepare the curve, dilute 25µL from the recombinant aliquot (20000pg/mL) in 225µL of Reagent Diluent (RD) 1% BSA-PBS to obtain the concentration of 2000pg/mL

First point of the curve: 2000pg/mL.

Add 100µL of IL-1 Rec on the 1st and 2nd well.

Add 100µL of RD 1% BSA-PBS on 2nd, 3rd, 4th, 5th, 6th, 7th and 8th well.

Serial dilution from the 2nd well to the last well.

8. Add 100µL of samples
9. Incubate for 3h at room temperature
10. Wash the plate 3x with PBS-T
11. Add 100µL of detection antibody (Work concentration =20ng/mL)

To reach the work concentration, dilute 50µL from the detection antibody aliquot at the concentration of 4µg/mL (with 0,5µg) in 10mL of RD to obtain a concentration of 20ng/mL

12. Incubate for 2h at room temperature
13. Wash the plate 3x with PBS-T
14. Add 100µL of Streptavidin 1:200
15. Incubate in dark chamber for 30 minutes at room temperature
16. Wash the plate 3x with PBS-T
17. Add 100µL of substrate

1 tablet of OPD (SigmaFast OPD, Sigma, Life Science, EUA) + 1 buffer + 20µL of distilled water → to do 2 plates

18. Incubate in dark chamber for 20 minutes at room temperature
19. Block with 50µL of H₂SO₄ 2N
20. Epoch reading at 492 nm (Test-dependent wavelength)(Epoch-BioTek, Winooski, VT, EUA)

ELISA TNF-alpha (DuoSet – R&D Systems (Minneapolis, MN, EUA))

1. Add 100µL of capture antibody (Work concentration= 4µg/mL) in 96-well plate.

To reach the work concentration, dilute 50µL from the capture antibody aliquot at the concentration of 800µg/mL (with 40µg) in 10mL of Carbonate-Bicarbonate Buffer to obtain a concentration of 4µg/mL.

2. Incubate the plate overnight in a wet chamber at 4°C
3. Wash the plate 3x with PBS-T (250µL/well)
4. Block with 200µL of PBS-T 5% skimmed milk powder
5. Incubate 3h at room temperature
6. Wash the plate 3x with PBS-T
7. To prepare the curve, dilute 25µL from the recombinant aliquot (10000pg/mL) in 225µL of Reagent Diluent (RD) 1% BSA-PBS to obtain the concentration of 1000pg/mL

First point of the curve: 2000pg/mL.

Add 100µL of IL-1 Rec on the 1st and 2nd well.

Add 100µL of RD 1% BSA-PBS on 2nd, 3rd, 4th, 5th, 6th, 7th and 8th well.

Serial dilution from the 2nd well to the last well.

8. Add 100µL of samples
9. Incubate for 3h at room temperature
10. Wash the plate 3x with PBS-T
11. Add 100µL of detection antibody (Work concentration =500ng/mL)

To reach the work concentration, dilute 100µL from the detection antibody aliquot at the concentration of 100µg/mL (with 5µg) in 10mL of RD to obtain a concentration of 500ng/mL

12. Incubate for 2h at room temperature
13. Wash the plate 3x with PBS-T
14. Add 100µL of Streptavidin 1:200
15. Incubate in dark chamber for 30 minutes at room temperature
16. Wash the plate 3x with PBS-T
17. Add 100µL of substrate

1 tablet of OPD (SigmaFast OPD, Sigma, Life Science, EUA) + 1 buffer + 20µL of distilled water → to do 2 plates

18. Incubate in dark chamber for 20 minutes at room temperature
19. Block with 50µL of H₂SO₄ 2N
20. Epoch reading at 492 nm (Test-dependent wavelength)(Epoch-BioTek, Winooski, VT, EUA)