

M9 base medium for bee gut bacteria

Description

This protocol is used to prepare stock solutions for a M9-based 'minimal' medium containing trace metals, vitamins, and casamino acids.

1. Preparation of HMB (hutner's mineral base) modified metals 44 (component 1)

- Weight on a balance the following:
 1. 1.095g $\text{ZnSO}_4 \times 7\text{H}_2\text{O}$, dissolve in **5ml** ddH₂O at pH=2
 2. 0.914g $\text{FeSO}_4 \times 7\text{H}_2\text{O}$, dissolve in **5ml** ddH₂O at pH=2
 3. 0.154g $\text{MnSO}_4 \times \text{H}_2\text{O}$, dissolve in **5ml** ddH₂O at pH=2
 4. 0.392g $\text{CuSO}_4 \times 5\text{H}_2\text{O}$, dissolve in **100ml** ddH₂O at pH=2
 5. 0.248g $\text{Co}(\text{NO}_3)_2 \times 6\text{H}_2\text{O}$, dissolve in **100ml** ddH₂O at pH=2
 6. 0.177g $\text{Na}_2\text{B}_4\text{O}_7 \times 10\text{H}_2\text{O}$, dissolve in **100ml** ddH₂O at pH=2
- Mix **1.**, **2.**, **3.**, and add 10ml of each of **4.**, **5.**, **6.** to obtain **component 1.**
- Adjust volume to 100ml and sterile-filter with 0.22 um filter
- Keep it at 4C, light protected

2. Preparation of additional salts (components 2, 3, 4, 5)

- **Component 2:** 1.72g of $\text{MgSO}_4 \times 7\text{H}_2\text{O}$, dissolve in **100ml** ddH₂O at pH=2, sterile-filter with 0.22 um filter
- **Component 3:** 3.33g of $\text{CaCl}_2 \times 2\text{H}_2\text{O}$, dissolve in **100ml** ddH₂O at pH=2, sterile-filter with 0.22 um filter
- **Component 4:** 0.99g of $\text{FeSO}_4 \times 7\text{H}_2\text{O}$, dissolve in **100ml** ddH₂O at pH=2, sterile-filter with 0.22 um filter
- **Component 5:** 0.974g of $(\text{NH}_4)_6\text{Mo}_7\text{O}_{24} \times 4\text{H}_2\text{O}$, dissolve in **100ml** ddH₂O at pH=2, sterile-filter with 0.22 um filter
- Keep at RT, light protected

3. Preparation of the vitamin 10x stock solution (component 6)

- 1. 5g of L-aspartate potassium, dissolve in **20ml** ddH₂O
- 2. 0.2g of Calcium pantothenate, dissolve in **10ml** ddH₂O
- 3. 0.01g of Thiamine-HCl, dissolve in **10ml** ddH₂O
- 4. 0.01g of biotin, dissolve in **10ml** ddH₂O
- Mix 59ml of ddH₂O, 20ml of **1.**, 10ml of **2.**, 10ml of **3.**, and 1ml of **4.**
- Sterile-filter the final 100 ml of vitamin stock solution (component 6)
- Component 6 can be kept for at least 6 months at 4°C in the fridge, protected from light

4. Preparation of 10x stock solution of casamino acids (component 7)

- 200g casamino acids (BD) , dissolve in **1000ml** ddH₂O, autoclave
- keep at RT

5. Preparation of 10x stock solution of M9 (component 8)

- Mix the following:
- 60g Na_2HPO_4
- 30g KH_2PO_4
- 5g NaCl
- 10g NH_4Cl
- Ad water to 1l, autoclave
- Keep at RT

6. Prepare 500ml of M9 base medium

- Prepare the final M9 base medium freshly on the day of use
- For 500ml mix the following:
- Component 1: 500 μl
- Component 2: 8390 μl
- Component 3: 1000 μl
- Component 4: 100 μl
- Component 5: 10 μl
- Component 6 (10x Vitamin stock): 50ml
- Component 7 (10x Casamino acids stock): 50ml
- Component 8 (10x M9 stock): 50ml
- ddH₂O: 340 ml

Changes/notes

Version	Changes	date
	.	