

Protein assessment by the Bradford method

Material:

1. Standard protein
2. Micro tubes
3. Bradford's ready-to-use solution

Procedures:

1. Prepare 3 to 5 dilutions of standard protein (BSA). The linearity for the BSA assay is 0.1 to 0.9 mg / mL. Pipette 2 μ L of each standard and sample, in triplicate, into the micro tubes.
2. Add 200 μ L of Bradford's ready-to-use solution and take the tubes for shaking in the vortex.
3. Incubate at room temperature for at least 5 minutes.
4. Read the absorbance in a spectrophotometer at a wavelength of 595 nm.
5. The protein concentration was calculated by using a standard curve dilution of 1 mg of BSA/mL with 1, 3, 5, 7, and 10 μ g.