

1. Follow Instructions in Step 64 of Potency Test: Glucose Stimulated Insulin Release Assay

<https://www.protocols.io/edit/potency-test-glucose-stimulated-insulin-release-as-bc3siyne>

UNOS Number: _____ Isolation Number: _____ Date of Assay: _____ Performed By: _____

	1	2	3	4	5	6	7	8	9	10	11	12
A	STD 0 0.00-0.00	STD 0 0.00-0.00	STD 0 0.00-0.00	KREBS Low Glucose	KREBS Low Glucose	KREBS Low Glucose						
B	STD 3 2.55-3.45	STD 3 2.55-3.45	STD 3 2.55-3.45	KREBS High Glucose	KREBS High Glucose	KREBS High Glucose						
C	STD 10 8.5-11.5	STD 10 8.5-11.5	STD 10 8.5-11.5									
D	STD 30 25.5-34.5	STD 30 25.5-34.5	STD 30 25.5-34.5									
E	STD 100 85-115	STD 100 85-115	STD 100 85-115									
F	STD 200 170-230	STD 200 170-230	STD 200 170-230									
G	LOW CONTROL	LOW CONTROL	LOW CONTROL									
H	HIGH CONTROL	HIGH CONTROL	HIGH CONTROL									
	Kit Lot#: _____			Temp. (Room): _____		Low Control Lot #: _____			High Control Lot#: _____			
	Exp. Date: _____			(acceptable range: 18-25°C)		Exp. Date: _____			Exp. Date: _____			