



Standard Operating Procedure for Potency Test: Glucose Stimulated Insulin Release Assay	<i>SOP #:</i>
<i>Details:</i> Attachment 3: Glucose Stimulated Insulin Release Worksheet (GSIR)	<i>Page Number:</i> Page 1 of 1

1. Glucose Stimulated Insulin Release Assay (Follow instructions in Step 66 and 77 of SOP)

<https://www.protocols.io/edit/potency-test-glucose-stimulated-insulin-release-as-bc3siyne>

UNOS Number: _____ Isolation Number: _____

2. Stimulation Test:

Performed by: _____ Date of Assay: _____

Low Glucose Start Time: _____ Low Glucose Stop Time: _____ Total Stimulation Time: _____

High Glucose Start Time: _____ High Glucose Stop Time: _____ Total Stimulation Time: _____

3. ELISA: (If using RIA, enter only who performed the assay and the date of assay.)

Performed by: _____ Date of Assay: _____ Shaker RPM: _____ Temperature: _____

Enzyme Conjugate Start Time: _____ Enzyme Conjugate Stop Time: _____ EC Incubation Time: _____

Substrate TMB Start Time: _____ Substrate TMB Stop Time: _____ Substrate TMB Incubation Time: _____

4. ELISA Results: (If using RIA, enter only Low/High Insulin Averages and Stimulation Index.)

Low Glucose Optical Density Average: _____ Low Glucose Insulin mU/L Average: _____

High Glucose Optical Density Average: _____ High Glucose Insulin mU/L Average: _____

Stimulation Index: _____ (Average High Insulin/Average Low Insulin)