

Supplemental Material

Table of Contents

Figure S1	2
Figure S2	3
Figure S3	4
Figure S4	5
Figure S5	6
Figure S6	7
Figure S7	8
Figure S8	9
Figure S9	10
Figure S10	12
Figure S11	13
Figure S12	14
Figure S13	15
Table S1	16
Table S2	18
Table S3	19
Table S4	20
Table S5	22
Table S6	24
Table S7	25
Table S8	26
Table S9	28
Table S10	29
Table S11	30
Table S12	31
Table S13.....	32
Supplementary Data Sheet 1. Protocol 1: <i>Pf</i> -SMART Amplification	35
Supplementary Data Sheet 2. Protocol 2: Library preparation	39
Data Availability	43




Figure S1. PCR estimated parasitemias based on real-time PCR. Panel A shows the distribution of estimated parasitemia in the entire study population. Panel B shows the estimated parasitemias of the 100 selected for this sequencing study. Panel C shows the estimated parasitemia after re-extraction. Red bars show mean and whiskers denote maximum and minimum parasitemia. Parasitemia is shown on a log₁₀ scale.




Figure S2. Sequencing depth for assay development controls. Plots depict sequencing depth on a log10 scale for 12 samples at each parasitemia. At each parasitemia level, individual points for values are plotted (left) and a box and whiskers plot with maximum and minimum values (right).




Figure S3. Number of successfully genotyped amplicons in controls at different parasitemias.

Successfully genotyped amplicons are defined as amplicons that retrieve at least one haplotype within a given replicate. At each parasitemia level, individual points for values are plotted (left) and a box and whiskers plot with maximum and minimum values (right).




Figure S4. Haplotype counts in monoclonal control samples at different parasitemias for diversity amplicons. This figure shows the relationship between parasitemia and the number of haplotypes detected for amplicons that are markers of diversity. There are 12 samples at each parasitemia, and because the samples are monoclonal, the expected number of haplotypes is one.




Figure S5. Haplotype counts in monoclonal control samples at different parasitemias for drug resistance amplicons. This figure shows the relationship between parasitemia and the number of haplotypes detected for amplicons containing known resistance loci. There are 12 samples at each parasitemia, and because the samples are monoclonal, the expected number of haplotypes is one.




Figure S6. Control samples with unexpected haplotypes. A total of 4 replicate/amplicon/parasitemia combinations (of 1,440 genotyped combinations) had unexpected haplotypes detected (red). All of the replicates with unexpected haplotypes occurred at a concentration of 1 parasite per microliter. The expected haplotype was detected in all samples as the majority variant.




Figure S7. Within sample allele frequency of known mixture of 3D7 and Dd2 (88:12). The amplicon specific within sample allele frequency for the Dd2 minor variant are shown for all 12 samples. Only amplicons that are variable between the two strains are shown. The expected allele frequency is 12 percent. Strain Dd2 has an imperfect second copy of MDR1 which alters the allele frequency at these amplicons (both the frequencies of individual copies and the sum of the frequencies of both copies are shown). At each amplicon, individual points for values are plotted (left) and a box and whiskers plot with maximum and minimum values (right).




Figure S8. Number of successfully genotyped amplicons per sample among 100 Cameroonian isolates.

Panel A shows the number of successful amplicons for each sample, where a successful amplicon is one containing at least one haplotype. The red bar is the mean for all 100 samples. Panel B shows the parasitemia distribution of samples with complete (24 amplicons) versus non-complete (<24 amplicons) genotyping on a log₁₀ scale. The line in the box displays the median parasitemia and the whiskers show the minimum and maximum parasitemia (ns, $p = 0.07$, t-test).




Figure S9. Within sample allele frequency comparison of replicates by amplicon for 100 Cameroonian isolates. The calculated within sample haplotype frequency (diversity markers) and antimalarial single nucleotide polymorphisms (SNP) are shown for each replicate. Black dots represent samples with two replicates. For samples where multiple haplotypes were detected for resistance genes, the lowest frequency haplotype was reported.




Figure S10. Read depth at each amplicon for 100 Cameroonian isolates. Each amplicon contains 100 Cameroonian isolates. Each isolate contains two replicates which are averaged to yield a single amplicon read count per isolate, illustrated by 100 points per amplicon. At each amplicon, individual points for values are plotted (left) and a box and whiskers plot with maximum and minimum values (right).





Figure S12. Principal Component Analysis of 50 Cameroonian Samples using *Pf*-SMART (A) and Molecular Inversion Probes (MIP)

A



B



Figure S13. Allele frequency in pooled samples. Weighted mean allele frequency for key resistance mutations shows consistency between pooled and individual samples. Exceptions occur for highly variable and low-frequency variants such as mutations K540E and A581G with frequencies varying from 0.0-0.7% and 0.1-51.5%, respectively, that are not found in pooled samples (A). Lower consistency between the number of total mutations found in pooled versus individual samples can be seen in variants that occur with a frequency lower than 40%, whereas more frequent mutations are successfully identified in both pooled and individual samples (B).

Table S1. Comparison of population characteristics

	Not sequenced (n=130)	Sequenced (n=100)	Total (N=230)	P-value
Age (years)	28.0 (21.0, 45.0)	25.0 (19.0, 35.0)	26.0 (20.0, 40.0)	0.28
Sex?				0.09
Male	53 (41%)	52 (52%)	105 (46%)	
Female	77 (59%)	48 (48%)	125 (54%)	
Education				0.35
Primary	55 (42%)	33 (33%)	88 (38%)	
Secondary	20 (15%)	17 (17%)	37 (16%)	
University	55 (42%)	50 (50%)	105 (46%)	
Occupation				0.60
Self-employed	36 (72%)	23 (82%)	59 (76%)	
Employer	5 (10%)	2 (7%)	7 (9%)	
Not working	9 (18%)	3 (11%)	12 (15%)	
Headache?				0.72
No	90 (69%)	67 (67%)	157 (68%)	
Yes	40 (31%)	33 (33%)	73 (32%)	
General fatigue?				0.56
No	93 (72%)	68 (68%)	161 (70%)	
Yes	37 (28%)	32 (32%)	69 (30%)	
Vomiting?				0.59
No	122 (94%)	92 (92%)	214 (93%)	
Yes	8 (6%)	8 (8%)	16 (7%)	
Cough?				0.59
No	122 (94%)	92 (92%)	214 (93%)	
Yes	8 (6%)	8 (8%)	16 (7%)	
Abdominal pain?				0.19
No	90 (69%)	61 (61%)	151 (66%)	
Yes	40 (31%)	39 (39%)	79 (34%)	
Gastric pain?				0.21
No	128 (98%)	100 (100%)	228 (99%)	
Yes	2 (2%)		2 (1%)	
Other symptoms?				0.47
No	122 (94%)	96 (96%)	218 (95%)	
Yes	8 (6%)	4 (4%)	12 (5%)	
Slept under a bednet in the last 24 hours?				0.11
No	55 (42%)	32 (32%)	87 (38%)	
Yes	75 (58%)	68 (68%)	143 (62%)	
Nights in the last 7 slept under a bednet?	7.0 (0.0, 7.0)	7.0 (0.0, 7.0)	7.0 (0.0, 7.0)	0.35
Number of people in the household?	2.5 (2.0, 5.0)	3.0 (2.0, 5.0)	3.0 (2.0, 5.0)	0.42
Number of bednets in the household?	1.0 (0.0, 2.0)	1.0 (0.0, 2.0)	1.0 (0.0, 2.0)	0.73
Water within a 2 minutes walk of the household?				0.17
No	102 (78%)	79 (79%)	181 (79%)	
Yes	21 (16%)	20 (20%)	41 (18%)	
Don't know	7 (5%)	1 (1%)	8 (3%)	0.17
Type of water source?				0.28
Stream	3 (2%)	6 (6%)	9 (4%)	
Pond/lake	9 (7%)	7 (7%)	16 (7%)	
Swamp/marsh	9 (7%)	7 (7%)	16 (7%)	
None	102 (78%)	79 (79%)	181 (79%)	
Don't know	7 (5%)	1 (1%)	8 (3%)	
Traveled in the last 3 months?				0.34
No	78 (84%)	53 (78%)	131 (81%)	

Yes	15 (16%)	15 (22%)	30 (19%)	
Is a member of the household working in the forest?				0.85
No	126 (99%)	97 (99%)	223 (99%)	
Yes	1 (1%)	1 (1%)	2 (1%)	
Animals in the house?				0.02
No	120 (94%)	84 (86%)	204 (91%)	
Yes	7 (6%)	14 (14%)	21 (9%)	
Pregnant?				0.96
No	122 (96%)	94 (96%)	216 (96%)	
Yes	5 (4%)	4 (4%)	9 (4%)	
Taken antimalarials in the last 14 days?				0.63
No	93 (73%)	70 (72%)	163 (73%)	
Yes	28 (22%)	24 (25%)	52 (23%)	
Don't know	4 (3%)	3 (3%)	7 (3%)	
Refused	2 (2%)		2 (1%)	
Which antimalarials?				0.51
Quinine	5 (15%)	4 (15%)	9 (15%)	
ACT	16 (47%)	11 (41%)	27 (44%)	
Don't know	11 (32%)	12 (44%)	23 (38%)	
Refused	2 (6%)		2 (3%)	
Source of antimalarials?				0.56
Drug seller	6 (22%)	8 (35%)	14 (28%)	
Clinic	14 (52%)	11 (48%)	25 (50%)	
Leftover	7 (26%)	4 (17%)	11 (22%)	
Lifelong malaria exposure?				0.34
No	9 (7%)	4 (4%)	13 (6%)	
Yes	118 (93%)	94 (96%)	212 (94%)	
Number of malaria exposures over lifetime?	7.0 (3.0, 7.0)	7.0 (6.5, 7.0)	7.0 (5.0, 7.0)	0.51
Number of malaria exposures in the last year?	0.0 (0.0, 1.0)	0.0 (0.0, 1.0)	0.0 (0.0, 1.0)	0.62
Has someone in the household been treated for malaria in the last year?				0.77
No	14 (11%)	12 (12%)	26 (12%)	
Yes	104 (82%)	76 (78%)	180 (80%)	
Don't know	8 (6%)	8 (8%)	16 (7%)	
Refused	1 (1%)	2 (2%)	3 (1%)	

P-value comparisons across groups for categorical variables are based on chi-square test of homogeneity; p-values for continuous variables are based on Kruskal-Wallis test for median.

Table S2. Within sample allele frequency of known mixture of 3D7 and Dd2 (88:12)^

	AMA1	Heome B	Heome C	Heome E	Heome F	Heome G	DHFR 51/59	DHFR 108	DHPS 436/437*	DHPS 613	MDR1 86 (HAP 1)	MDR1 86 (HAP2)	MDR1 86 (TOTAL)	CRT
Average	0.119	0.094	0.105	0.097	0.087	0.107	0.096	0.088	0.011	0.093	0.118	0.068	0.185	0.085
STDEV	0.042	0.030	0.040	0.035	0.036	0.047	0.036	0.030	0.001	0.023	0.032	0.017	0.038	0.039

*: when detected (5/12 reps)

^: Not all amplicons are variable between the strains. Those with no polymorphisms are not reported.

Table S3. Individual characteristics and sequencing results for 100 participants sequenced with *Pf*-SMARRT

See the uploaded Excel file.

- Tab 1: Individual metadata and averaged frequencies for antimalarial resistance polymorphisms and diversity amplicons in each individual
- Tab 2: Data dictionary for metadata
- Tab 3: Underlying sequencing depth (count) data for sequencing replicates for each individual used to calculate within sample allele frequencies

Table S4. Haplotypes detected for AMA in 100 Cameroonian participants

Haplotype Name	Haplotype Sequence
ama1.00	ATTTGGTAAAGGTATAATTATTGAGAATTCAAAACTACTTTTTTAACACCGGTAGCTACGGAAAAATCAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAAATCCTCCTAT GTCACCAATGACATTAAATGGTATGAGAGATTTATATAAAAAATAATGAATATGTAAAAAATTTAGATGAATTGACTTTA
ama1.01	ATTTGGTAAAGGTATAATTATTGAGAATTCAAATACTACTTTTTTAACACCGGTAGCTACGGAAAAATCAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAAATCCTCTTAT ATCACCAATGACATTAGATCATATGAGAGATTCCTATAAAAAATAATGAATATGTAAAAAATTTAGATGAATTGACTTTA
ama1.02	ATTTGGTAAAGGTATAATTATTGAGAATTCAAATACTACTTTTTTAACACCGGTAGCTACGGGAAATCAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAAATCCTCTTAT ATCACCAATGACATTAAATGGTATGAGAGATTTTATAAAAAATAATGAATATGTAAAAAATTTAGATGAATTGACTTTA
ama1.03	ATTTGGTAAAGGTATAATTATTGAGAATTCAAATACTACTTTTTTAACACCGGTAGCTACGGAAAAATCAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAAAACCTCTTAT GTCACCAATGACATTAGATCAAATGAGACATTTTATAAAGATAATAAATATGTAAAAAATTTAGATGAATTGACTTTA
ama1.04	ATTTGGTAAAGGTATAATTATTGAGAATTCAAATACTACTTTTTTAACACCGGTAGCTACGGGAAATCAATATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAGAACCTCATAT GTCACCAATGACATTAGATGAAATGAGACATTTTATAAAGATAATAAATATGTAAAAAATTTAGATGAATTGACTTTA
ama1.05	ATTTGGTAAAGGTATAATTATTGAGAATTCAAATACTACTTTTTTAACACCGGTAGCTACGGGAAATCAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAGAACCTCTTAT ATCACCAATGACATTAGATGATATGAGAGATTTTATAAAAAATAATGAATATGTAAAAAATTTAGATGAATTGACTTTA
ama1.06	ATTTGGTAAAGGTATAATTATTGAGAATTCAAATACTACTTTTTTAACACCGGTAGCTACGGAAAAATCAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAGAACCTCTTAT GTCACCAATGACATTAGATCAAATGAGACATTTTATAAAGATAATAAATATGTAAAAAATTTAGATGAATTGACTTTA
ama1.07	ATTTGGTAAAGGTATAATTATTGAGAATTCAAATACTACTTTTTTAACACCGGTAGCTACGGAAAAATCAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAAAACCTCTTAT GTCACCAATGACATTAGATCAAATGAGAGATTTTATAAAAAATAATGAATATGTAAAAAATTTAGATGAATTGACTTTA
ama1.08	ATTTGGTAAAGGTATAATTATTGAGAATTCAAATACTACTTTTTTAACACCGGTAGCTACGGAAAAATCAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAAAACCTCTTAT GTCACCAATGACATTAGATGAAATGAGACATTTTATAAAGATAATAAATATGTAAAAAATTTAGATGAATTGACTTTA
ama1.09	ATTTGGTAAAGGTATAATTATTGAGAATTCAAATACTACTTTTTTAACACCGGTAGCTACGGGAAATCAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAGAACCTCTTAT ATCACCAATGACATTAAATGGTATGAGAGATTTTATAAAAAATAATGAATATGTAAAAAATTTAGATGAATTGACTTTA
ama1.10	ATTTGGTAAAGGTATAATTATTGAGAATTCAAATACTACTTTTTTAACACCGGTAGCTACGGGAAATCAATATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAGAACCTCTTAT GTCACCAATGACATTAGATGAAATGAGACATTTTATAAAGATAATAAATATGTAAAAAATTTAGATGAATTGACTTTA
ama1.11	ATTTGGTAAAGGTATAATTATTGAGAATTCAAATACTACTTTTTTAACACCGGTAGCTACGGGAAACAAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAAATCCTCTTAT ATCACCAATGACATTAAATGGTATGAAAGATTTTATAAAGATAATGAAGATGTAAAAAATTTAGATGAATTGACTTTA
ama1.12	ATTTGGTAAAGGTATAATTATTGAGAATTCAAAACTACTTTTTTAACACCGGTAGCTACGGAAAAATCAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAAAACCTCTTAT GTCACCAATGACATTAGATGATATGAGAGATCCTTATAAAAAATAATGAATATGTAAAAAATTTAGATGAATTGACTTTA
ama1.13	ATTTGGTAAAGGTATAATTATTGAGAATTCAAAACTACTTTTTTAACACCGGTAGCTACGGAAAAATCAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAGAACCTCTTA TGTCACCAATGACATTAGATGATATGAGACGTTTTTATAAAGATAATGAATATGTAAAAAATTTAGATGAATTGACTTTA
ama1.14	ATTTGGTAAAGGTATAATTATTGAGAATTCAAATACTACTTTTTTAACACCGGTAGCTACGGGAAACAAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAAATCCTCTTAT ATCACCAATGACATTAGATCATATGAGAGATTTTATAAAAAAATGAATATGTAAAAAATTTAGATGAATTGACTTTA
ama1.15	ATTTGGTAAAGGTATAATTATTGAGAATTCAAATACTACTTTTTTAACACCGGTAGCTACGGGAAATCAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAAAACCTCTTAT GTCACCAATGACATTAGATGATATGAGACTTTTGTATAAAGATAATGAAGATGTAAAAAATTTAGATGAATTGACTTTA
ama1.16	ATTTGGTAAAGGTATAATTATTGAGAATTCAAAACTACTTTTTTAACACCGGTAGCTACGGAAAAATCAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAAAACCTCTTAT GTCACCAATGACATTAGATGAAATGAGACATTTTATAAAGATAATAAATATGTAAAAAATTTAGATGAATTGACTTTA
ama1.17	ATTTGGTAAAGGTATAATTATTGAGAATTCAAATACTACTTTTTTAACACCGGTAGCTACGGAAAAATCAAGATTTAAAAGATGGAGGTTTTGCTTTTCCTCCAACAGAACCTCTTAT GTCACCAATGACATTAGATCGTATGAGAGATTTTATAAAAAATAATGAAGATGTAAAAAATTTAGATGAATTGACTTTA

Table S5. Haplotypes detected for Heome A in 100 Cameroonian participants

Haplotype Name	Haplotype Sequence
heome-a.00	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTTT GTTTCATCTTTGAATTATCCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.01	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTTT GTTTCATCTTTGAATTATCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.02	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTCT GTTTCATCTTTGAATTATCCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.03	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTTT GTTTCATCTTTGAATTATCCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.04	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTTTG TTCATCTTTGAATTATCCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.05	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTAACTTTG TTCATCTTTGAATTATCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.06	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTTT GTTTCATCTTTGAATTATCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.07	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTTTT TTCATCTTTGAATTATCCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.08	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTCT GTTTCATCTTTGAATTATCCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.09	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTTTG TTCATCTTTGAATTATCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.10	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTCT GTTTCATCTTTGAATTATCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.11	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTTT GTTTCATCTTTGAATTAACCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.12	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTCT GTTTCATCTTTGAATTATCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.13	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTTT GTTAATCTTTGAATTATCCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.14	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTCT GTTTCATCTTTGAATTATCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.15	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTTTG TTCATCTTTGAATTATCCGTTCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.16	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTTTT TTCATCTTTGAATTATCCATCTGTTTGTATGTCTCAAGCAAGGTT
heome-a.17	TTATTAATGTCATTTTGTTTTCTTTATTATTATTATTTTATTTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTATTTTGACTTTG TTCATCTTTGAATTAACCGTTCTGTTTGTATGTCTCAAGCAAGGTT

heome-a.18	TTATTAATGTCATTTTGTTCCTTATTATTATTTTATTTTCCCTGGTTATCTTCCTTATTACCTAGCTTTTTTTGTTGTCATTATTTCTGTCTGATTATTTTGACTTT GTTTCATCTTTGAATTATCCGTTCTGTTTGAATGTCTCAAGCAAGGTT
------------	---

Table S6. Haplotypes detected for Heome B in 100 Cameroonian participants

Haplotype Name	Haplotype Sequence
heome-b.00	AATGTACTTATCCCCTGATATGATTTTCCAATATCATCAAAATCAGTAATCTCTAAATCATTTATATGATTATATGTATAAATCTTTTATCTTGTTCAAATAATTTCTGTAATAAAATA TATTCTCTTTATCATCATCAATTATTTTCATATAATCC
heome-b.01	AATGTACTTATCCCCTGATATGATTTCTCCAGTATCATCAAAATCAGTAATCTCTAAATCATTTATATGATTATATGTATAAATCTTTTATCTTGTTCAAATAATTTCTGTAATAAAATA TATTCTCTTTATCATCATCAATTATTTTCATATAATCC
heome-b.02	AATGTACTTATCCCCTGATATGATTCTCCAATATCATCAAAATCAGTAATCTCTAAATCATTTATATGATTATATGTATAAATCTTTTATCTTGTTCAAATAATTTCTGTAATAAAATA TATTCTCTTTATCATCATCAATTATTTTCATATAATCC
heome-b.03	AATGTACTTATCCCCTGATATGATTGTTCCAATATCATCAAAATCAGTAATCTCTAAATCATTTATATGATTATATGTATAAATCTTTTATCTTGTTCAAATAATTTCTGTAATAAAATA TATTCTCTTTATCATCATCAATTATTTTCATATAATCC
heome-b.04	AATGTACTTATCCCCTGATATGATTTCCCAATATCATCAAAATCAGTAATCTCTAAATCATTTATATGATTATATGTATAAATCTTTTATCTTGTTCAAATAATTTCTGTAATAAAATA TATTCTCTTTATCATCATCAATTATTTTCATATAATCC
heome-b.05	AATGTACTTATCCCCTGATATGATTCTCCAATATCATCAAAATAGTAATCTCTAAATCATTTATATGATTATATGTATAAATCTTTTATCTTGTTCAAATAATTTCTGTAATAAAATA TATTCTCTTTATCATCATCAATTATTTTCATATAATCC
heome-b.06	AATGTACTTATCCCCTGATATGATTTTCCAATATCATCAAAATAGTAATCTCTAAATCATTTATATGATTATATGTATAAATCTTTTATCTTGTTCAAATAATTTCTGTAATAAAATA TATTCTCTTTATCATCATCAATTATTTTCATATAATCC
heome-b.07	AATGTACTTATCCCCTGATATGATTTCTCCAATATCATCAAAATCAGTAATCTCTAAATCATTTATATGATTATATGTATAAATCTTTTATCTTGTTCAAATAATTTCTGTAATAAAATA TATTCTCTTTATCATCATCAATTATTTTCATATAATCC
heome-b.08	AATGTACTTATCCCCTGATATGATTTCTCCATTATCATCAAAATCAGTAATCTCTAAATCATTTATATGATTATATGTATAAATCTTTTATCTTGTTCAAATAATTTCTGTAATAAAATA TATTCTCTTTATCATCATCAATTATTTTCATATAATCC
heome-b.09	AATGTACTTATCCCCTGATATGAGTTCTCCAATATCATCAAAATCAGTAATCTCTAAATCATTTATATGATTATATGTATAAATCTTTTATCTTGTTCAAATAATTTCTGTAATAAAAT ATATTCTCTTTATCATCATCAATTATTTTCATATAATCC

Table S7. Haplotypes detected for Heome C in 100 Cameroonian participants

Haplotype Name	Haplotype Sequence
heome-c.0	TACCCACGGGACATTTTAATTTATACCTATTGGAACATTTTTATTAGGTCTTTCTTTATGATTCTTCTTTAAATAGGTTTCATCATTATTTGAGTAACTAACAGATCATTATCGTTCT CAAAATTATTATAACCATCAATATTTTATTTCCATATACATCTATAAT
heome-c.1	TACCCACGGGACATTTTAATTTATACCTATTGGAACATTTTTATTAGGTCTTTCTTTATGATTCTTCTTTAAATAGGTTTCATCATTATTTGAGTAACTAACAGATTATTATCGTTCT CAAAATTATTATAACCATCAATATTTTATTTCCATATACATCTATAAT
heome-c.2	TACCCACGGGACATTTTAATTTATACCTATTGGAACATTTTTATTAGGTCTTTCTTTATGATTCTTCTTTAAATAGGTTTCATCATTATTTGAGTAACTAACAGATCATTATCGTTCT CAAAATTATTATAACCATCAATATTTTATTTCCATATACATCTATAAT
heome-c.3	TACCCACGGGACATTTTAATTTATACCTATTGGAACATTTTTATTAGGTCTTTCTTTATGATTCTTCTTTAAATAGGTTTCATCATTATTTGAGTAACTAACAGATCATTATGGTTC TCAAAATTATTATAACCATCAATATTTTATTTCCATATACATCTATAAT
heome-c.4	TACCCACGGGACATTTTAATTTATACCTATTGGAACATTTTTATTAGGTCTTTCTTTATGATTCTTCTTTAAATAGGTTTCATCATTATTTGAGTAACTAACAGATCATTATCCTTCT CAAAATTATTATAACCATCAATATTTTATTTCCATATACATCTATAAT

Table S8. Haplotypes detected for Heome D in 100 Cameroonians participants

Haplotype Name	Haplotype Sequence
heome-d.00	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGAATACACCAG
heome-d.01	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCATGAATACACCAG
heome-d.02	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGAATACACCAG
heome-d.03	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAACATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGAATACACCAG
heome-d.04	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAACATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGAATACACCAG
heome-d.05	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAACGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGAATACACCAG
heome-d.06	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATACATCTTGAATACACCAG
heome-d.07	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGTATACACCAG
heome-d.08	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAACATTTAT TTTTGAGAAAGATTTCATATATTCATATACATCTTGAATACACCAG
heome-d.09	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCGTGAATACACCAG
heome-d.10	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAAAATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGAATACACCAG
heome-d.11	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACCTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGAATACACCAG
heome-d.12	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGAATACACCAG
heome-d.13	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTGTTATTTTTATTTTAAAAAAGGATATATCATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGAATACACCAG
heome-d.14	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAAAAAGATTTCATATATTCATATACATCTTGAATACACCAG
heome-d.15	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACCTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAAAAAGATTTCATATATTCATATATATCTTGAATACACCAG
heome-d.16	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTCGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGAATACCCAG
heome-d.17	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTAAATTGATAGTTTGCAACAATTATACGTTTGTTATTTTTATTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCATGAATACACCAG

heome-d.18	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTTAAATTGATAGTTTGCAAACAATTATACGTTTGTTATTTTTATTTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTAAGAAAGATTTCATATATTCATATATATCTTGAATACACCAG
heome-d.19	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTTAAATTGATAGTTTGCAAACAATTATACGTTTGTTATTTTTATTTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTAAGAAAGATTTCATATATTCATATATATCATGAATACACCAG
heome-d.20	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTTAAATTGATAGTTTGCAAACAATTATACGTTTGTTATTTTTATTTTTGAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGGATACACCAG
heome-d.21	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTTAAATTGATAGTTTGCAAACAATTATACGTTTGTTATTTTTATTTTTAAAAAAGGATATATGATCATCTAAACATTTAT TTTTAAGAAAGATTTCATATATTCATATATATCTTGAATACACCAG
heome-d.22	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTTAAATTGATAGTTTGCAAACAATTATACGTTTGTTATTTTTATTTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCATGAATACACCAG
heome-d.23	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTTAAATTGATAGTTTGCAAACAATTATACGTTTCGTTATTTTTATTTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGAATACACCAG
heome-d.24	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTTAAATTGATAGTTTGTAACAATTATACGTTTGTTATTTTTATTTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGGATACACCAG
heome-d.25	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTTAAATTGATAGTTTGCAAACAATTATACGTTTCGTTATTTTTATTTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGGATACACCAG
heome-d.26	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTTAAATTGATAGTTTGCAAACAATTATACCTTGTTATTTTTATTTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCATGAATACACCAG
heome-d.27	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTTAAATTGATAGTTTGCAAACAATTATACGTTTGTTATTTTTATTTTTAAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCTTGGATACACCAG
heome-d.28	CATGATTCCTCAAAAGTTTAACATTCCACATAGTATCTGTTAAATTGATAGTTTGCAAACAATTATACGTTTGTTATTTTTATTTTTGAAAAAGGATATATGATCATCTAAAGATTTAT TTTTGAGAAAGATTTCATATATTCATATATATCGTGAATACACCAG

Table S9. Haplotypes detected for Heome E in 100 Cameroonian participants

Haplotype Name	Haplotype Sequence
heome-e.0	GTATCTTTCCATATATAATAACGCATGAATTCTTTAAATATATCTTTACTTTATCTGTAAGTTTCATGGAATCTATTATTTTCATGTGCATATTTATTTTTGCCTCGTTTATATATTAC CTGCTTCTAATTTACATAAACATTTTCTATTAAGTTTGA
heome-e.1	GTATCTTTCCATATATAATAACGCATGAATTCTTTAAATATATCTTTACTTTATCTGTAAGTTTCATGGAATCTATTATTTTCATGTGCATATTTATTTTTGCCTCGTTTATATATTAC CTGCTTTTAATTTACATAAACATTTTCTATTAAGTTTGA
heome-e.2	GTATCTTTCCATATATAATAACGCATGAATTCTTTAAATATATCTTTACTTTATCTGTAAGTTTCATGGAATCTATTATTTTCATGTGCATATTTATTTTTGCCTCGTTTATATATTAC CTGCTTCTAATTTACATAAACATTTTCTATTAAGTTTGA
heome-e.3	GTATCTTTCCATATATAATAACGCATGAATTCTTTAAATATATCTTTACTTTATCTGTAAGTTTCATGGAATCTATTATTTCTTCTGCATATTTATTTTTGCCTCGTTTATATATTAC CTGCTTCTAATTTACATAAACATTTTCTATTAAGTTTGA
heome-e.4	GTATCTTTCCATATATAATAACGCATGAATTCTTTAAATATATCTTTACTTTATCTGTAAGTTTCATGGAATCTATTATTTTCATGTGCATATTTATTTTTGCCTCGTTTATATATTAC CTGCTTTTAATTTACATAAACATTTTCTATTAAGTTTGA
heome-e.5	GTATCTTTCCATATATAATAACGCATGAATTCTTTAAATTTATCTTTACTTTATCTGTAAGTTTCATGGAATCTATTATTTTCATGTGCATATTTATTTTTGCCTCGTTTATATATTAC CTGCTTCTAATTTACATAAACATTTTCTATTAAGTTTGA
heome-e.6	GTATCTTTCCATATATAATAACGCATGAATTCTTTAAATATATCTTTACTTTATCTGTAAGTTTCATGGAATCTATTATTTCTTCTGCATATTTATTTTTGCCTCGTTTATATATTAC CTGCTTTTAATTTACATAAACATTTTCTATTAAGTTTGA
heome-e.7	GTATCTTTCCATATATAATAACGCATGAATTCTTTAAATATATCTTTACTTTATCTGTAAGTTTCATGGAATCTATTATTTTCATGTGCATATTTATTTTTGCCTCGTTTATATATTAC CTGCTTTTAATTTACATAAACATTTTCTATTAAGTTTGA

Table S10. Haplotypes detected for Heome F in 100 Cameroonian participants

Haplotype Name	Haplotype Sequence
heome-f.00	TTTGGGTAGATATCGTTATAAGGGCTATAATGTATTGACTTCTCCGGGTTTGTACATTATCATTACAACCAATATTATTAAGGTTAATATTATTTTATCATTATTATTGTTGTTTTCAT TTGTGTGGTTATATAATAACAATGATGTGTCGCTATTGGATCATTT
heome-f.01	TTTGGGTAGATATCGTTATAAGGGCTATAATGTATTGACTTCTCCGGGTTGGTTACATTATCATTACAACCAATATTATTAAGGTTAATATTATTTTATCATTATTATTGTTGTTTTCAT TTGTGTGGTTATATAATAACAATGATGTGTCGCTATTGGATCATTT
heome-f.02	TTTGGGTAGATATCGTTATAAGGGCTATAATGTATTGACTTCTCCGGGTTGGTTACATTATCATTACAACCAATATTATTAAGGTTTCATATTATTTTATCATTATTATTGTTGTTTTCAT TTGTGTGGTTATATAATAACAATGATGTGTCGCTATTGGATCATTT
heome-f.03	TTTGGGTAGATATCGTTATAAGGGCTATAATGTATTGACTTCTCCGGGTTTGTACATTATCATTACAACCAATATTATTAAGGTTTCATATTATTTTATCATTATTATTGTTGTTTTCAT TTGTGTGGTTATATAATAACAATGATGTGTCGCTATTGGATCATTT
heome-f.04	TTTGGGTAGATATCGTTATAAGGGCTATAATGTATTGACTTCTCCGGGTTTGTACATTATCATTACAACCAATATTATTAAGGTTAATACTATTTTATCATTATTATTGTTGTTTTCAT TTGTGTGGTTATATAATAACAATGATGTGTCGCTATTGGATCATTT
heome-f.05	TTTGGGTAGATATCGTTATAAGGGCTATAATGTATTGACTTCTCCGGGTTGGTTACATTATCATTACAACCAATATTATTAAGGTTAATACTATTTTATCATTATTATTGTTGTTTTCAT TTGTGTGGTTATATAATAACAATGATGTGTCGCTATTGGATCATTT
heome-f.06	TTTGGGTAGATATCGTTATAAGGGCTATAATGTATTGACTTCTCCGGGTTTGTACATTATCATTACAACCAATATTATTAAGGTTTCATATTATTTTATCATTATTATTCTTGTGTTTTCAT TTGTGTGGTTATATAATAACAATGATGTGTCGCTATTGGATCATTT
heome-f.07	TTTGGGTAGATATCGTTATAAGGGCTATAATGTATTGACTTCTCCGGGTTTGTACATTATCATTACAACCAATATTATTAAGGTTAATATTATTTTACCATTATTATTGTTGTTTTCAT TTGTGTGGTTATATAATAACAATGATGTGTCGCTATTGGATCATTT
heome-f.08	TTTGGGTAGATATCGTTATAAGGGCTATAATGTATTGACTTCTCCGGGTTTGTACATTATCATTACAACCAATATTATTAAGGTTAATATTATTTTATCATTATTATTCTTGTGTTTTCAT TTGTGTGGTTATATAATAACAATGATGTGTCGCTATTGGATCATTT
heome-f.09	TTTGGGTAGATATCGTTATAAGGGCTATAATGTATTGACTTCTCCGGGTTTGTACATTATCATTACAACCAATATTATTAAGGTTAATATTATTTTATCATTATTATTGTTGTTTTCAT TTGTGTGGTTATATAATAACAATGATGTGTCGCTATTAGATCATTT
heome-f.10	TTTGGGTAGATATCGTTATAAGGGATATAATGTATTGACTTCTCCGGGTTGGTTACATTATCATTACAACCAATATTATTAAGGTTTCATATTATTTTATCATTATTATTGTTGTTTTCAT TTGTGTGGTTATATAATAACAATGATGTGTCGCTATTGGATCATTT
heome-f.11	TTTGGGTAGATATCGTTATAAGGGCTATAATGTATTGACTTCTCCGGGTTGGTTACATTATCATTACAACCAATATTATTAAGGTTAATATTATTTTACCATTATTATTGTTGTTTTCAT TTGTGTGGTTATATAATAACAATGATGTGTCGCTATTGGATCATTT
heome-f.12	TTTGGGTAGATATCGTTATAAGGGCTATAATGTATTGACTTCTCCGGGTTGGTTACATTATCATTACAACCAATATTATTAAGGTTAATATTATTTTATCATTATTATTCTTGTGTTTTCAT TTGTGTGGTTATATAATAACAATGATGTGTCGCTATTGGATCATTT

Table S11. Haplotypes detected for Heome G in 100 Cameroonian participants

Haplotype Name	Haplotype Sequence
heome-g.0	AATGATGTAGATGCTAACACAACATCCTGGTCAGTATCTTGTGATAAATAATTTAATTTGTCCGCTAAAAAATAAAGTATAATTTACAGCTTTTTTTATTTGTATCTGTGTTTTTCAT GTTATCATAAAAAGACAAAGATGAAGTAGCACATGAGTTAATTAT
heome-g.1	AATGATGTAGATGCTAATAGCACATCCTGGTCAGTATCTTGTGATAAATAATTTAATTTATCCGCTAAAAAATAAAGTATAATTTACAGCTTTTTTTATTTGTATCTGTGTTTTTCAT GTTATCATAAAAAGACAAAGATGAAGTAGCACATGAGTTAATTAT
heome-g.2	AATGATGTAGATGCTAACACAACATCCTGGTCAGTATCTTGTGATAAATAATTTAATTTATCCGCTAAAAAATAAAGTATAATTTACAGCTTTTTTTATTTGTATCTGTGTTTTTCAT GTTATCATAAAAAGACAAAGATGAAGTAGCACATGAGTTAATTAT
heome-g.3	AATGATGTAGATGCTAATAGCACATCCTGGTCAGTATCTTGTGATAAATAATTTAATTTGTCCGCTAAAAAATAAAGTATAATTTACAGCTTTTTTTATTTGTATCTGTGTTTTTCAT GTTATCATAAAAAGACAAAGATGAAGTAGCACATGAGTTAATTAT
heome-g.4	AATGATGTAGATGCTAACACAACATCCTGGTCAGTATCTTGTGATAAATAATTTAATTTGTCCGCTAAAAAATAAAGTATAATTTACAGCTTTTTTTATTTGTATCTGTGTTTTTCAT GTTATCGTAAAAAGACAAAGATGAAGTAGCACATGAGTTAATTAT
heome-g.5	AATGATGTAGATGCTAACACAACATCCTGGTCAGTATCTTGTGATAAATAATTTAATTTATCCGATAAAAAATAAAGTATAATTTACAGCTTTTTTTATTTGTATCTGTGTTTTTCAT GTTATCATAAAAAGACAAAGATGAAGTAGCACATGAGTTAATTAT

Table S12. Haplotypes detected for Heome H in 100 Cameroonian participants

Haplotype Name	Haplotype Sequence
heome-h.0	AAGAGTGTTAATAATACTACAGCTAGTAATATTTCCAAATCCAAAAATGAAC TTGAAAAAAGTCATCAGTTAATAATAATGCTTATATTAAAATTATCGAAAATTATGATTTATTATGG TTAAAATTAATAGAAATTGAATTATGTTGTAATAATAAAAAATTTAAATCCTATAA
heome-h.1	AAGAGTGTTAATAATACTACAGCTAGTCATATTTCCAAATCCAAAAATGAAC TTGAAAAAAGTCATCAGCTAATAATAATGCTTATATTAAAATAATCGAAAATTATGATTTATTATG GTAAAATTAATAGAAATTGAATTATGTTGTAATAATAAAAAATTTAAATCCTATAA
heome-h.2	AAGAGTGTTAATAATACTACAGCTAGTCATATTTCCAAATCCAAAAATGAAC TTGAAAAAAGTCATCAGCTAATAATAATGCTTATATTAAAATTATCGAAAATTATGATTTATTATG GTAAAATTAATAGAAATTGAATTATGTTGTAATAATAAAAAATTTAAATCCTATAA
heome-h.3	AAGAGTGTTAATAATACTACAGCTAGTCATATTTCCAAATCCAAAAATGAAC TTGAAAAAAGTCATCAGTTAATAATAATGCTTATATTAAAATTATCGAAAATTATGATTTATTATGG TTAAAATTAATAGAAATTGAATTATGTTGTAATAATAAAAAATTTAAATCCTATAA
heome-h.4	AAGAGTGTTAATAATACTACAGCTAGTAATATTTCCAAATCCAAAAATGAAC TTGAAAAAAGTCATCAGCTAATAATAATGCTTATATTAAAATAATCGAAAATTATGATTTATTATGG TTAAAATTAATAGAAATTGAATTATGTTGTAATAATAAAAAATTTAAATCCTATAA

Table S13. ITrue Index ordering

Primer Name	Primer Sequence
iTru5_01_A	AATGATACGGCGACCACCGAGATCTACACACCGACAAACACTCTTTCCCTA*C
iTru5_01_B	AATGATACGGCGACCACCGAGATCTACACAGTGGCAAACACTCTTTCCCTA*C
iTru5_01_C	AATGATACGGCGACCACCGAGATCTACACCACAGACTACACTCTTTCCCTA*C
iTru5_01_D	AATGATACGGCGACCACCGAGATCTACACCGACACTTACACTCTTTCCCTA*C
iTru5_01_E	AATGATACGGCGACCACCGAGATCTACACGACTTGTGACACTCTTTCCCTA*C
iTru5_01_F	AATGATACGGCGACCACCGAGATCTACACGTGAGACTACACTCTTTCCCTA*C
iTru5_01_G	AATGATACGGCGACCACCGAGATCTACACGTTCCATGACACTCTTTCCCTA*C
iTru5_01_H	AATGATACGGCGACCACCGAGATCTACACTAGCTGAGACACTCTTTCCCTA*C
iTru5_02_A	AATGATACGGCGACCACCGAGATCTACACCTTCGCAAACACTCTTTCCCTA*C
iTru5_02_B	AATGATACGGCGACCACCGAGATCTACACGTGGTATGACACTCTTTCCCTA*C
iTru5_02_C	AATGATACGGCGACCACCGAGATCTACACCACTGTAGACACTCTTTCCCTA*C
iTru5_02_D	AATGATACGGCGACCACCGAGATCTACACAGACGCTAACACTCTTTCCCTA*C
iTru5_02_E	AATGATACGGCGACCACCGAGATCTACACCAACTCCAACACTCTTTCCCTA*C
iTru5_02_F	AATGATACGGCGACCACCGAGATCTACACAACACGCTACACTCTTTCCCTA*C
iTru5_02_G	AATGATACGGCGACCACCGAGATCTACACTGGATGGTACACTCTTTCCCTA*C
iTru5_02_H	AATGATACGGCGACCACCGAGATCTACACTTCAAGCACACTCTTTCCCTA*C
iTru5_03_A	AATGATACGGCGACCACCGAGATCTACACAACACCACACACTCTTTCCCTA*C
iTru5_03_B	AATGATACGGCGACCACCGAGATCTACACTGAGCTGTACACTCTTTCCCTA*C
iTru5_03_C	AATGATACGGCGACCACCGAGATCTACACCACAGGAAACACTCTTTCCCTA*C
iTru5_03_D	AATGATACGGCGACCACCGAGATCTACACTGACAACCACACTCTTTCCCTA*C
iTru5_03_E	AATGATACGGCGACCACCGAGATCTACACTGTTCCGTACACTCTTTCCCTA*C
iTru5_03_F	AATGATACGGCGACCACCGAGATCTACACCCTAGAGAACACTCTTTCCCTA*C
iTru5_03_G	AATGATACGGCGACCACCGAGATCTACACGCATAACGACACTCTTTCCCTA*C
iTru5_03_H	AATGATACGGCGACCACCGAGATCTACACCAAGTGCTTACACTCTTTCCCTA*C
iTru5_04_A	AATGATACGGCGACCACCGAGATCTACACCGTATCTCACACTCTTTCCCTA*C
iTru5_04_B	AATGATACGGCGACCACCGAGATCTACACCGTCAAGAACACTCTTTCCCTA*C
iTru5_04_C	AATGATACGGCGACCACCGAGATCTACACCCATGAACACACTCTTTCCCTA*C
iTru5_04_D	AATGATACGGCGACCACCGAGATCTACACGGTACTTCACACTCTTTCCCTA*C
iTru5_04_E	AATGATACGGCGACCACCGAGATCTACACACCGCTATACACTCTTTCCCTA*C
iTru5_04_F	AATGATACGGCGACCACCGAGATCTACACTTCCAGGTACACTCTTTCCCTA*C
iTru5_04_G	AATGATACGGCGACCACCGAGATCTACACTCGAACCTACACTCTTTCCCTA*C
iTru5_04_H	AATGATACGGCGACCACCGAGATCTACACTAGTGCCAAACACTCTTTCCCTA*C
iTru5_05_A	AATGATACGGCGACCACCGAGATCTACACGGTACGAAACACTCTTTCCCTA*C
iTru5_05_B	AATGATACGGCGACCACCGAGATCTACACAAGCATCGACACTCTTTCCCTA*C
iTru5_05_C	AATGATACGGCGACCACCGAGATCTACACGCCAATACACACTCTTTCCCTA*C

iTru5_05_D	AATGATACGGCGACCACCGAGATCTACACCTGTATGCACACTCTTTCCCTA*C
iTru5_05_E	AATGATACGGCGACCACCGAGATCTACACCTTAGGACACACTCTTTCCCTA*C
iTru5_05_F	AATGATACGGCGACCACCGAGATCTACACTCAGCCTTACACTCTTTCCCTA*C
iTru5_05_G	AATGATACGGCGACCACCGAGATCTACACACATGCCAACACTCTTTCCCTA*C
iTru5_05_H	AATGATACGGCGACCACCGAGATCTACACGATGGAGTACACTCTTTCCCTA*C
iTru5_06_A	AATGATACGGCGACCACCGAGATCTACACCGATCGATACACTCTTTCCCTA*C
iTru5_06_B	AATGATACGGCGACCACCGAGATCTACACTACTCCAGACACTCTTTCCCTA*C
iTru5_06_C	AATGATACGGCGACCACCGAGATCTACACAGCTACCAACACTCTTTCCCTA*C
iTru5_06_D	AATGATACGGCGACCACCGAGATCTACACTCGACAAGACACTCTTTCCCTA*C
iTru5_06_E	AATGATACGGCGACCACCGAGATCTACACTATGACCGACACTCTTTCCCTA*C
iTru5_06_F	AATGATACGGCGACCACCGAGATCTACACAGCCAACACTACTCTTTCCCTA*C
iTru5_06_G	AATGATACGGCGACCACCGAGATCTACACGATCTTGACACTCTTTCCCTA*C
iTru5_06_H	AATGATACGGCGACCACCGAGATCTACACCCTCGTTAACACTCTTTCCCTA*C
iTru7_101_01	CAAGCAGAAGACGGCATAACGAGATGGTAACGTGTGACTGGAGTTCA*G
iTru7_101_02	CAAGCAGAAGACGGCATAACGAGATCAACACAGGTGACTGGAGTTCA*G
iTru7_101_03	CAAGCAGAAGACGGCATAACGAGATACACCTCAGTGACTGGAGTTCA*G
iTru7_101_04	CAAGCAGAAGACGGCATAACGAGATCATGGATCGTGACTGGAGTTCA*G
iTru7_101_05	CAAGCAGAAGACGGCATAACGAGATTGATAGGCGTGACTGGAGTTCA*G
iTru7_101_06	CAAGCAGAAGACGGCATAACGAGATCGGTTGTTGTGACTGGAGTTCA*G
iTru7_101_07	CAAGCAGAAGACGGCATAACGAGATCAACGAGTGTGACTGGAGTTCA*G
iTru7_101_08	CAAGCAGAAGACGGCATAACGAGATACCATAGGGTGACTGGAGTTCA*G
iTru7_101_09	CAAGCAGAAGACGGCATAACGAGATGGTGTACAGTGACTGGAGTTCA*G
iTru7_101_10	CAAGCAGAAGACGGCATAACGAGATCAGCATACTGACTGGAGTTCA*G
iTru7_101_11	CAAGCAGAAGACGGCATAACGAGATGGACATCAGTGACTGGAGTTCA*G
iTru7_101_12	CAAGCAGAAGACGGCATAACGAGATAGAAGGACGTGACTGGAGTTCA*G
iTru7_102_01	CAAGCAGAAGACGGCATAACGAGATCGCCTTATGTGACTGGAGTTCA*G
iTru7_102_02	CAAGCAGAAGACGGCATAACGAGATCAGGTAAGGTGACTGGAGTTCA*G
iTru7_102_03	CAAGCAGAAGACGGCATAACGAGATTTGCAACGGTGACTGGAGTTCA*G
iTru7_102_04	CAAGCAGAAGACGGCATAACGAGATGCTGAATCGTGACTGGAGTTCA*G
iTru7_102_05	CAAGCAGAAGACGGCATAACGAGATGAACGTGAGTGACTGGAGTTCA*G
iTru7_102_06	CAAGCAGAAGACGGCATAACGAGATAACGCACAGTGACTGGAGTTCA*G
iTru7_102_07	CAAGCAGAAGACGGCATAACGAGATCGCAACTAGTGACTGGAGTTCA*G
iTru7_102_08	CAAGCAGAAGACGGCATAACGAGATTGGCTCTTGTGACTGGAGTTCA*G
iTru7_102_09	CAAGCAGAAGACGGCATAACGAGATTGAGCTGTGTGACTGGAGTTCA*G
iTru7_102_10	CAAGCAGAAGACGGCATAACGAGATGCCTTAACGTGACTGGAGTTCA*G
iTru7_102_11	CAAGCAGAAGACGGCATAACGAGATTGTGGCTTGTGACTGGAGTTCA*G
iTru7_102_12	CAAGCAGAAGACGGCATAACGAGATAACCGTGTGTGACTGGAGTTCA*G

iTru7_103_01	CAAGCAGAAGACGGCATAACGAGATAATCGCTGGTGACTGGAGTTCA*G
iTru7_103_02	CAAGCAGAAGACGGCATAACGAGATGGTCACTAGTGACTGGAGTTCA*G
iTru7_103_03	CAAGCAGAAGACGGCATAACGAGATTAGTCTCGGTGACTGGAGTTCA*G
iTru7_103_04	CAAGCAGAAGACGGCATAACGAGATACCATGTCGTGACTGGAGTTCA*G
iTru7_103_05	CAAGCAGAAGACGGCATAACGAGATAGACATGCGTGACTGGAGTTCA*G
iTru7_103_06	CAAGCAGAAGACGGCATAACGAGATGATGGAGTGTGACTGGAGTTCA*G
iTru7_103_07	CAAGCAGAAGACGGCATAACGAGATCAGTCACAGTGACTGGAGTTCA*G
iTru7_103_08	CAAGCAGAAGACGGCATAACGAGATGTTCTTCGGTGACTGGAGTTCA*G
iTru7_103_09	CAAGCAGAAGACGGCATAACGAGATAAGACACCGTGACTGGAGTTCA*G
iTru7_103_10	CAAGCAGAAGACGGCATAACGAGATGCCTTCTTGACTGGAGTTCA*G
iTru7_103_11	CAAGCAGAAGACGGCATAACGAGATTCGAACCTGTGACTGGAGTTCA*G
iTru7_103_12	CAAGCAGAAGACGGCATAACGAGATGGAACATGGTGACTGGAGTTCA*G
iTru7_104_01	CAAGCAGAAGACGGCATAACGAGATTATGGCACGTGACTGGAGTTCA*G
iTru7_104_02	CAAGCAGAAGACGGCATAACGAGATCTACAAGGGTGACTGGAGTTCA*G
iTru7_104_03	CAAGCAGAAGACGGCATAACGAGATAATCCAGCGTGACTGGAGTTCA*G
iTru7_104_04	CAAGCAGAAGACGGCATAACGAGATCCTCGTTAGTGACTGGAGTTCA*G
iTru7_104_05	CAAGCAGAAGACGGCATAACGAGATGCAACCATGTGACTGGAGTTCA*G
iTru7_104_06	CAAGCAGAAGACGGCATAACGAGATGGTATAGGGTGACTGGAGTTCA*G
iTru7_104_07	CAAGCAGAAGACGGCATAACGAGATCGACCTAAGTGACTGGAGTTCA*G
iTru7_104_08	CAAGCAGAAGACGGCATAACGAGATGATCTTGCGTGACTGGAGTTCA*G
iTru7_104_09	CAAGCAGAAGACGGCATAACGAGATAAGGCTCTGTGACTGGAGTTCA*G
iTru7_104_10	CAAGCAGAAGACGGCATAACGAGATTCCATTGCGTGACTGGAGTTCA*G
iTru7_104_11	CAAGCAGAAGACGGCATAACGAGATTACTCCAGGTGACTGGAGTTCA*G
iTru7_104_12	CAAGCAGAAGACGGCATAACGAGATCGATGTTGCGTGACTGGAGTTCA*G

Supplementary Data Sheet 1. Protocol 1. *Pf*-SMARRT Amplification

Plasmodium falciparum Simultaneous Multiplex Antimalarial Resistance and Relatedness Testing v1.9

Infectious Disease Epidemiology and Ecology Lab University of North Carolina at Chapel Hill, Chapel Hill, NC, USA

www.med.unc.edu/infdis/ideel

References: PMID:29246158, PMID:15132750, PMID:15273102,
<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC11370457/>

Introduction: This is a multiplex PCR protocol for simultaneously amplifying *AMA1*, *Heome*, *CRT*, *MDR1*, *DHFR* and *K13* from Chelex-extracted Dried Blood Spot (DBS) samples.

The PCR primers were adopted from:

- *AMA1* - PMID:29246158
- *CRT* - designed for this protocol
- *MDR1* - PMID:15132750
 - *MDR184* designed for this protocol
- *DHFR* & *DHPS* - PMID:15273102
- *K13* – designed for this protocol
- *HeOME*: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC11370457/>

1. **Materials:**

1.1 List of primers and expected amplicon lengths based on lab-cultured 3D7 isolate.

Gene	Primer Name	Primer Sequence (5'-3')
<i>AMA1</i>	<i>AMA1_F</i>	CCATCAGGGAAATGTCCAGT
	<i>AMA1_R</i>	TTTCCTGCATGTCTTGAACA
<i>CRT</i>	<i>PFCRT_F_v2</i>	GGTGGAGGTTCTTGTCTTGG
	<i>PFCRT_R_v2</i>	AGTTGTGAGTTTCGGATGTTACA
<i>MDR1</i>	<i>MDR1_86_F</i>	TGTATGTGCTGTATTATCAGGAGGAAC
	<i>MDR1_86_R</i>	AATTGTACTAAACCTATAGATACTAATGATAATATTATAGG
	<i>MDR1_184_F_v2</i>	GTGAGTTCAGGAATTGGTACGA
	<i>MDR1_184_R_v2</i>	GCCTCTTCTATAATGGACATGGT
	<i>MDR1_1034_F</i>	AAAAAGAAGAATTATTGTAAATGCAGCTT
	<i>MDR1_1034_R</i>	GGATCCAAACCAATAGGCAAAA
<i>DHFR</i>	<i>DHFR_51_59_F</i>	TGAGGTTTTTAATAACTACACATTTAGAGGTCT
	<i>DHFR_51_59_R</i>	TATCATTTACATTATCCACAGTTTCTTTGTT
	<i>DHFR_108_F</i>	TGGATAATGTAAATGATATGCCTAATTCTAA
	<i>DHFR_108_R</i>	AATCTTCTTTTTTTAAGGTTCTAGACAATATAACA
<i>DHPS</i>	<i>DHPS_437_F</i>	TGAAATGATAAATGAAGGTGCTAGTGT
	<i>DHPS_437_R</i>	AATACAGGTACTACTAAATCTCTTTCATAATTTTT
	<i>DHPS_540_F</i>	AATGCATAAAAGAGGAAATCCACAT
	<i>DHPS_540_R</i>	TCGCAAATCCTAATCCAATATCAA
	<i>DHPS_581_F</i>	CCTCGTTATAGGATACTATTTGATATTGGAT
	<i>DHPS_581_R</i>	TGGGCAATAAATCTTTTTCTTGAATA
	<i>DHPS_613_F</i>	TGGATTAGGATTTGCGAAGAAAC
	<i>DHPS_613_R</i>	GTTGTGTATTTATTACAACATTTTGATCATTC
<i>K13</i>	<i>K13-A-432-466_F</i>	GAAAGTGAAGCCTTGTTGAAAGAAG
	<i>K13-A-432-466_R</i>	GTACACATACGCCAGCATTGTTG
	<i>K13-B-461-531_F2</i>	TGATGGTGTAGAATATTTAAATTCGATG
	<i>K13-B-461-531_R</i>	CTACCATTTGACGTAACACCACA

	K13-C-513-570_F	TTTGAAACTGAGGTGTATGATCG
	K13-C-513-570_R	GCTGATGATCTAGGGGTATTCAA
	K13-F-567-630_F	TGGCACCTTTGAATACCCCT
	K13-F-567-630_R	AGGTAATTAAGCTGCTCCTGA
	K13-G-660-709_F	TGGCAATTTCTAAATGGTGTACC
	K13-G-660-709_R	GCCAAGCTGCCATTCATTG
Heome	HeOME-A F	TTTTAGTTTCGGTATTTTGTGTTCCCTCTT
	HeOME-A R	AAGAAATTTATCAGAGTTACAAAAGGGAAATC
	HeOME-B F	TTTATCCTTATCATTATTTCCATCATTCTGG
	HeOME-B R	AAAAATAAAAGGAACAATGTAATGGTTGAAAA
	HeOME-C F	TCCCGAAAACATATCACTAGATCCAT
	HeOME-C R	AAAAATTAAACATGATGCCACATTTTAGTAGT
	HeOME-D F	AGCTATCATTACATGCTGACACAATAT
	HeOME-D R	GATTAGTTGTGGAGATGATAAACTATCAAATT
	HeOME-E F	AATTTTCTTATATAACCTAAGTTGATGACTTGG
	HeOME-E R	AGAACAGATGAAGTAACTACTCGATTAAATGA
	HeOME-F F	ATCTTTTTCGTTGTATGTGCATAATCA
	HeOME-F R	AAACATAATTCTAATGATATTGACCTTGTCGA
	HeOME-G F	ACATTACACAAAATAGAAAAATCTTCATTTTTC
	HeOME-G R	ATCAATAATCAAAATCATGATAACAACCAATT
	HeOME-H F	AATAACTTAAATAAAAAATATGGACGGCTCC
	HeOME-H R	GACATTCTTTCAATGCTTCCGAAA

	Chr	fullStart	fullStop	Size	He
AMA1	11	1294287	1294523	236	
DHFR-108	4	748343	748493	150	
DHFR-51-59	4	748173	748361	188	
DHPS-436-437	8	549631	549742	111	
DHPS-540	8	549958	550119	161	
DHPS-581	8	550076	550216	140	
DHPS-613	8	550102	550254	152	
Heome-A	4	115454	115684	230	0.842943
Heome-B	5	213026	213251	225	0.837785
Heome-C	14	771406	771633	227	0.789254
Heome-D	13	2545064	2545291	227	0.77364
Heome-E	8	1344617	1344845	228	0.715463
Heome-F	11	1816151	1816378	227	0.703947
Heome-G	13	815861	816091	230	0.607549
Heome-H	11	408542	408771	229	0.605662
K13-A	13	1725576	1725728	152	
K13-B	13	1725381	1725645	264	
K13-C	13	1725264	1725482	218	
K13-F	13	1725105	1725295	190	
K13-G	13	1724870	1725020	150	
MDR1-1034	5	960948	961039	91	
MDR1-184	5	958372	958599	227	
MDR1-86	5	958071	958206	135	
CRT	7	403503	403693	190	

1.2 Primer pool ratios

- Mix the primers above (**both forward and reverse**) into pools 1 and 2 in the volumes below.

Primer	Pool 1 Volume	Primer	Pool 2 Volume
AMA1	20	PFCRT_v2	12
MDR1-86	20	MDR1-184_v2	12
MDR1-1034	15	K13-B	18
K13-A	20	K13-F	12
K13-C	25	DHPS-540	24
K13-G	15	DHPS-613	18
DHPS-436-437	15	DHFR-51-59	24
DHPS-581	15	Heome-A	24
DHFR-108	25	Heome-C	18
Heome-B	25	Heome-E	18
Heome-D	20	Heome-G	36
Heome-F	15	Heome-H	12
Total Pool 1 Volume	460μL	Total Pool 2 Volume	456μL

1.3 Other reagents

- Nuclease-free water
- [QIAGEN Multiplex PCR kit](#) (Cat.206145 or Cat.206143)
- ZR-96 DNA Clean & Concentrator-5 (Cat.4024)**
- Qubit or PicoGreen supplies for quantifying product (Thermo Fisher, Cat.Q32854 or Cat.Q33232)
- IDEEL iTru library preparation protocol

1.4 Consumables/Supplies

- Nuclease-free PCR microtubes/strips/plates.
- Nuclease-free 1.5uL Eppendorf tubes (for master mixes).
- Cold plates for microtubes and 1.5uL tubes.
- Consumables/Supplies

1.5 Experiment Setup

- Thaw reagents and samples' DNA on ice.
- Clean pipettes and working bench tops with 70% ethanol
- UV-treat the working surface, all equipment, and tubes/strips.
- Briefly vortex and centrifuge each reagent before use.

2. Methods:

- Dilute the 100 μ M stock primers down to 10 μ M working solutions for each forward and reverse primer.
- Prepare the two primer pools as indicated in **section 1.1** by mixing all forward and reverse primers in the ratios shown in **section 1.2**.
- Prepare the PCR master mix based on the table below. Since there are two primer pools, the recipe below applies to both pools 1 and 2. Additionally, two or more PCR replicates are required per primer pool to resolve PCR and sequencing errors. The resulting setup has four PCR reactions per sample, i.e., two PCR replicates for primer pool 1 and two PCR replicates for primer pool 2.

Reagents for a 25µl reaction	1x Reaction
Qiagen PCR Master Mix	12.5 µL
Primer Pool Mix*	2.5 µL
Template DNA*	2.0 µL
PCR grade water	8.0 µL
Total Volume:	25.0 µL

*The “primer pool mix” represents one for pool 1 that is separate from pool 2

*Template DNA volume can be adjusted depending on the concentration of DNA

- For a sample called “001”, a simple naming convention across the two PCR replicates can be:
 - 001-p1-A & 001-p1-B** - representing the two replicates for primer pool 1.
 - 001-p2-A & 001-p2-B** - representing the two replicates for primer pool 2.
- Combine the master mix and sample DNA, seal caps/strips/plates, and place in a thermocycler using the following conditions:

Step	Temperature	Duration
Initial activation step	95°C	15 min
45 Cycles	Denaturation 94°C	30 sec
	Annealing 60°C	90 sec
	Extension 72°C	90 sec
Final Extension	72°C	10 min

- Spot check PCR products from both pools by gel electrophoresis against a 50 bp or 100 bp ladder, using a 2% agarose gel. Expected amplicon sizes are listed above.

3. PCR purification

- During the optimisation of this assay, the **ZR-96 DNA Clean & Concentrator-5 (Cat.4024)** was used. Refer to insert protocol for additional guidance.

4. Product quantification and mixing of P1 and P2

- If desired, assay the purified products using capillary electrophoresis to ensure that you have the desired PCR fragments (*approx.* 100-300bp) before proceeding to library preparation.
- For small numbers of samples, quantify the DNA using Qubit reagent and mix equal amounts of P1 and P2 product prior to library preparation.
- For larger numbers of samples, quantification by FI on a plate reader is recommended.
 - Make a dilution series of 6 controls containing DNA in TE-Buffer (0.1ng to 100ng)
 - Dilute PicoGreen reagent 1:200 according to manufacturer's protocol in TE buffer.
 - Mix 48µL diluted PicoGreen with 2 µL DNA or control in a black flat-well microplate.
 - Shimmy plate *gently* against palm to level the fluid in the bottom of the well. Cover plate if needed with optical film.
 - Incubate plate for 2 minutes at room temperature in the multimode plate reader before quantifying the fluorescence signal with 485Ex - 535Em filters.
- Use calculated concentration to mix equal ng of P1 and P2 reaction products for each sample replicate to go into library preparation

5. Library preparation

- Refer to the IDEEL iTru library preparation protocol.

Supplementary Data Sheet 2. Protocol 2. Library Preparation

iTru Protocol Adapted from Published Materials

Infectious Disease Epidemiology and Ecology Lab University of North Carolina at Chapel Hill, Chapel Hill, NC,
USA

www.med.unc.edu/infdis/ideel

Reference: <https://peerj.com/articles/7755/#supplemental-information>

1.1 Overview:

This protocol is intended for the preparation of libraries from purified *Pf*-SMARTT amplicons. The original iTru protocol followed methods for Kapa Hyper Prep Kits using half volumes. This protocol halves the volume again to give ¼ reactions. We're using NEB Modules instead of the actual Kapa Reagents referenced in Adapterama, but the same "Sera-Pure" SpeedBeads. We also omitted the dual size selection. The recipe for SpeedBead Mix is within the supplementary material of this paper.

1.2 Consumables

- a. P20 Pipette Tips
- b. P200 Pipette Tips
- c. Reagent Reservoir (sterile)
- d. 96-well PCR plates
- e. 50mL Conical Tube (sterile)
- f. Foil plate seals
- g. Optical plate seals

1.3 Reagents

Library Preparation

- a. New England Biotechnology Modules:
 - i. End Repair (NEB, Cat.E6050L)
 - ii. dA Tailing (NEB, Cat.E6053L)
 - iii. T4 Quick DNA Ligase (NEB, Cat.E6056L)
- b. 100% Anhydrous Ethanol
- c. Molecular Grade Water
- d. 1X TE Buffer
- e. Elution Buffer (Qiagen, Cat.19086)
- f. SpeedBead Mix (see Reference)
- g. KAPA HiFi HotStart ReadyMix 2X (KAPA Biosystems, Cat. KK2602)
- h. 5uM Stubby iTru Adapter (see Reference)
- i. i5 Index Primer Set (**Table S14**)
- j. i7 Index Primer Set (**Table S14**)

Library Quantification

Qubit HS dsDNA Kit (Thermo Fisher, Cat.Q32854), or
Quant-iT PicoGreen (Thermo Fisher, Cat.Q33232)

1.4 Other Supplies

DynaMag 96-well Magnetic Plate (Thermo Fisher, 12331D)

Multi-mode Plate Reader (FI Unit)
Thermocycler
P10, P20 and P200 12-tip multi-channel pipettor

2.0 Protocol

Notes:

- I. Begin from amplified, purified, and combined equimolar *Pf*-SMART pools for each sample replicate.
 - II. Allow SpeedBeads to equalize to room temperature for 30 minutes before using.
 - III. Keep all reagents and Master Mixes in a cold block or on ice.
 - IV. Do not vortex NEB enzymes. Gently flick to mix, then spin briefly.
 - V. Combine Master Mixes immediately before use and return reagents to the freezer promptly.
 - VI. A list of iTru i5 and i7 Dual Indexing Primers, as well as instructions for synthesizing Stubby Y-yoke Adapters, can be found in the original Adapterama paper within their published Supplemental Materials or an extended list.
-

2.1 End Repair

1. Insert 12.5µL of purified amplicons in 5 µL of End Repair Master Mix.
 - a. 2.5µL buffer, 1.25µL enzyme, & 1.25µL water (1/4 NEB module recommendations).
 2. Cover with foil seal, spin plate briefly, then Incubate in a Thermocycler at 20°C for 30 minutes.
 3. Add 50µL of SpeedBeads to the 17.5µL volume, mix by gently pipetting, and incubate at room temperature for 5-15 minutes.
 - Magnet capture (DynaMag 96 well magnet plate).
 - Wash twice with 80µL of 80% EtOH.
-

2.2 A-Tailing

4. Aspirate final wash with 200µL pipette tip, then re-aspirate final ethanol drop in bottom of well with 20µL tip.

If working in 96 well plates, and ethanol is aspirated entirely, the first row or column is likely dry enough for immediate addition of A-tailing Master Mix. If beads appear cracked, they have been allowed to dry for too long.

5. To dry beads, add 12.5µL A-tailing Master Mix.
 - a. 1.25µL buffer, 0.75µL enzyme*, & 10.5µL water.
 6. Incubate in a Thermocycler at 30°C for 30 minutes.
 7. Add 22.5µL of PEG/NaCl to 12.5µL volume, mix by pipetting, and incubate at room temperature for 5-15 minutes. This step re-binds DNA to beads.
 - Magnet capture & wash twice with 80µL of 80% EtOH.
-

2.3 Stubby Adapter Ligation

8. To dry beads add 1.25 µL of 5 µM iTru Y-yoke adapters. Dispense adapters directly onto bead pellets with multichannel. Do not attempt to mix with pipette.

If using the DynaMag 96-well magnet, it is easiest to add adapters in columns of 8 samples

9. Add 11.2 µL of Ligation master mix to the adapter-moistened bead. Pipette gently to mix.
 - a. 2.5 µL buffer, 1.25 µL enzyme*, & 7.5 µL water.
10. Incubate in a Thermocycler at 20°C for 15 minutes.
11. Add 12.5 µL of PEG/NaCl to the 12.5 µL volume, mix, incubate at room temperature for 5-15 minutes.
 - Magnet capture & EtOH wash twice with 80 µL of 80% EtOH.

~Here is where you could incorporate size selection if desired~

12. To dry beads add 12.5 µL of 1X TE Buffer. Incubate at room temperature for 5 minutes. This is the pre-PCR library prep.

There is no need to capture beads or aspirate supernatants. Adapter ligated sample with beads goes into dual-indexing PCR. Keep pre-PCR Library Prep at 4°C while preparing the Dual-Index PCR Plate.

2.4 Dual-Index PCR Library Amplification

13. To a new tube add 12.5 µL of KAPA Hotstart Ready Mix, 1.25 µL of i5 primer, 1.25 µL of i7 primer.
14. Mix 10 µL of the pre-PCR library prep (with speedbeads). Cover with foil seal and spin briefly.
15. Place in thermocycler for 8-12 cycles using the following conditions:

Step	Temperature	Duration
Initial activation step	98°C	45 sec
45 Cycles	Denaturation	98°C
	Annealing	60°C
	Extension	72°C
Final Extension	72°C	1 min

-
16. Use Qubit or PicoGreen to quantify the PCR product. Capture beads with a magnet and aspirate libraries in supernatant to quantify. The DNA is fully in the TE Buffer at this point. Consider running a HS tape if there is not a significant increase in DNA concentration.

17. If pooling samples for sequencing, combine libraries in equal nanogram amounts by dividing the calculated ng/ μ L value by the ng of each library desired in the pool.

Ex. To pool 200 ng of a sample with a concentration of 45ng/ μ L, divide $200/45 = 4.444\sim\mu$ L

18. After pooling libraries, purify with Speedbeads

- a. Add 32 μ L of Speedbeads to a 25 μ L volume of pooled libraries. Mix gently by pipetting and incubate at room temperature for 5 minutes.
- b. Magnet capture and wash twice with 80 μ L of 80% EtOH
- c. To dry beads add 25uL 1X TE or Elution Buffer

19. Libraries are ready to be diluted and loaded onto an Illumina sequencing platform

Data Availability: Sequencing data are available on the sequence read archive (PRJNA1171111). Metadata are available in Table S3 and upon request to Dr. Innocent M. Ali. Custom scripts are available at https://github.com/bailey-lab/pfsmarrt_cameroon_7-6-24.