

## Media for *V. natriegens*

### LB2:

Standard LB media (10.0 g/L Tryptone, 5.0 g/L Yeast Extract, 10.0 g/L NaCl) with additional 1.5% NaCl

### LBv2 Salts:

Standard LB media (10.0 g/L Tryptone, 5.0 g/L Yeast Extract, 10.0 g/L NaCl) with additional with additional v2 salts (204 mM NaCl, 4.2 mM KCl, and 23.14 mM MgCl<sub>2</sub>)

### BHlv2:

Brain heart infusion Broth (37 g/L Teknova Brain Heart Infusion Broth Dry Media (Cat. No. B9500) with additional v2 salts (204 mM NaCl, 4.2 mM KCl, and 23.14 mM MgCl<sub>2</sub>)  
note: BHI can not be autoclaved together with the salts, they have to be autoclave separate and later on added.

Matthew T Weinstock, Eric D Heseck, Christopher M Wilson, Daniel G Gibson  
*Vibrio natriegens* as a fast-growing host for molecular biology  
Nature Methods volume 13, pages 849–851 (2016)